CEREC Ortho SW 2.1

Upgrade Description (EN)

Installation notes:

- Perform a backup of all patient data before starting the installation process.
- Do not perform an installation on several acquisition units or computers at the same time if these are within the same network.
- Before starting the installation process, please check the patient database directory in the software (Configuration/Settings/Patient Database). The path should always point to a directory, never to the root of a drive, i. e. never to "D:\", but rather to "D:\Data".

Hardware Support:

CEREC Ortho SW 2.1 supports Omnicam (with guided scan) and Primescan (no guided scan required). CEREC Ortho SW 2.1 now also supports CEREC AC/AF (Connect) units of the previous generation with Omnicam. This requires an installed Windows 10 upgrade.

Changes:

- Hardware support for first-generation CEREC AC units with Windows 10 Upgrade
- Virtual Articulation functionality added to the MODEL Phase.
- Improved scanning performance
- Faster and more accurate tooth segmentation proposals
- Palatal Height and depth of curve of Spee measurement added to Model Analysis steps.
- Improved tooth positioning in Treatment Simulation

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