CEREC Tessera™

Advanced Lithium Disilicate

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Processing Guide for CEREC® Software Version 5.2.9 and upwards

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
CEREC Tessera™ Advanced Lithium Disilicate is a restorative material specifically designed to fabricate restorations using an in-office CAD/CAM workflow. CEREC Tessera™ Restorations are matrix-fired in just 4 minutes 30 seconds using the CEREC SpeedFire Furnace for increased efficiency.

Preparation Guidelines


Crowns fabricated from CEREC Tessera™ blocks may be conventionally cemented using resin-modified glass ionomer cements or may be bonded using adhesive resin cements or self-adhesive universal resin cements such as those found within the Dentsply Sirona Calibra® family of cements. To ensure optimal restoration success, follow the tooth preparation guidelines below. In all cases, inlays, onlays or veneers must be placed using a **fully adhesive bonding protocol**.

Minimum Wall Thickness (Conventional Cementation with RMGI Cements*)

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
*Taper between 4° and 8°
minimum coronal length: 4.0 mm
Round internal line angles

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
*Taper between 4° and 8°
minimum coronal length: 4.0 mm
Round internal line angles


Minimum Wall Thickness (Adhesive Bonding)

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

Inlays

Onlays

Veneers

Crowns

- Axial reduction of 1.0 to 1.5 mm.
- Wall angles 4 – 8 degrees with the long axis of the tooth.
- Centric and dynamic occlusion, reduce incisally/occlusally by 1.5 mm.
- Lingual shoulders extended at least 1.0 mm into the proximal contacts surfaces.
- Shoulder preparation without a bevel with all angles rounded.
- Smooth all surfaces.

Inlays/Onlays


- Preparation should be free of undercuts, exhibiting a draw.
- All internal line angles should be rounded.
- Minimum 1.0 mm preparation depth in the central fossa.
- Ensure margins are placed clear of occlusal contact points.
- Use adhesive bonding.

Veneers

- Reduce to 0.6 mm for the labial surface and 0.4 mm in the gingival area.
- Reduce the labiolingual incisal angle by 1.0 to 1.5 mm.
- Locate margins in enamel.
- Use chamfer or rounded-shoulder preparation.
- Locate proximal extensions proximally to conceal preparation margins and to avoid proximal gingival undercuts.
- Use adhesive bonding.

Glazing (optional) and Firing

- While the restoration is milling, pre-heat the CEREC SpeedFire furnace by pressing the “pre-heat” button (lower left section of the control screen). The oven will begin pre-heating to 400°C (taking approximately 1 minute).
- Prior to glazing, the restoration must be clean and free of oil and other surface contaminants. Clean the surface of the restoration with a steam cleaner or clean thoroughly for 30 seconds with soap and water. Any contamination after cleaning must be prevented.

Options for Glazing:

Spray or Paste Glaze only

- The Dentsply Sirona Universal Spray Glaze or Dentsply Sirona Universal Overglaze are recommended.
- Apply a thin, even coat of glaze to all sides of the restoration.
- Ensure glaze does not get on the intaglio surface by blocking out with moldable silicone or similar material.
- Air dry for 10 seconds.
- Apply a second coat if needed.

Stains and Spray Glaze

- Apply stain to restoration as desired.
- Apply a thin, even coat of glaze to all sides of the restoration.
- Ensure glaze does not get on the intaglio surface by blocking out with moldable silicone or similar material.
- Air dry for 10 seconds.
- Apply a second coat if needed.

Note: glazes and stains by third-party manufacturers may require other firing temperature and cycles than those recommended for CEREC Tessera CAD/CAM blocks

Firing Directions:

- To achieve the mechanical properties of CEREC Tessera, the firing cycle is mandatory – with or without glaze.
- For glaze firing, choose firing program “Glaze”. For no-glaze firing (dry firing), polishing of the surface is required before firing to obtain a glossy appearance of the crown. Then choose firing program Universal*.
- Ensure the restoration is centered in the CEREC SpeedFire stage.
- Firing of two (2) restorations is possible when selecting program Universal* with Multijob.

For Firing Pad Firing:

- Place the stained and/or glazed restoration on the firing pad, select the job from the control menu of the CEREC SpeedFire and press “Start”.
- If processing an anterior restoration, place the restoration with the lingual side down on the firing pad. With premolars, place the restoration with the interproximal in contact with the firing pad.
- If any imperfections exist on the restoration after firing, they can be corrected by polishing.

For Pin Supported Firing:

- **Use only Dentsply Sirona Investment Pins (REF# 5365490111).**
- Do not use other support devices or refractory putties or processing failure may occur.
- Do not place investment pins on the firing tray, as the height of the firing chamber is too low resulting in a risk of sticking of the crown in the furnace and its failure.
- Do not use investment pins with Dentsply Sirona Firing Pad.

*Designation is effective May 2024; prior to May 2024, the programs are still called Glaze Spray for glaze firing and Glaze (Stain & Paste) for no-glaze (dry) firing.

General firing parameters for third-party ovens:

CEREC Tessera™ firing programs will have to be added to specific oven programs.

General Firing Recommendation

	Start Temperature	Closing time	Preheating time	Heating rate	Final temperature	Holding time	VAC
	°C	min	min	°C/min	°C	min	min
Paint-On Glaze (only)	400	2:00	2:00	55	760	2:00	–
2nd & Subsequent Glaze Firing – if needed	400	2:00	2:00	55	760	2:00	–
Spray Glaze	400	1:00	1:00	55	760	2:00	–
Stain and Spray Glaze	400	1:00	1:00	55	760	2:00	–

Programat CS2, EP 5000/5010

	Standby Temperature B	Closing time S	Heating rate t ↑	Firing temperature T	Holding time H	Vacuum on/off	Long-term Cooling L
	°C	min	°C/min	°C	min	Vac 1 (°C) Vac 2 (°C)	°C
Paint-On Glaze (only)	403	6:00	55	760	2:00	off	0
2nd & Subsequent Glaze Firing – if needed	403	6:00	55	760	2:00	off	0
Spray Glaze	403	2:00	55	760	2:00	off	0
Stain and Spray Glaze	403	3:00	55	760	2:00	off	0


Programat CS6

	Pre-drying	Program Levels	Heating rate t	Final temp. T °C	Holding time H (min)	Quick opening	
Universal Stain & Glaze / Universal Spray Glaze	→off	1↗ 1↘	55 C°/min 70 C°/min	760 651	2:00 0:00	→OFF	V2 0 °C V1 0 °C

Note: Additional firing programs are available. Consult your local Dentsply Sirona representative for a complete list or visit www.dentsplysirona.com/CERECtessera for other furnace firing programs.

CEREC Tessera™ blocks can be adhesively bonded or conventionally cemented based on the prep guidelines.

Adhesive Cementation with Calibra® Ceram Cement and Prime&Bond Adhesive*

1. Apply and agitate Prime&Bond adhesive for 20 seconds, then evaporate solvent.** (Phosphoric Etching is optional).
2. Carefully etch the internal surfaces of the restoration with 5% hydrofluoric acid for 30 seconds, then rinse thoroughly. Apply Calibra® Silane Coupling Agent to the dried intaglio surface and leave undisturbed for at least 60 seconds. Remove excess silane with compressed air. The surface appears "moist"; this is normal.
3. Apply Calibra® Ceram Adhesive Resin Cement.
4. Seat restoration. Briefly light cure cement at the margins by moving the curing light tip in a constant motion around the margins for no more than 5 seconds per surface to initiate the gel phase. Clean up excess cement.
5. For final polymerization, light cure all areas of the restoration for 20 seconds total (buccal, lingual, and occlusal).
6. Polish the margins to ensure a smooth surface and impart a high luster.

*Prime&Bond elect®, Prime&Bond Universal™ or Prime&Bond active® Adhesive
** For restorations thicker than 2.5 mm, cure adhesive prior to seating the restoration.

Conventional Cementation

1. Clean, rinse and dry. Do not desiccate. Isolate with cotton rolls.
2. Etch intaglio surfaces for 30 seconds with 5% HF acid.
3. **Application**
Apply cement.
4. **Seating**
Seat the restoration and stabilize it under pressure, to avoid any movements, until the cement reaches a rubber-like consistency.
5. **Excess Cement Cleanup**
Remove excess material when the cement has reached the rubber-like consistency, approximately 2 min. after seating crown.

Note: Do not use any pre-treatment agents, eg. silanes, primers, or bonding agents when using Calibra® Bio cement. Use of any such pre-treatments will destroy the bond between cement and restoration.

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Programat and Indenco are not products of Dentsply Sirona and are trademarks of their respective companies.
DeguDent GmbH