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Biosil[®] f

IM05

R only

CE 0124

Made in Germany

Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso
Instrução de uso
Návod k použití
Brugsanvisningen
Οδηγία χρήσης
Használati utasítás
Naudojimo instrukcija
Lietošanas instrukcija
Bruksanvisning
Gebruiksaanwijzing
Instrukcja obsługi
Bruksanvisning

DeguDent
A Dentsply Company

Biosil® f

D	Gebrauchsanweisung	3
GB	Instructions for use	8
F	Mode d'emploi	13
I	Istruzioni per l'uso	18
E	Instrucciones para el uso	23
P	Instrução de uso	28
CZ	Návod k použití	33
DK	Brugsanvisningen	38
GR	Οδηγία χρήσης	43
H	Használati utasítás	48
LT	Naudojimo instrukcija	53
LV	Lietošanas instrukcija	58
N	Bruksanvisning	63
NL	Gebruiksaanwijzing	68
PL	Instrukcja obsługi	73
S	Bruksanvisning	78

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1 Indications for use

Biosil® f is a CoCr alloy for the fabrication of partial model cast dentures.

For dental use only.

1.1 Technical data

Composition in mass %

Co	Cr	Mo	Si	Mn	C
64.8	28.5	5.3	0.5	0.5	0.4

Biosil® f is free from nickel and beryllium and may contain Fe.

Color	white
Melting range	1320 – 1380 °C
Vickers hardness	400 HV 10
0.2 % yield strength	700 MPa
Tensile strength	900 MPa
Modulus of elasticity	220 GPa
Elongation at fracture	5 %
Density	8.4 g/cm ³
Preheating temperature	1000 °C
Casting temperature	1500 °C

It is in accordance with the international standard ISO 6871-1.

2 Contraindications

Do not use in case of hypersensitivity to one or more of the metal constituents of the alloy.

3 Warnings

In the mechanical processing of dental alloys, always work with local aspiration plus face or respiratory protection.

Cobalt - vapour is released during processing.

May cause cancer.

May release metal vapours when melting.

Do not breathe Dusts and vapours.

Ensure suitable suction/aeration at the work place and with operational machinery.

4 Precautions

Avoid occlusal and approximal contact of different alloy types.

5 Adverse reaction

Allergies to metals contained in the alloy as well as electrochemically conditioned dysesthesia are possible. In individual cases systemic side effects of metals contained in the alloy have been reported.

6 Step-by-Step instructions

6.1 Sprues

2 – 4 sprues – depending on the size of the frame work should be attached to each object. In order to facilitate flow of way and controlled solidification, wax sprue wire with a diameter of 3 – 4 mm should be attached to the thickest points of the casting **without** tapering. The length of the sprues must be determined individually by the selection and the position of the sprue former.

6.2 Investing

Apply a thin coat of Waxit® wetting agent to the wax object or spray with Waxit® and then blow it dry carefully (do not use compressed air). Invest the pattern with Optivest® or Biosint® investment. Observe the working instructions of the investment that is used. A fine investment can be used. However, it is not required if a fine-grain investment compound like e.g. Optivest® is used.

6.3 Preheating

Preheat slowly with the suction system activated. Adhere to the heating rates and holding times given by the manufacturer of the respective investment. After the final temperature of 1000 °C has been reached, preheat the casting ring for 60 min. In the case of particularly filigree constructions the preheating temperature must be raised to 1050 °C. The holding time must be extended correspondingly, if a larger number of casting rings is used.

6.4 Casting

Due to their shape the ingots can be placed easily into the ceramic crucibles. Use ceramic crucibles only.

- a) HF-casting unit (e.g. Degutron®): Do not use graphite or carbon crucibles. After a time of further heating of 5 – 7 s, depending on the quantity of alloy, casting is carried out before the oxide skin tears up.
- b) Flame-melting (e.g. Motorcast® compact): melt with the reducing flame of a propane-oxygen mixture. After exceeding the liquidus temperature and a period of continued heating – 5 – 15 s depending on the alloy quantity – begin casting.
- c) Arc (e.g. Degumat®): melt with the arc on step 3 – 4. Begin the casting process after a period of further heating – 2 – 7 s depending on the alloy quantity. Since melting in the Degumat® unit is carried out under protective gas (argon), no oxide film is formed on the melt.

6.5 Devesting

After casting let the casting ring cool down at room temperature for at least 20 min. Quenching in water is not recommended because of the risk of deformation of the casting. After devesting blast with blasting agent. In order not to destroy the passivation layer on the surface, CoCr alloys **should not be treated with stripping agent**.

6.6 Finishing

All finishing instruments that are suitable for cast dentures can be used. Use covering varnish to protect filigree elements, e.g. clasps, against abrasion in case of electro-polishing.

6.7 Soldering

Use flux Oxynon® for soldering of Biosil® f to precious metal alloys. The solder that is selected depends on the precious metal alloy to be soldered. For soldering Biosil f we recommend the Degudent®-Lot (solder) N1W.

6.8 Reusability

The fabrication of high-quality dentures requires the use of extremely pure materials. Therefore, to avoid risks only original Biosil® f ingots should be used.

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DeguDent GmbH