

IMPORTANT INFORMATION

- Not for use on roots that have immature or open/perforated apices
- Not for use with chloroform or other solvents
- Do not operate without flowing irrigant
- Stop flow of irrigant before removing needle from canal or changing canals
- Do not apply excessive force on the syringe
- Do not over-insert or bind tip into the canal
- Working times given in this technique card are for NaOCl (6%)
- Single patient use only

DENTSPLY
TULSA DENTAL
SPECIALTIES

Manufactured By:

DENTSPLY Tulsa Dental Specialties

DENTSPLY International, Inc.

608 Rolling Hills Dr.

Johnson City, TN 37604

1-800-662-1202


1-800-597-2779 (fax)

www.tulsadentalspecialties.com



Do Not Reuse

For Dental Use Only



Non-Returnable If Opened

Rx only



CONSULT
INSTRUCTIONS
FOR USE


© 2012 DENTSPLY International, Inc. PUPFTC Rev. 2 6/12

TECHNIQUE CARD



ProUltra[®] PIEZOFLOW[™]
ULTRASONIC IRRIGATION NEEDLE

DENTSPLY
TULSA DENTAL
SPECIALTIES



Use personal protective equipment and fully instrument and use a rubber dam before irrigation – Some additional orifice shaping may be required (see DFU).

- Alcohol wipe ProUltra® PiezoFlow® Ultrasonic Irrigation Needle, wrench tighten onto ultrasonic handpiece
- Set working length by trial fitting the needle into canal
 - **Set stopper 1 mm short of binding – no deeper than 75% of working length – needle should never be taken to apex**
- Attach minimum 20 ml syringe full of irrigant to PiezoFlow® Luer Lock
- Set older ProUltra® Piezo Booster to power level #5 – do not adjust
- Set ProUltra® Piezo Ultrasonic to power level #9 – do not adjust
- **Place inactive ProUltra® PiezoFlow® Ultrasonic Irrigation Needle into canal**
- **Begin irrigant delivery**, after irrigant has filled canal:
 - **Activate ultrasonic booster**
 - **Begin continuous suction** at access opening
 - **Move needle up and down** within canal – gently work in a straight line from top of canal to stopper setting
 - Irrigate for **1 minute with flow of 15 ml/min** per canal (NaOCl 6%)
- **Stop irrigant flow** before changing canals or removing needle