Duceralloy C

(D) Gebrauchsanweisung	(I) Istruzioni per l'uso	(SLO) Navodila za uporabo
(DK) Brugsanvisning	(PL) Instrukcja stosowania	(P) Instruções de uso
(GB) Instructions for use	(E) Instrucciones de uso	(CZ) Navod k použití
(FIN) Käyttöohjeet	(RUS) Правила пользования	(TR) Kullanım Talimatı
(F) Mode d'emploi	(S) Brugsanvisning	(H) Használati utasítás
(GR) Οδηγίες Χρήσης	(SK) Návod na používanie	(UAE) إرشادات الاستخدام

DeguDent
A Dentsply Company

(D)	Gebrauchsanweisung	1	(RUS)	Правила пользования	53
(DK)	Brugsanvisning	7	(S)	Brugsanvisning	59
(GB)	Instructions for use	13	(SK)	Návod na používanie	64
(FIN)	Käyttöohjeet	19	(SLO)	Navodila za uporabo	69
(F)	Mode d'emploi	24	(SLO)	Instruções de uso	74
(F)	Οδηγίες Χρήσης	30	(P)	Navod k použití	80
(GR)	Istruzioni per l'uso	36	(CZ)	Kullanım Talimatı	86
(I)	Instrukcja stosowania	42	(TR)	Használati utasítás	91
(PL)	Instrucciones de uso	47	(H)		96
(E)			(UAE)	إرشادات الاستخدام	

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CoCr Dental Casting Alloy for the Ceramic Bonding Technique

Duceralloy C is a highly corrosion-resistant non-noble ceramic-bonding alloy based on CoCr. It can be veneered with all dental ceramics (e.g. Duceram® Plus, Duceram® Kiss, Duceram® love) in the CTE range of $14.3 \mu\text{m/m} \cdot \text{K}$ (25–600 °C).

Duceralloy C meets the requirements of the EN ISO 9693 Standard.

Material Specifications: Composition:

Co	59.4 %
Cr	24.5 %
W	10.0 %
Nb	2.0 %
V	2.0 %
Si	1.0 %
Mo	1.0 %
Fe	0.1 %

Technical Data

Type:	Non-precious ceramic bonding alloy
Vickers hardness:	330 HV
0.2% yield strength:	630 MPa
Tensile strength:	790 MPa
Modulus of elasticity:	200 GPa
Elongation at fracture:	3 %
Density:	8.8 g/cm ³
Solidus point:	1,270 °C
Liquidus point:	1,370 °C
WAK value: (25 – 500 °C)	14.0 µm/m·K
(25 – 600 °C)	14.3 µm/m·K

Contraindications: Do not use in known insensitivity to one or several metals contained in the alloy.

For dental use only.

Keep under lock and out of the reach of children.

Warning! Duceralloy C contains cobalt and chromium. There are maximum allowable concentrations (MAC values) for the workplace for certain compounds containing these metals. All finishing with rotary instruments should be undertaken with a vacuum system in place. In addition, we recommend wearing a mouth guard and face shield.

This alloy should not be delivered to patients or processed by persons with a known hypersensitivity to one of its components.

Adverse effects: Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.

Interactions: Avoid occlusal and proximal contact with different alloy types.

Processing:

Modelling: Before modelling, the dies should be coated with die spacer. The crown should be modelled to reduced anatomical shape. Care should be taken in the modelling process to ensure that no sharp angles or edges are created that might later result in cracks in the ceramic veneer.

When modelling and finishing the frameworks, make sure to maintain a minimum wall thickness of 0.3 mm for single crowns and 0.5 mm for abutment crowns.

Sprueing: If single-crown sprueing is preferred to the beam-casting method, the sprues should have a diameter of 3–4 mm. When using the beam-casting method the following sprue dimensions should be observed.

Sprue from the cone	3–4 mm diameter
Transverse sprue (beam)	4–5 mm diameter
Connection to crown	3–3.5 mm in diameter and 5–8 mm in length

Investing: Any phosphate-bonded, graphite-free investment compound can be used that is suitable for alloys that melt at high temperatures (e.g. Deguvest SR). When using metallic muffle rings, an adequate – if necessary multi-layered – muffle lining is required.

Follow the Instructions for Use for the investment compound employed.

Wax elimination/preheating: At 300 °C and 950 °C depending on flask size:

Flask size	1 x 20 min
	3 x 30 min
	6 x 45 min
	9 x 60 min

Required metal quantity: The quantity of metal required for casting is calculated according to the following rule of thumb: “Wax weight of the mould multiplied by the alloy density (8.8 g/cm³) + approx. 10 g”.

Melting procedure: Duceralloy C can be melted using the high-frequency procedure, induction casting or with an open flame.

High-frequency procedure: The casting is initiated 3–5 seconds after complete melting.

Induction casting (vacuum-die casting): The casting is initiated 3–5 seconds after the oxide film is broken.

Open-flame casting: To melt Duceralloy C, the burner must be used such that the heat required for melting is uniformly distributed. Use the maximum burner setting.

Follow the Instructions for Use for the burner used.

The alloy will form an oxide skin on melting. It must be further heated until the metal under the oxide skin moves visibly due to the pressure of the flame. The casting is initiated after 3–5 seconds.

Caution! Do not use any melting/flux powder.

Crucible materials: Use only ceramic melting crucibles to cast Duceralloy C. Before casting, the crucible (without metal) should be preheated in an oven.

Finishing: Carbide cutters are suitable for finishing. Diamond tools should not be used. Always ensure that the tools are sharp. Always use high cutting speeds and low pressure.

Ceramic veneering: After removal from the mould, the surface of the Duceralloy C is roughened by intensive sandblasting with Al_2O_3 (250 μm) at a pressure of 3–4 bar.

Ceramic bonding is performed in accordance with the Instructions for Use for the ceramic veneering material.

Duceram Plus, Duceram Kiss and Duceram love ceramics are ideally suited.

Observe the ceramic manufacturers' recommendations for non-precious alloys (see "Safe veneering guideline for Duceram® Kiss/love and non-precious alloys").

Oxide firing: There is no fundamental need for oxidizing.

Polishing: To reduce oxidation on the uncovered alloy fractions, it is advisable to pre-polish them before ceramic firing. After firing, the alloy surface is first polished to a satin lustre using a rubber stone or wheel and then polished to a high gloss using polishing or diamond paste. Polishing the

metal surface will increase its corrosion resistance. It reduces the release of alloy components and further improves the biocompatibility of the alloy. Residual polish can be removed with a steam jet cleaner or in an ultrasound bath.

Soldering:

Before firing:	Degudent solder U1W (1,120 °C)
Flux:	Oxynon; DS1

Reusing old material: Neither sprues or cones nor failed castings may be reused, as this may result in oxides being embedded in the casting object, compromising its quality.

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