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CoCr Dental Casting Alloy for the Burn-On Technique

StarLoy soft is a very corrosion-stable, non-ferrous burn-on alloy based on CoCrW. It possesses excellent mechanical properties and is appropriate for use as a burn-on alloy. It can be veneered with the dental ceramics StarLight Ceram and Duceram® Kiss. Comply with the recommendations issued by the ceramics manufacturer.

Contraindications: Do not use in known hypersensitivity to one or several metals contained in the alloy.

Adverse effects: Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.

Interactions: Avoid occlusal and approximal contacts of different alloy types.

Avoid inhaling dust and vapours while in contact with dental alloys. Use suitable vacuum appliances for protection against dust and vapours. Use a facemask or respiratory protection.

For additional information concerning handling, dust protection and application notes, please consult our brochure "Precious-metal dental alloys – Processing instructions".

For dental use only.

Keep locked up and out of the reach of children.

StarLoy soft meets the requirements of the EN ISO 9693 Standard.

Material Specifications: Composition:

Co	54,1 %
Cr	20,0 %
W	16,4 %
Nb	0,2 %
Fe	7,5 %
Si	1,5 %
Mn	0,3 %

Technical Data:

Type:	NEM burn-on alloy
Vickers hardness:	280 HV10
Yield strength:	550 MPa
Tensile strength:	710 MPa
Modulus of elasticity:	200 GPa
Elongation at break:	12 %
Density:	9.1 g/cm ³
Solidus point:	1,390 °C
Liquidus temperature:	1,410 °C
WAK value: (25 – 500 °C)	14.6 µm/m·K
(25 – 600 °C)	14.9 µm/m·K

Modelling: Before modelling, the stubs should be coated with stub varnish. The shape of the crowns should represent the tooth but be reduced in size. While modelling care should be taken that no sharp corners or edges are formed. These could later result in cracks in the veneer ceramic. For reasons of mechanical stability and leaking of the casting, the crown wall thickness should not be less than 0.4 mm in metal-ceramic work.

Instructions for using modelling plastics (e.g. Pattern Resin) with metal:

When producing caps for secondary crown patterns in the cone and telescopic area, the surface quality of the primary parts has a crucial effect on the possibility of removing the plastic cap.

Here it applies that the smoother the primary metal, the more easily the plastic cap comes away.

Parallel walls respectively tapering of the primary part are prerequisite for removing the plastic cap.

Pinning: If the beam-casting method is not given preference, but instead single crown pinning is implemented, the casting channel (sprue) should have a diameter of 3.5 – 4 mm. In the beam casting method the following dimensions should be observed in the sprue systems.

Sprue from the cone	4 mm diameter
Transverse sprue (beam)	4 – 5 mm diameter
Connection to crown	3.5 mm in diameter and 5 – 8 mm in length

Investment: Any phosphate-bonded, graphite-free investment material that is suitable for alloys that melt at high temperatures can be used. If metallic muffle rings are used, an adequate, if necessary a multilayer, muffle lining is required.

For the investing process, the investment material's instructions for use are to be followed.

Wax Expulsion/Preheating: At 300 °C and 950 °C depending on the size of the casting flask:

Casting flask size	1 x 20 min
	3 x 30 min
	6 x 45 min
	9 x 60 min

Required Quantity of Metal: The quantity of metal required for casting is calculated according to the following rule of thumb: "Wax weight of the mould multiplied by the alloy density (9.1 g/cm³) + 10 g".

Melting procedure: StarLoy soft can be melted using the high-frequency procedure or the induction procedure or with an open flame.

High-frequency procedure: The casting is initiated 3–5 seconds after complete melting.

Induction casting (vacuum die casting): The casting is initiated 3–5 seconds after the oxide film is broken.

Open-flame casting: To melt StarLoy soft, the burner must be used such that the heat required for melting is uniformly distributed. The maximum burner setting is used.

Please comply with the burner manufacturer's instructions for use.

The alloy forms an oxide skin on melting. It must be further heated until the metal under the oxide skin moves visibly due to the flame pressure. Casting is initiated after 3–5 seconds.

Attention! Do not use any melting/flux powder.

Crucible materials: Only ceramic melting crucibles may be used to pour StarLoy soft. Before the casting process the crucible should be preheated in the oven.

Elaboration: Carbide milling cutters are suitable for machining. Diamond tools should not be used. In all cases ensure that the tools are sharp. In principle the material should always be machined at high cutting speed and with low pressure. The minimum wall thickness for StarLoy soft is 0.3 mm and wall should not be thinner than this. After removal from the mould, the surface of the StarLoy soft is roughened by intensive sandblasting with Al₂O₃ (250 µm) at a pressure of 3–4 bar.

Ceramic Veneer: Burning-on the ceramic is performed in accordance with the instructions for using the ceramic material used for the veneer. StarLight ceram and Duceram Kiss ceramics are ideally suitable. Note that the temperature must be increased by approx. 10 °C for dentin firing. To avoid discoloration, steam-clean the casting after each firing.

Oxide Burning: There is not a fundamental necessity for oxidising.

Cooling after Burning: It is generally advisable to allow slow cooling. Please comply with the ceramic manufacturer's recommendation on firing.

Polishing: To reduce oxidation on the uncovered alloy fractions, it is advisable to prepolysh them before the ceramic burning. After burning, the alloy surface is polished to a satin lustre and subsequently polished to high-gloss with polish or diamond paste using a rubber polisher. By polishing the metal surface its corrosion-resistance is increased. In this manner the release of alloy components is reduced and consequently the alloy's biocompatibility is again improved. The polish residues can subsequently be removed with a steam jet cleaner or in an ultrasound bath.

Soldering:

Before burning:	Degudent-Lot U1W (1120 °C)
Soldering flux:	Oxynon

Laser welding: Laser-welding wire is available in the size 0.50 mm.

Reusability: The fabrication of high-quality dentures requires the use of extremely pure materials. Therefore, to avoid risks only original StarLoy soft ingots should be used.

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