User instructions

WISIL® • WISIL® M

EN

Cobalt based dental casting alloy for removable restorations.

Nickel and beryllium free.

ISO 22674

DeguDent GmbH

Rodenbacher Chaussee 4

63457 Hanau-Wolfgang

GERMANY

Tel.: +49/6181/59-50

Made in Germany
For dental use only

Chemical composition in mass-percentage	Wisil	Wisil M	Physical and alloy properties (target values)	Wisil	Wisil M
Co	64	63.1	Type	V extra hard	V extra hard
Cr	28	28	Colour	White	White
Mo	5.1	6	Density	8,4 g/cm³	8,4 g/cm³
Mn	1.0	1.0	Vickers hardness*	390 (HV5)	410 (HV5)
Si, W, C	Trace	Trace	0.2% Yield strength*	580 MPa	650 MPa
			Elongation*	7 %	5 %
			Tensile strength*	840 MPa	900 MPa
			Modulus of elasticity*	220 GPa	230 GPa
			Solidus temperature	1,355 °C/2,470 °F	1,335 °C/2,435 °F
			Liquidus temperature	1,375 °C/2,505 °F	1,365 °C/2,490 °F
			Casting temperature	1,535 °C/2,795 °F	1,530 °C/2,786 °F
Rx only					

* as cast

Indications for use:
Type V, cobalt based dental casting alloy for the manufacturing of removable dental restorations.

Contraindications:
Do not use in case of hypersensitivity to one or more of the metal constituents of the alloy.

Warnings:
Prolonged exposure to alloy dust and/or fumes may lead to lung irritation and/or pulmonary complications. Use appropriate engineering controls to limit exposure. For excessive inhalation of dust or fumes, seek medical advice.

This alloy contains chromium. Some compounds of this element are potential carcinogens. Advise your health care provider of exposure to this element.

Adverse effects:
In exceptional cases, certain components of the alloy may cause an allergic reaction or sick feeling due to electrochemical processes. Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications.

Interactions:
In exceptional cases, occlusal or approximal contact with different alloys may cause a sick feeling due to electrochemical processes.

Processing instructions

1. Sprues
Apply 2 to 4 sprues to each wax pattern, depending on the size of the object.
Sprue diameter: 3-3,5 mm.
Connect at the pattern's thickest point expanding in delta fashion.
Arrange sprues in line with the direction of metal flow. Keep sprues short to minimize distance between cone and cast. Attach sprue former at least 1 cm above pattern.

2. Painting with liquid investment
Apply a thin, even coat of liquid investment to the pattern, avoiding bubbles. Invest model immediately, do not allow liquid investment to dry.

3. Investing
PH3 is ideally suited to reach optimal fit.
Pour into casting ring, avoiding bubbles.
Mixture ratio powder : liquid
Models (with duplicating compound) 100 g : 13 ml
(with duplicating silicone) 100 g : 16 ml
Moulds 100 g : 16 ml
hardening time: 30-45 min
[see Processing instructions PH3] Alternatively to PH3, we also offers the HARTEX binder investment material for the WISIL system.

4. Preheating
Preheat casting rings in preheating furnace according to PH3 processing instructions:
300 °C (572 °F) holding time 30 min
600 °C (1,112 °F) holding time 30 min
1,000 °C (1,832 °F) holding time 30-45 min
For plates and intricate designs preheat at 1,050 °C (1,922 °F) holding time 30-45 min

5. Required metal quantities
Simple bar: 3 ingots = 25 g
Skeletal plate: 3-4 ingots = 25-34 g
Complete plate: 5-6 ingots = 42-51 g

6. Casting
a) Motorized or spring-loaded centrifugal caster for open-flame melting
Ignite the acetylene-oxygen or propane-oxygen gas mixture. Fully open both valves and adjust the flame by throttling the intake of acetylene or propane such that blue cones approx. 3 mm long are visible at the holes of the multi-orifice head. Pressure values for the flame: acetylene = 0,7 bar, oxygen =1,9 bar or propane = 0,3 bar, oxygen = 1,5 bar. This gives a neutral flame which is absolutely essential for melting WISIL. Place WISIL ingots into preheated crucible. Heat ingots evenly until they form a molten bath by keeping the torch about 40 mm above the crucible whilst making rotating movements. As soon as the oxide shade breaks, initiate centrifugal movement.

b) Atmospheric high-frequency melting and centrifugal casting unit
Do not use graphite crucibles. As soon as ingots are molten and shaded areas have disappeared from the melt, wait 4 s and then initiate centrifugal casting.

c) Vacuum high-frequency melting and centrifugal casting unit
This unit can be used for melting and casting all WISIL alloys under vacuum conditions.
As soon as ingots are molten and shaded areas have disappeared from the melt, wait 4 s and then initiate centrifugal casting.

7. Disembedding
After casting, allow casting ring to cool for about 20 min. Devest under water.
Use sandblasting to remove any investment residue.

8. Finishing
WISIL frameworks can be easily finished using rotary instruments for CoCrMo removable partial dentures; use dust collector. Clean after finishing by blaster using corundum and then glassbeads.

9. Electrolytic polishing
Immerse dry WISIL cast into the electrolytic bath and polish for 5 min; rinse, dry and inspect. Repeat procedure if necessary to improve lustre.

10. Soldering
Flame soldering: Clean surfaces, remove oxides. Secure parts to be soldered. Cover surfaces with flux. Use CoCr solder in conjunction with soldering powder; suitable for all repair work.
Alternatively, gold solders for CoCr can be used [observe manufacturer's instructions].

11. Polishing
After electrolytic polishing, remove surface with rubber polisher.
Polish with brush and paste. Use wool buffs for high lustre.
Clean with steam jet or in an ultrasonic cleaning bath.

12. Reusability
The fabrication of high-quality dentures requires the use of extremely pure materials. Therefore, to avoid risks only original Wisil/Wisil M ingots should be used.

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