EN

Working procedure for CARRARA PdF

- 1. Indication.** Rich yellow dental alloy for metal-ceramic, inlays, crowns, bridges and telescopic technique.
 - 2. Contraindication.** If patients are known to be allergic to any of the components, the material should not be applied.
 - 3. Sprueing.** According to the indirect method: 3 mm Ø for crowns, 3,5 mm Ø for pontics, runner bar 4-4,5 mm Ø.
 - 4. Investing.** Invest in a reliable phosphate-bonded investment, such as Carrara Universal DL Investment.
 - 5. Burn out.** Preheat at 800°C/1470°F for 30-120 min. depending on the size of the muffle ring.
 - 6. Casting.** Casting temperature: 1210°C/2210°F in graphite crucible; or with the flame in a ceramic crucible with for example Elephant Proflux. Use at least 1/3 new metal at each casting. Clean the metal well before re-use, for example by means of sandblasting and steam-cleaning.
 - 7. Finishing.** Use carbide burs and/or ceramic bonded stones. The final trimming must be carried out in one direction only by means of fine crossteethed carbide burs! Sandblast with 50-125 µm aluminum oxide (Alublast), (pressure max. 3 bar).
 - 8. Cleaning.** With distilled water ultrasonic or by means of a steam cleaner.
 - 9a. Oxidation.** 3 minutes at 860°C/1560°F without vacuum.
 - 9b. Oxidoreduction.** Sandblast with aluminum oxide 125 µm (pressure max. 2 bar) or pickle for max. 1 min. in for example clean Pickling Agent. Afterwards clean again (see point 8).
 - 10. Bake.** With Carrara Interaction, cooling stage; quick.
 - 11. Pre-soldering.** Carrara 950°C/1740°F with for example Elephant Paste Flux.
 - 12. Post-soldering.** Carrara 750°C/1380°F (at 820°C/1508°F) with for example Elephant Paste Flux.
- For a reliable quality we advise to use Carrara PdF in combination with Carrara ceramic.**

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