User instructions

SUPRANIUM®

EN Nickel based dental metal-ceramic alloy
Beryllium and gallium free

ISO 22674
ISO 9693

Made in Germany
For dental use only

DeguDent GmbH
Rodenbacher Chaussee 4
63457 Hanau-Wolfgang
GERMANY
Tel.: +49/6181/59-50

Chemical composition in mass percentage	Physical and alloy properties (target values)			
Ni 61.0	Type	: IV extra hard	Density	: 8,4 g/cm ³
Cr 21.5	Colour	: Grey	Melting Range	: 1,300-1,360 °C / 2,372-2,480 °F
Mo 9.0	Vickers hardness*	: 220 (HV10)	Casting temperature	: 1,545 °C / 2,813 °F
Nb 4.0	Modulus of elasticity*	: 200 GPa	CTE 25-500 °C	: 13,9 µm/m-K
Co 2.0	Yield strength*	: 400 MPa (Rp 0.2)		
Fe 1.5	Tensile strength*	: 710 MPa		
Si, Mn < 1.0	Elongation*	: 30 %		* after firing

Indication

Nickel based dental alloy for fixed metal-to-ceramic dental restorations.

Contra-indication

Supranium contains nickel, it should not be used by or for people with a known nickel allergy, or with an allergy or sensitivity to one of the other components.

Adverse effects

Certain components of the alloy may cause an allergic reaction or sick feeling, due to electro-chemical reactions. Nickel is classified as a suspect carcinogen (R40). For processing precautions see MSDS.

Interactions

In exceptional cases occlusal or approximal contact with different alloys may cause a sick feeling due to electro-chemical processes or a slight sensation of electrical pulses, due to the difference in electric potential.

Processing instructions

1. Waxing up

To ensure problem-free casting of the sub-structure, do not model crown walls thinner than 0,4 – 0,5 mm to ensure that wall thickness after finishing is at least 0,25 mm. Form crowns and pontics in line with anatomical tooth shapes to permit uniform ceramic coating [incisal max. 2 mm, occlusal max. 1,5 mm]. Avoid sharp edges and receding areas in the pattern; the transition from metal to ceramic should be smooth. The lingual metal end should be located below the incisal zone to avoid spalling caused by unwanted pressure/tensile stresses.

2. Spruing

Direct method: Crowns and pontics must be provided with a sprue of 3 mm diameter; use sprues with reservoir (diameter 5, 6 or 7 mm) for full-cast crowns and pontics to avoid shrinkage cavities during solidification. Distance between reservoir and cast: 1 to 2 mm.

Indirect method: For casting with transverse sprue, the following dimensions apply: sprues to pattern = diam. 3 mm, transverse sprue = diam. 5 mm, feed channels to the transverse sprue = diam. 5 mm. Distance from transverse sprue to cast = 1 to 2 mm. Indirect method for pontics and full-cast crowns: particularly short distance from sprue button to cast.

Arrange sprues according to the direction of metal flow. Attach wax marking pin to indicate the position of the cast on the base of the casting ring. Sprues should be as short as possible. Do not place the object to be cast in the hottest point in the casting ring.

3. Stress-relieving the pattern

Use stress-relieving agent indicated by the manufacturer.

4. Investing

A fine-grained, phosphate-bonded, graphite and gypsum-free investment for crown and bridge work, is ideally suited to producing accurately fitting casts. We recommend, Carrara Universal Dustless [Cat. no. 5325460500].

5. Preheating

Follow the instructions from the investment manufacturer. For Carrara Universal Dustless, we recommend bench setting the moulds for one hour. Place the moulds in a cold furnace and heat them to 982 °C (1800 °F). Hold for one hour.

6. Required metal quantities

The number of SUPRANIUM® ingots required for a specific casting operation is determined by multiplying the weight of the wax pattern (g/including sprues) by the density of SUPRANIUM® (8,4 g/cm³). A SUPRANIUM® ingot weighs 6 g.

7. Insert the casting ring

Insert the heated casting ring in the centrifugal casting arm of the casting machine such that the marking point is pointing in the opposite direction of the rotational direction of the casting arm.

8. Casting

a) With motorized centrifugal caster for open-flame melting and casting

Ignite oxy-acetylene or oxy-propane mixture. Fully open both valves and adjust the flame by reducing the flow of acetylene or propane until approx. 3 mm long blue cones become visible at the multiple-orifice tip. Pressure values: acetylene = 0,7 bar, oxygen = 1,9 bar or propane = 0,3 bar, oxygen = 1,5 bar. This gives a neutral flame which is essential for melting SUPRANIUM®. Place SUPRANIUM® ingots in the preheated crucible. Keeping the torch about 40 mm above the top of the crucible and rotating the tip, heat the ingots evenly until they form a molten bath. Initiate centrifugal action immediately after molten bath has formed.

b) With atmospheric high-frequency melting and centrifugal casting unit
Do not use graphite crucibles. Initiate centrifugal action as soon as the ingots have melted and the shadow has disappeared.

c) With vacuum high-frequency melting and centrifugal unit

This type of unit can be used for melting and casting SUPRANIUM under vacuum conditions. Initiate the centrifugal action as soon as the ingots have melted and the shadow has disappeared.

9. Removal

After casting, allow the casting ring to cool down in air. Carefully remove the cast from the casting ring. Carefully remove residual investment compound. Clean the cast by blasting with aluminium oxide abrasive (90 – 125 µm). Do not blast crown margins directly.

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