Celtra® Ceram

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Thank you for choosing Celtra® Ceram veneering porcelain. For more than 100 years, Dentsply Sirona has been shaping the dental industry with innovative products designed to establish a new standard in esthetics and quality in fabricating all-ceramic dental prostheses. We are pleased to introduce this new product system with the highest level of performance our customers demand.

Indications

Celtra Ceram is suitable for veneering all-ceramic frameworks and is designed exclusively for dental use only by trained professionals. Specifically, it is a low-fusing, leucite-reinforced feldspathic ceramic optimized for veneering and characterizing Celtra Zirconia-reinforced Lithium Silicate (ZLS) frameworks produced using heat-pressing technology in a dental laboratory. With a CTE 25-500°C of $9.0 \times 10^{-6} \text{K}^{-1}$ and a firing temperature of 770°C (1st dentine), Celtra Ceram is also suitable for lithium disilicate frameworks.

- > Celtra Press Zirconia-reinforced Lithium Silicate (ZLS) frameworks:
CTE 25-500°C: $9.7 \times 10^{-6} \text{K}^{-1}$
- > Lithium disilicate frameworks: CTE 100-500°C: $10.5 \times 10^{-6} \text{K}^{-1}$

Contraindications

Only the indications listed above are suitable. The following are not acceptable for use of Celtra Ceram:

- > Celtra Duo blocks (CTE $11.6 \times 10^{-6} \text{K}^{-1}$ (25-500°C))
- > Titanium or alloy frameworks with a similar CTE
- > Alumina frameworks
- > Any other veneering ceramic system
- > Bruxism or other parafunctions
- > Insufficient occlusal distance

Warnings

Undesired side effects of these medical products are extremely rare with proper processing and use. In case of skin sensitization or rash, discontinue use and seek medical attention. Immunoreactions, such as allergies and/or local irritations (taste or oral mucosa), cannot completely be excluded on principle. For patients with hypersensitivity to any of the ingredients, this medical device must not be used or only used under strict supervision of a treating physician/dentist. Similarly, known cross-reactions or interactions of this medical product with other materials existing in the mouth must be considered by the physician/dentist.

Precautions

- > Do not inhale abrasive dusts
- > Do not ingest any of the materials (paste/powder/fluid)
- > Frameworks that show evidence of cracking or surface voids must be discarded and not processed further; do not attempt to repair any cracks during the veneering process or additional firings
- > Consider the following during firing to avoid potential for cracking:
 - Do not completely fill restorations with refractory putty
 - Use only ceramic pins, if possible
 - If wire pins are used, cover the wire with refractory putty but do not completely fill restoration
- > Frameworks below minimum thickness are not recommended
- > For minimum thickness recommendations for restorations, please consult table on page 8/10

Adverse reactions

No adverse reactions have been reported for Celtra Ceram. When working with these materials, make sure to comply with the Instructions for Use and the pertinent Safety Data Sheets (SDS).

Technical data

- > Coefficient of Thermal Expansion (CTE): $9.0 \pm 0.5 \times 10^{-6} \text{K}^{-1}$ (25-500°C), $T_g = 520^\circ\text{C}$
- > Device classification: Dental ceramics, type 1, class 1 per ISO 6872: 2015
- > Flexural strength: 108 MPa
- > Chemical solubility: $28 \mu\text{g}/\text{cm}^2$

Handling, transport and storage

- > Powders: Avoid direct exposure to sunlight and high heat, protect from moisture, when not in use keep lid tightly closed

Symbols on product labels

	Product code or reorder number
	Batch or Lot number
	Expiration date
	Keep dry
	Date of manufacture
	Consult instructions for use

Compatible materials

Multiple compatible materials and systems are recommended for use with the Celtra system and referenced herein. For more information related to these products, including ordering information, please visit celtra-dentsplysirona.com

Compatible liquids

For best results, the following accessory liquids are recommended:

- > Dentsply Sirona Modeling Liquid DU
- > Dentsply Sirona Modeling Liquid U
- > Ducera Liquid SD
- > Dentsply Sirona Stain and Glaze Liquid

If longer working times are desired, the following liquid can be used as well:

- > Ducera Liquid Form

Isolation

- > Ducera Isolating Fluid SEP
- > Dentsply Sirona Die Release

System Overview

Please remember when selecting shades that Celtra Press ingots are made to the corresponding shade and will match the tooth shade exactly. The pressed restoration will therefore correspond to the dentin shade and will require characterization only in the incisal area.

Cut-Back Technique

Enamels (E)

E1-Extra Light, E2-Light, E3-Medium

Enamels are optimized in translucency and color range to perfectly mimic natural effects in the incisal region. **6**, **8**, **9**, **12** & **13**

Enamel Transparent (E4)

Transparent powder without opalescence that can be used in a thin application to create depth and can be mixed with any other powders to increase translucency. **14**

Enamel White (E5)

Whitish effect powder to enhance occlusal cusps, palatal/lingual ridges or the incisal of anterior areas - can be diluted by using Enamel Opal Transparent EO4.

Enamels Opal (EO)

E01-Extra Light, E02-Light, E03-Medium

Opalescent enamel inspired by nature to expand design options in the incisal area without affecting core color. Enamel Opal increases vitality without being too transparent (no greying effect) yet maintains a transparent visual character. A true multifunctional material to quickly achieve high esthetic results.

Enamel Opal Transparent (EO4)

A strong opalescent, almost transparent multi-functional powder with broad application within the system. EO4 can be used in pure form and/or mixed with all other powders to increase opalescence.

Enamel Opal HT (EO5)

A translucent enamel to enhance accents of yellow and orange opalescent effects. Especially formulated for highly translucent core materials, yet versatile to encourage creativity.

Enamel Opal LT (EO6)

A opalescent powder with higher value and higher opacity than EO5 intended for use with lighter and bleach shade substructures, yet similar to EO5 in versatility to promote creativity.

Enamel Effect (EE)

Enamel Effects of varying hues that can be applied in incisal areas to enhance color depth and introduce natural features of adjacent dentition. **2**, **3**, **4**, **7**, **10**, **11** & **13**

Enamel Effect Sunrise (EE1) / Enamel Effect Sunset (EE3)

Opalescent effect powder can be used for yellow/orange (Sunrise) **7** & **13**, as well as orange/reddish (Sunset) accents in the enamel areas. They are well suited for increasing the chroma level at the 2nd or 3rd dentine firings. Sunset is mainly used for A-shades, whereas Sunrise is used mainly for B-shades. Both powders can be diluted using EO4 Enamel Opal Transparent to reduce the intensity of yellow and reddish colors.

Enamel Effect Violet (EE2) / Enamel Effect Sky (EE5)

Opalescent effect powders for discreet accent rod formations **2** & **11** and for creating strong illusions of depth and translucency in incisal areas. They can be diluted using EO4 Enamel Opal Transparent to reduce the intensity of the violet or sky.

Enamel Effect Fog (EE4)

Opalescent effect powder for grayish incisal areas **4** & **11** – can be diluted by using EO4 Enamel Opal Transparent to reduce the gray intensity.

Enamel Effect Ivory (EE6)

Whitish opalescent effect powder **3** & **10** for palatal/lingual ridges in the anterior area and to enhance occlusal cusps in the posterior region as well – can be diluted by using EO4 Enamel Opal Transparent to reduce the intensity of milky-ivory color.

Layering Technique

Dentins (D)

Dentins are optically balanced in hue, chroma and value, with out-of-the-bottle shade accuracy per the Vita® Classical standard. All 16 A-D shade designations¹ are available in individual bottles. Additionally, bleach shades are available in the BL1-4 standard.

Power Dentins (PD)

Power Dentins are highly chromatic, fluorescent and intensive powders for individual chromatic adjustments. The Power Dentins can be used in all areas to increase the chroma in the cervical, palatal or occlusal fossa areas, and are especially useful for lithium disilicate frameworks lacking chroma. The powders can be used pure or in varying mixtures depending upon desired effect and intensity.

Dentin Effects (DE)

Dentins are intensive powders in a range of cool and warm hues which can be used as rod formations, accents and chromatic or value adjustments plus varied effects. EO4 Enamel Opal Transparent can be used to dilute and reduce the intensity.

Dentin Gingiva (DG)

Gingiva powders in multiple shades mimicking soft tissue areas. Dentin Gingiva has same translucency as dentin and is fired during dentin applications at 770°C, which is especially helpful in veneering implant-supported frameworks.

Add-On Correction (C) & Add-On Gingiva (G)

Add-On Correction porcelain can be used for making final adjustments. Similarly, Add-On Gingiva can be used for final corrections in gingival areas. Both are low fusing ceramics with a firing temperature of 750°C or 760°C (see firing charts), and can be used either with or after glaze firing.

¹The A-D designations correspond to Vita® Classical Shade Guide which is meant to be a guide not a match.

Vita is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co.

Stains & Overglaze

Dentsply Sirona Universal Stains and Glaze is especially developed with broad compatibility with a wide array of substructures and veneering systems, including Celtra Press and Ceram restorations. More information about this universal system is available at celtra-dentsplysirona.com.

Note: For prescribed staining recipes of full contour Celtra restorations in all 16 Vita® Classical shades, please see the Celtra Staining Guide at celtra-dentsplysirona.com.

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***Power Dentin (mixing) example:** Shade B2 requires 3 parts PD1 (3X) and 1 part PD4 (1X) mixed together and applied where needed. Please consider the Dentsply Sirona Portioner to improve mixing ratio accuracy.

Framework preparation

Celtra Press frameworks should be pressed in anatomical form to support the ceramic being applied to ensure an even layered thickness of the veneering material. Fine corrections, e.g. adjustment in the marginal area, are carried out by using fine diamonds or silicone wheels.

Celtra Press – minimum framework wall thickness / veneer thickness (mm)

Technology	Area		Inlays	Onlays / Tabletops	Veneers	Anterior crowns	Posterior crowns	Anterior bridges	Posterior bridges
								Connector cross-section 16 mm ²	
Staining technique	Framework wall thickness (fully contoured)	Full-arch	1,0 ≥ isthmus width	1,5	0,6	1,2	1,5	1,2	1,5
		incisal/occlusal	1,5	1,5	0,6	1,5	1,5	1,5	1,5
Cut-back	Framework wall thickness	Full-arch	-	-	0,6	1,2	1,5	1,2	1,5
		incisal/occlusal	-	-	0,4	0,8	0,8	0,8	0,8
	Veneering (thickness)	-	-	0,4	0,7	0,7	0,7	0,7	
Layering technique	Framework wall thickness	full-arch/incisal/occlusal	-	-	-	0,8	0,8	0,8	0,8
	Veneering (thickness)	-	-	-	0,4-0,7	0,7	0,7	0,7	

*Connector for 3-unit bridge should have a minimum of a cross-sectional area of 16 mm²

PowerFire (Celtra Press frameworks only)

Use 50-micron sized aluminum oxide at 20 psi pressure and lightly blast the exterior surface of the ceramic restoration. Care must be taken not to harm the margins.

Use a steam cleaner to clean the surfaces or put the restoration in distilled water and place in an ultrasonic cleaner for 10 minutes.

PowerFire is a firing program that is carried out before the first ceramic firing of the veneering porcelain. PowerFire increases the flexural strength of the Celtra Press restoration to >500 MPa. After PowerFire, **1** blasting must be avoided as it will reduce the strength of the restoration.

Helpful Tips:

- To avoid tensions in the Celtra Ceram veneering ceramic, it is necessary to round off any sharp corners and edges of the frame.
- Irrigation will generally not be required when working on the material with diamond cutters.
- Water cooling is recommended, however, when using high-speed rotary instruments.
- To avoid overheating of the framework material, do not use a high grinding pressure.

Composite dies

The **Dentsply Sirona Die Material** is designed to mimic the actual shade of the patient's prepared tooth. When this material is placed inside the pressed Celtra crown, it will assist in accurate shade reproductions.

The dentist should take a prep shade of the tooth being restored for the laboratories reference with the Dentsply Sirona Prep Guide.

If the dentist did not take the shade of the prepared tooth, the composite die material guide below may be used to verify the final shade. Select the proper composite die material from the appropriate table.

- 1) Apply the Dentsply Sirona Prosthetics Die Release to the inside of the ceramic restoration and allow it to dry.
- 2) Place Dentsply Sirona Die Release in the inside of the restoration, then, a small amount of the Dentsply Sirona composite die material on the inside of the restoration. Pack the material to remove any voids. Immediately push a dowel pin into the uncured composite die material. Remove any excess composite from the margin area.
- 3) Light cure the composite for 1-2 minutes using a hand held light curing unit or the Triad 2000 curing unit from Dentsply Sirona.
- 4) Remove the composite die material from the restoration and carefully clean it using a steam cleaner or in distilled water in an ultrasonic cleaner for 10 minutes.

Due to the high translucency of Celtra Press, the influence of the die shade on the shade of the restoration must be taken into account. The aesthetic result is also influenced by the color of the adhesive material. Using the supplied light-curing die material, the dental technician has the ability to map the shade information supplied by the dentist to a control die to replicate information about the oral situation in the shade reproduction. The aim is to simulate the shade of the prepared tooth (follow the working instructions).

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Celtra Ceram processing

Note: When firing a Celtra® Press All-Ceramic restoration it is important to **use only ceramic/porcelain** type pins/pegs or to place the object directly on the firing pad, to prevent issues with Celtra® Press during porcelain and glaze firings. When other than recommended type of pins/pegs are used, while initial results with some firing pins may appear acceptable, internal stress can compromise long term success. Do not fill the entire restoration with refractory putty. Doing so could cause cracking of restoration.

Framework-porcelain ratio standards for Celtra Press restorations.*

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*The thickness of the veneering must not exceed 2.0 mm at any point.

- > The minimum framework wall thickness will always be based on the total thickness of the restoration.
- > The thickness ratio of the framework wall to the ceramic layer must be at least 1 : 1 to ensure framework stability and esthetics.

Cut-back Technique

The cut-back framework is complemented in the incisal or occlusal areas using incisal materials. This produces high-quality esthetic restorations in only a few steps. The exact tooth shade is ensured after applying the incisal material.

Additional dentine materials need not be applied to Celtra Press restorations.

Layering Technique

By applying the individual layering technique, one can use the broad range of dentin and enamel powders **2**, **3**, **4**, **10** & **11** giving you the ability to create natural and highly sophisticated restorations.

If the structure requires dentin proceed with the build-up. Aesthetic effects can be achieved in the enamel area using enamel opal **5**, **11** & **15** and enamel effect powders **3**, **4**, **10** & **11**. For cervical, occlusal fossa and body areas use EE1 Sunrise **7** & **13** or EE3 Sunset to enhance chromatic effects. Opal effect powders can be applied as an overlay or as needed **5**, **11** & **15** to further enhance individualization effects.

Three unit bridge after first firing **16**.

Staining / Overglazing

The Dentsply Sirona Universal Stain and Glaze system (available separately) is recommended for finishing full contour and veneered restorations. This system is especially matched to the Celtra Press framework and Celtra Ceram veneering system.

Place a small amount of the enamel stain or glaze on the palette. If necessary, mix the stain and glaze with the Dentsply Sirona Stain and Glaze Liquid to achieve a creamy consistency and apply the mixture to the porcelain surface.

Additional Universal Stains may be applied for individual characterizations **17** & **18**.

Remove the die material from the restoration. To ensure proper fit, remove excess glaze from the interior of the crown as well as the internal margin areas.

Fire the crown according to recommended firing cycles **19**.

Note: If a higher sheen is desirable, either raise the high firing temperature 10°C or use an additional 30 second hold time at the high temperature.

General firing recommendations - Cut-back and Layering Technique

PowerFire: Framework only

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
0:00	1:00	400	1:00	55	760	0	0	0:00	2:00	0:00

First firing: Dentin & Enamel

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	770	400	770	1:00	1:00	5:00

Second firing: Dentin & Enamel

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	760	400	760	1:00	1:00	5:00

Glaze Firing

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	750	0	0	0	2:00	5:00

Add-on (with and after glaze firing)

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	750	400	750	1:00	1:00	5:00

Full Contour Technique

Power firing incl. Glaze (Celtra Press framework only) - 1st firing

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	760	0	0	0	2:00	5:00

Glaze - 2nd firing

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	750	0	0	0	2:00	5:00

Add-on with 1st glaze firing

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	760	400	760	1:00	1:00	5:00

Add-on after glaze firing

Drying	Closing	Start Temp	Pre-heating	Heating Rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Hold Time	Hold Time	Cool
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	750	400	750	1:00	1:00	5:00

Note:

1. Slow cooling is mandatory; this includes correction firings of restorations after try-in.
- 2 Firing temperatures must be adapted to the number of units fired in the same cycle.
 - a. 5 to 9 units require an increase by 5 °C to 10 °C;
 - b. 10 or more units require an increase by 10 °C to 20 °C.

The values indicated here are recommended values and serve only as a guidance. Deviations of the firing results are possible. The firing results depend on the respective output per oven and are due to the manufacturer and to age. Therefore the recommended values must be individually adapted at each firing. We recommend a test fire to control the oven. All indications have been carefully elaborated and tested by us, but they are passed on without any guarantee.

For up-to-date firing recommendations please visit celtra-dentsplysirona.com.

Cementation

Preparation of the Celtra restoration

- > Clean the restoration with a steam cleaner, in ultrasonic bath or with alcohol.
- > Apply 5% - 9% hydrofluoric acid etching gel (available separately, see manufacturer's complete Directions for Use) to the interior of the restoration only and allow to soak for 30 seconds.
- > CAUTION: Follow manufacturer's precautions. Do not allow tissue or eyes to come into contact with the acid! Protective eye glasses and gloves should be worn.
- > Remove the hydrofluoric acid as per the manufacturer's instructions.
- > Dry the restoration in an air stream. It is recommended to silanize the etched surfaces immediately.
- > At chairside, apply silane only to those surfaces required for adhesive cementing.
- > Allow to soak for 60 seconds. If the silane layer is no longer liquid, add more silane. Blow-dry in a powerful air stream. (Recommended material: Calibra® Silane Coupling Agent, available separately, see complete Directions for Use).

Cementing

Depending on the indication for Celtra Press restorations a self-adhesive or full-adhesive cementation can be chosen. Crowns can be cemented with Ketac™ Cem* glass-ionomer cement. Compatible time-proven adhesive cementing materials are available as part of the Dentsply Sirona range of products. Cements are available separately.

	Self-adhesive	Fully adhesive	Glass-ionomer ¹
Inlays	R	HR	-
Onlays	R	HR	-
Veneers	-	HR	-
Crowns	HR	HR	R
Bridge	HR	HR	R

* Ketac™ is a registered trademark of 3M ESPE

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DENTSPLY International Inc.
York, PA 17401 U.S.A.
DeguDent GmbH