For better dentistry

DENTSPLY
PROTHETICS


Manufacturer:
DENTSPLY International Inc.
570 West College Ave.
York, PA 17405-0872, USA
Tel: 717-699-4190
Fax: 717-849-4543
www.ceramco.com or www.dentsply.com

CE REP
DeguDent GmbH
Postfach 1364
63403 Hanau
Germany
+49/6181/5950

Distributed by:
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3
905-851-6060

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R^{only}
CE
0086

Photograph courtesy of :
Functional Esthetics, Farmington, MO
Restoration designed by: Douglas J. Frye C.D.T.

For better dentistry

DENTSPLY
PROTHETICS

radico[®]
provisional & diagnostic resin

Directions For Use
Mode d'emploi
Istruzioni per l'uso
Gebrauchsanleitung
Instrucciones de uso

provisional & diagnostic resin

Directions for Use

Caution: Federal law restricts this device to sale by or on the order of a Dentist.

DESCRIPTION

The Radica Provisional System includes visible light-cured, resin-based composites and accessories used to create provisional restorations for patients awaiting final prosthodontic restoration.

INDICATIONS FOR USE:

1. Indicated for fabrication of all provisional dental restorations for up to 6 months.
2. Indicated for fabrication of the following provisional dental restorations (without metal) for greater than 6 months up to 2 years.
 - Inlays/Onlays
 - Veneers
 - Crowns
 - Bridges through the second bicuspid (1 pontic)
3. Indicated for the build up of a tooth on denture base materials.

CONTRAINDICATIONS:

1. Radica is contraindicated for patients with a history of sensitivity to acrylate resins.
2. Radica is not indicated for use as a fixed definitive restoration.

WARNINGS:

1. Radica contains components which may cause eye and skin irritation. May cause skin sensitization (allergic contact dermatitis) after prolonged contact. Wash thoroughly with soap and water after direct skin contact. If dermatitis or other symptoms persist, seek medical advice.
2. Sealer (liquid) is highly flammable and is an eye and skin irritant. Do not allow contact of liquid sealer with eyes due to risk of serious damage. Irritating to respiratory tract and skin. May cause allergic skin reaction. May cause sensitization after repeated skin contact, especially in susceptible individuals. Use in a ventilated workspace to avoid inhalation. High vapor concentrations of sealant/stain vapors in a non-ventilated space can induce headache, drowsiness, and irritation of eyes and respiratory system. Excessive long-term exposure may be associated with other more serious health effects. Monitor air quality per OSHA standards.

First aid for:

- **Eye Contact:** Flush eyes promptly with copious amounts of water for 15 minutes, and consult a physician.
- **Inhalation:** Move subject to fresh air. Give oxygen or artificial air as required.
- **Ingestion:** If conscious and alert, rinse mouth and drink 2-4 cupfuls of milk or water. Do not induce vomiting. Get medical attention.
- **Skin Contact:** Wash skin with plenty of soap and water for at least 15 minutes. Get medical aid if symptoms persist.

PRECAUTIONS:

1. Use the Radica System in a well-ventilated workspace.
2. This product is light sensitive! Minimize exposure to light and keep the lid on the container to maintain its effective use.
3. Store in a cool, dry area; 60°- 80°F (16° - 27°C) is recommended.
4. Keep sealer bottles tightly closed when not in use. Material is volatile and flammable. Keep away from heat, sparks and flames.
5. Avoid inhaling grinding dust. Use suction and adequate industrial ventilation.
6. Restorations with mesial-distal spans greater than 38mm cannot be processed in the Triad® 2000 curing unit. Restoration/model with a height greater than 40 mm should not be processed in the Enterra®VLC Curing Unit.
7. When making adjustments to repair or refinish restorations by a dentist intraorally, be sure to fully remove the cured sealer layer. Repolish using manual polishing techniques after repair.

ADVERSE REACTIONS:

The Radica System contains components which may cause skin dryness, sensitization (allergic contact dermatitis) or other allergic reactions in susceptible persons.


EQUIPMENT:

Radica Syringe Heater

The Radica resin must be heated prior to application. Plug in the Syringe Heater and turn on the unit. The heater is preset to 140°- 147°F (60°-64°C). Allow the heater to warm up for about 30 minutes. Approximately 10 minutes prior to use, place the necessary dentin and enamel syringes in the heater slots. The Radica syringes may remain in the heater during the day, but should not remain heated for an extended period of time.

Electric Wax Pencil

The electric wax pencil should be used for manipulating the resin. The special nickel-plated tips were designed to prevent any discoloration of the resin. Do not use other electric wax pencils unless the tips are chrome or nickel-plated steel. Ensure the wax tips remain clean. Do not use hot instruments that have been warmed with an external heat source.

Light Curing Units

Be sure to use the Enterra® VLC Curing Unit, Triad® 2000 Curing Unit (maximum of 38mm span in the Triad 2000), or the Eclipse® Processing Unit.

To ensure an optimum cure, follow these guidelines:

1. When using the Enterra, the model + restoration should not exceed 40mm in height as indicated by the perforated maximum height line on the Enterra unit.

Enterra Unit

*Perforated
Maximum
Height Line*

2. When using the Triad 2000 VLC Unit, appliance must be placed on the center of the turntable. Position the turntable height so the bottom edge of the device lines up with the lowest edge of the Lamp Access Door.

Triad 2000

*Triad 2000
Note: center and proper
height for curing*

INSTRUCTIONS:

MATRIX PREPARATION

1. Ensure that the teeth that are to be prepared are clean and free of any debris.
2. Fill in any edentulous areas requiring a provisional with a denture tooth.
3. Use the Radica Matrix Putty and mix the base and the catalyst in a one-to-one ratio until thoroughly mixed.

Filling in the edentulous space

4. Apply the putty to the area of the model requiring the provisional restorations. Be sure to also capture the adjacent teeth, and extend the putty down at least 2 mm past the marginal area.
5. Allow the Matrix Putty to harden. Remove it from the model.
6. Using a sharp instrument, remove the Matrix Putty to within 2 mm of the margin area.

Applying the Matrix Putty


Trimming the Matrix Putty

MODEL/DIE PREPARATION

Note: Proper die reduction is important. Over reduction will result in a provisional that is too thick, and will need to be hollow ground prior to insertion. Minimum connector size is 4mm x 5mm.

1. Select a Midwest® HP 1558 carbide bur that has a diameter of 1.0 mm. Make depth cuts in the die stone as a guide.
2. Use a die saw or a diamond disk and open the adjacent contacts 1 mm. Be careful not to damage the adjacent tooth. Use the carbide bur to reduce the required stone die by approximately 1.0mm.

Midwest 1558 bur with a diameter of 1.0 mm

Do depth cuts of 1.0 mm

Use a saw on the interproximals

1.0 mm reduction of the crowns

3. Apply a thin coat of the Radica Model Separator to the prepared areas of the model.

Applying the Model Separator

RADICA® DENTIN & ENAMEL SHADE CHART

Dentin Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4
Enamel Shade	Light	Light	Light	Light	Light	Light	Light	Light	Light
Dentin Shade	C1	C2	C3	C4	D2	D3	D4	i2	RD-001
Enamel Shade	Light	Light	Light	Medium	Light	Light	Light	Light	Light

PROVISIONAL FABRICATION

1. Express a small amount of the pre-heated enamel resin into the incisal area of the matrix. Use the electric spatula to spread out the resin and remove any excess. Excess resin will prevent the matrix from fully seating. Allow the resin to air cool about one minute.
2. There are four effect dentins: red-orange, yellow-orange, deep red tissue tint and tissue tint. They can be used to increase the chroma of a tooth in a specific area such as the gingival or occlusal or fill in for tissue voids around the gingival.
3. Quickly express enough dentin to fill the remaining voids in the matrix. Immediately place the matrix onto the prepared model, allowing the excess material to escape around the margins.
4. Allow the material to set approximately 2-3 minutes.
5. Remove the matrix putty from the model. If necessary, additions may be made by using the electric spatula and applying some enamel or dentin resin.

Placing the enamel

Smoothing out the enamel

Effect dentin placement

Filling with dentin

Placing the matrix

Removing the matrix

6. Excess resin may be easily removed using an instrument. Connector should be U-shaped (not V-shaped or pointed) and the size should be at least 4 X 5 mm.
7. Use a brush to apply a thin layer of the sealer. Do not point sealer onto the model while applying it to the Radica crown or bridge. If overlapped onto the model, the restoration will be difficult to remove from the model after curing.
Note: Immediately clean the brush with acetone, or place the brush in a dark area.
8. Following the placement guidelines, put the restoration/model into the Enterra, Triad 2000, or the Eclipse Processing Unit, and do the INITIAL CURE according to Table 1.

Adding on to the resin

Clean up the excess

Excess resin removed

Applying the sealer

Curing the resin

9. For crowns, proceed to FINISHING THE PROVISIONAL, STEP 1.
10. After curing, allow the model/restoration to return to room temperature before removing the cured restorations. This may be accelerated by placing the model in cold/ice water.
11. Use a brush to apply a **thin layer** of the sealer to the bottom of the pontic.
12. Place the bridge, upside down, on the Enterra or Eclipse Processing Unit turntable or on the center of the Triad 2000 turntable and cure according to the PONTIC CURE step on Table 1.

*Removing the bridge
from the model*

Applying the sealer

FINISHING THE PROVISIONAL

1. Use an instrument to remove the provisional from the model and do any additional contouring.
2. To clean the provisional, blast all internal and external surfaces with aluminum oxide at 20 psi. Steam clean all surfaces; then dry.
3. There are five stains: brown, blue, white, orange, and olive. If stains are required, place a small amount of the stains on a palette and mix with a small amount of the sealer PRIOR to applying sealer to the restoration.
4. Use a brush to apply a thin layer of the sealer to the external surfaces. Immediately clean the brush with acetone. If stains are going to be used, IMMEDIATELY apply stains to the provisional using a fine brush.
5. Cure the sealer in the Enterra, Eclipse Processing Unit or the Triad 2000 unit according to the SEALER CURE on Table 1.
6. Radica provisionals may also be mechanically polished using buffing wheels and abrasives.

Blast with aluminum oxide

Mixing the stains

Applying the stains

Applying the sealer

Completed provisionals

REPAIRS / ADD-ONS / TISSUE TINT

1. Roughen the area intended for repair, add-on or tissue tint using a carbide bur. Extend the roughening 1 mm beyond the area of repair, add-on or tissue tint.
2. If necessary, apply Model Separator to the repair area of the model and place the restoration back on the model.
3. Wet the roughened surface by brushing a thin layer of sealer onto the prepared resin surface.
4. Extrude a small amount of dentin or enamel resin (as needed) onto a palette; then use the electric spatula to apply the resin to the repair/add-on area. For the tissue tint, use the electric spatula to apply the resin to the desired area.
5. Allow to cool for 2 minutes.

- 6. Use a brush to apply a thin layer of the sealer to the repair/add-on. Immediately clean the brush with acetone, or place the brush in a dark area.
- 7. Cure in the Enterra, Triad 2000 or the Eclipse Processing Unit using the INITIAL CURE on Table 1.
- 8. Finish according to previous techniques.

Roughen the surface

Apply the sealer

Add the resin and cure

RADICA DIAGNOSTIC

Radica can also be used as a diagnostic and provisional material at the same time. This will allow for a very esthetic diagnostic presentation, and the case will already be prepared for provisionals.

- 1. Prepare the model using a uniform 1.0mm reduction.
- 2. Wax the case to full contour, optimizing the contours.
- 3. Make a mould using the Matrix Putty.
- 4. Fabricate the Radica provisional as previously indicated. Follow all steps as listed in the Provisional Fabrication section.
- 5. Cure the diagnostic restoration before delivery to your account.

Working Model

1.0 mm depth cuts

Uniform 1.0 mm reduction

Waxed to full contour

Matrix Putty
8

*Completed
Diagnostic/Provisional*

BUILD UP OF A TOOTH ON A DENTURE BASE

NOTE: Radica can be bonded to conventional acrylics or Eclipse resin using the following procedure. **To prevent any distortion of the denture base resin, be sure to have the appliance on the model when curing.**

1. After the denture base is fully cured, use a course carbide bur to roughen the surface. Steam clean or scrub under running water to remove any debris. Allow to dry thoroughly.

2. Apply a **thin coat** of the VLC Sealer. Allow to bench set for 2 minutes. Do not cure.

3. Apply the Radica to the desired area using the electric wax pencil. Allow to fully cool.

4. Apply a **thin coat** of the VLC Sealer and cure for 5 min (Initial) in the Enterra unit, or 10 min in the Triad 2000, or 3 min (Initial) in the Eclipse Processing Unit.

5. Finish and polish the denture and Radica restoration using a rag wheel with pumice and polishing compound.

DELIVERING THE RADICA RESTORATION

Ensure all internal areas are clean. If necessary blast with aluminum oxide at 20 psi. Also clean model to remove all traces of the separator before reseating for delivery to dentist.

Prior to cleaning

After cleaning

TABLE 1. CURE TIMES FOR DENTSPLY CURING UNITS

	Enterra® Curing Unit	Triad® 2000 Curing Unit	Eclipse® Processing Unit
Initial cure	5 min	10 min	3 min
Pontic cure	1.5 min	2 min	1/2 min
VLC Sealer	2 min	2 min	1 min

Note: Do not use early versions of the Triad Curing Unit.

INSTRUCTIONS FOR THE ENTERRA CURING UNIT OPERATION:


START UP

Plugging in the unit lights the display window. The number of the software version will be illuminated. The machine begins self testing before use. This happens very quickly. The name "Eclipse" appears automatically.

PRODUCTION SELECTION

Choose the F3 button on the right under the display—this button scrolls through the available programs. **When "Radica" appears**, press the center F2 button. This locks in the menu for the Radica programs.

PROGRAM SELECTION

Either of the outside buttons will continue to scroll through the Radica program selections.

Choose **the "5:00" minute program** for the **initial cure**.

Place work inside and close the door manually.

Press the "Start" button.

Use **the "1.5" minute program** for the **pontic cure**.

Use **the "2:00" minute program** for the **sealer cure**.

TIPS, TRICKS & TROUBLESHOOTING

<i>Can a handheld tack light be used with the Radica Provisional?</i>	• Yes – but only to "tack" the resin into position. You still must follow with a full curing cycle.
<i>Cannot remove the cured provisional from the model.</i>	• Did not use enough model separator. • Undercuts on the model. • Sealer was overlapped onto the model during application.
<i>Gray streaks/discoloration in the completed restoration.</i>	• Used an electric wax pencil with tips that were not nickel or chrome coated. • Used a warm instrument from an external heat source.
<i>Voids or porosity in my cured resin.</i>	• Air was trapped in the resin when expressing the resin into the silicone matrix. • Surface was not properly prepared prior to doing an add-on procedure.
<i>Shades are off.</i>	• Too much or too little enamel used in relation to the dentin. • Used the wrong enamel. • Used excess effect dentin. • Improper positioning of the model in the curing unit when curing (Triad 2000 only). • Did not use an approved curing unit. • Dentist used the wrong reline resin.
<i>Cannot express the resin from the syringe.</i>	• The syringe was not in the syringe warmer long enough. • The syringe was not fully inserted in the syringe warmer. • The tip of the syringe was left off and the material polymerized.
<i>Final provisional did not appear to have the proper strength.</i>	• Did not use the correct curing cycle. • Did not use the recommended curing unit(s). • Did not use the recommended connector size. • Improper positioning of the model in the curing unit. • Walls of the restoration are too thin. • Did not use enough separator on model.
<i>Surface of the provisional is rough after curing.</i>	• VLC sealer started to cure prior to applying to the surface. • Surface not smooth prior to the sealer/glaze application. • Application brush was contaminated. • Continued brush strokes after the sealer started to dry.
<i>Surface has lost its detail.</i>	• VLC Sealer was applied too thick.
<i>Observed blue tint on restoration due to Matrix Putty.</i>	• Insufficient mixing of base/catalyst.
<i>What is the best choice of Reline Material?</i>	• Choice of reline material and cement will affect longevity. Best results are also obtained with frequent recalls and supervision through the treatment sequence. • Best results are also obtained with more Radica and less reline. Preferably all Radica. • Reline Materials: Triad® VLC Provisional Material, Biolon® C&B Resin (DENTSPLY Trubyte), Temporary Bridge Resin, Integrity® Temporary Crown and Bridge Material.

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