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Instructions for use
Mode d'emploi
Istruzioni per l'uso
Instrucciones para el uso

DENTSPLY

ExtraStar

D	Gebrauchsanweisung	4
GB	Instructions for use.	8
F	Mode d'emploi	12
I	Istruzioni per l'uso	16
E	Instrucciones para el uso	20

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ExtraStar is a phosphate-bonded, rapid burnout investment, for pressable ceramics. ExtraStar has no plaster or graphite. ExtraStar is mixed with ExtraStar liquid.

1 Safety Tips

Do not mix with other liquids or materials during use, because this can lead to health risks or impair the quality of the pressed object.

- This product contains quartz and cristobalite, some in the form of fine dust that may penetrate the pulmonary alveoli. Inhalation of the dust constitutes a potential carcinogenic hazard. Ensure that an adequate suction line/ventilation is available in the workspace and near the machines.
- Inhalation the dust of the product may trigger silicosis.
- Employers must provide Class P2 respiratory protection.
- Heating the investment compound above 200°C will result in the formation of ammonia. Ammonia irritates the respiratory organs.
- May cause mechanical irritation of the skin and mucosa of the eye and respiratory passages.
- Irritation and caustic reactions to the skin and mucosa must be expected on local contact.

1.1 Information for the processing of mixing liquids for dental investment compounds

Following contact with skin:

Upon skin contact, wash with plenty of water.

Following eye contact:

Rinse thoroughly with plenty of water keeping eyelid open.

In case of persistent discomfort: Consult an ophthalmologist.

Mixing liquids are alkaline.

Further information is contained in the safety data sheet.

2 Processing Recommendations

Processing Parameters	Recommended
Material Temperature	21 °C / 70 °F
Mixing Ratio	100 g : 25 ml
Vacuum Mixing Time	1 min
Working Time	5 – 7 min

3 Investment storage

Most of you are probably aware of the problems which can occur when investment is subjected to variations in environment temperature (e.g. summer and winter) and the effect this can have on your castings. That is why it is very important to find the right place to store your supply of investment. The investment should be stored in a cool and dry place, between +5 °C to +30 °C at 30 % to 90 % relative humidity. If you have a place in your laboratory where these conditions can be found throughout the year, then that's the place to store your investment!

If it is not possible for you to store your investment at this temperature, then you should at least make sure that the work bench where investing is done fits these requirements (use a temperature control cupboard or similar).

Please note the symbols on product labels:

REF Product code

LOT Batch or lot number

 Use before

 Consult instructions for use

Avoid any nasty surprises by observing the following:

- Regularly check the temperature and humidity in the embedding region - **application temperature: +20 °C to +23 °C at 30 % to 90 % relative humidity.**
- Avoid using any investment that has just been brought out of an especially warm or cold place. (e.g. storage rooms, silos, cupboards etc.)!
- Fill up stocks at your workbench before your supply runs out. Give the investment time to adjust to the environment!
- Avoid using investment that has just been delivered to the laboratory. Again, give the investment time to acquire the right temperature.
- If space allows, stock up on investment (especially LIQUID) before the cold winter months arrive. This way, you forego the risk of exposing the liquid to frost during delivery

4 Expansion Control

ExtraStar is mixed with ExtraStar liquid. Expansion can be controlled by diluting this liquid with distilled water, thus achieving an optimal fit for any geometry.

Specific concentration *recommendations* to be individually adapted:

– Incisor restoration:	approx. 80 % → 20: 5 (ml)
– Premolar restoration:	approx. 75 % → 19: 6 (ml)
– Inlay, triple surface:	approx. 85 % → 21: 4 (ml)
– Inlay, double surface:	approx. 90 % → 22: 3 (ml)
– Inlay, single surface:	approx. 90 % → 22: 3 (ml)
– Veneer:	approx. 85 % → 21: 4 (ml)

Please note that many processing parameters (e.g. temperature, modelling material etc.) can influence the fit and deviations from the aforementioned recommendations are therefore possible.

5 Mixing the ExtraStar Investment

Premix powder and liquid in a clean mixing bowl of a mixing unit for about 30 seconds (including period of sprinkling by hand) with a clean spatula. Let the investment stand under vacuum for about 30 sec.

Next, mix under vacuum for about 60 sec.

When mixing ExtraStar Investment, please observe the following:

The mixing ratio and mixing time are important for consistent results. We recommend that the mixing ratios and times provided be complied with exactly. High vacuum in the mixing bowl determines quality; only start the mixing procedure after waiting for a full vacuum display.

5.1 Investing

Lubricate carefully base former and ring leveller with vaseline. Place ring former on sprue base former. While gently vibrating ring assembly, pour the vacuum-mixed investment into the ring, allowing the investment to flow bubble free into and over wax patterns. Place the ring leveller device into the ring former. The excess investment will be forced up and through the opening in the leveller device. Placement of the leveller insures correct ring height, and proper alignment of the plunger and the ring. After **10 minutes**, the ring may be removed from the ring former. Be sure to remove the uneven nub from end of ring with a knife, making sure you have a flat and even surface !

6 Burnout

ExtraStar is suitable for so-called rapid heating techniques. In order to achieve consistent results in this heating procedure, please observe the following:

After **15 min.**, measured from the moment the powder and liquid are mixed together, the ring must be placed immediately into the furnace preheated to **850 °C**. Lean the ring against the back or side walls of the oven at a **45°** angle.

In order to achieve even heating of the ring, the holding time of **60 min.** must be adhered to at the final temperature.

7 Pressing

Pressing is done in accordance with the processing instructions of the particular press ceramic manufacturer.

8 Devesting

After pressing, let the ring cool in the open air at room temperature. Stand the cooled ring on its end on the plunger. Place a plunger next to it and lay the tip of a pencil over the edge slightly. Rotate the ring, engaging the pencil lead to the side of the ring, scoring a line all around. Use a cut-off disk at the pencil line, and cut off plunger from ring. Also, cut the investment off of the plunger.

Sandblast the ring with **4 bar** pressure until the units become visible, then drop the pressure to about **2 bars**. Use only 50 micron glass bead grit (**50 µm**).

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