Lucitone[®] Clear

Dental Resin

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Resina dentale

Dentalkunststoff

Resina dental

Directions for Use

Mode d'emploi

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Indications for Use:

LUCITONE CLEAR Dental Resin is a powder and liquid system based on polymethyl methacrylate, intended for use in the fabrication of clear dental appliances such as occlusal splints.

Contraindications:

1. LUCITONE CLEAR Dental Resin is contraindicated for patients and users with a history of allergic reaction to methyl methacrylate monomer.

Warnings:

1. Wear protective gloves and protective clothing. Wash thoroughly with soap and water after contact. If dermatitis or other symptoms persist, seek medical advice.
2. Avoid inhalation or ingestion. High vapor concentration can induce headache, nausea, and irritation of eyes and respiratory system. Work in a well ventilated area. Liquid contact with eyes may cause possible corneal damage. Wear eye protection. Excessive long-term exposure may be associated with other more serious health effects. Monitor air quality per OSHA standards.

Inhalation: Move subject to fresh air. Give oxygen or artificial respiration as required.

Ingestion: Contact your regional Poison Control Center immediately.

Eye Contact: Flush eyes promptly with copious amounts of water for 15 minutes; consult a physician. Wash skin with soap and water.

Precautions:

1. When grinding prosthodontic resins, proper ventilation, masks, and vacuum systems should be used.
2. Store at 60°-80°F (16-27°C), away from moisture and direct sunlight. The liquid contains unsaturated monomers which may polymerize prematurely if stored at excessively high temperature and/or in sunlight. The powder contains an organic polymerization initiator which may degrade if stored at excessively high temperature.
3. Lucitone Liquid and LUCITONE CLEAR Dental Resin uncured resin are hazardous materials. Dispose of in accordance with Federal, State, and Local regulations.
4. Lucitone Liquid contains methyl methacrylate monomer, a flammable liquid with a flash point of 50°F (10°C). Keep away from heat, sparks, and open flame.
5. Use liquid in a well-ventilated area. Replace cap when not in use.
6. Do not pack denture resin after work time has elapsed or material has become rubbery.
7. **DO NOT USE** the long cure cycle [9 hours at 163°F (73°C)] for processing LUCITONE CLEAR Dental Resin. The long cure cycle will not produce optimum clarity of LUCITONE CLEAR Dental Resin.

Adverse Reactions:

1. Corneal damage, headache, nausea, and vomiting may occur with exposure to methyl methacrylate monomer.

(See Warnings and Precautions statements.)

appropriate engineering controls are not used.

- Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals. Residual monomer in fully cured materials can be minimized by soaking the cured prosthesis in warm water for several days.
- Particulates will be generated when grinding acrylic resins. Eye, skin, and respiratory irritation may occur if

STEP BY STEP INSTRUCTIONS:

PROPERTIES AND RECOMMENDED PROCEDURES

The following table describes typical properties of LUCITONE CLEAR Dental Resin, plus recommended temperatures and other important processing information. Further detail of this information is found later in the text.

ISO 20795-1.....	Type I Class 1
Storage temperature for powder and liquid.....	60°-80°F (16°-27°C)
Powder/Liquid ratio.....	21 g (32 cc)/10 ml
Mixing time (time to wet all particles).....	15-30 seconds
Time to reach packing plasticity @ 73° ±2°F (23° ±1°C).....	9±2 minutes
Working time.....	8±4 minutes
Material used to prepare mould.....	gypsum
Temperature of mould when packing.....	approx. 110°F (43°C)
Recommended cure time and temperature	
1st stage.....	1 ½ hours @ 163°F (73°C)
2nd stage.....	½ hour @ 212°F (100°C)
Method of cooling flask, time and temperature	
1st stage.....	½ hour in air @ 60°-80°F (16°-27°C)
2nd stage.....	¼ hour in water @ 60°-80°F (16°-27°C)

FLASKING

Use conventional dental compression moulding or the Success® Injection System, and gypsum materials for flasking.

Use clean boiling water for the final flush. DO NOT USE WAX SOLVENTS. Be sure case is free of wax and grease.

WAX ELIMINATION

Soften wax in boiling water for approximately 6 minutes. Separate flask and remove wax by flushing with a solution of boiling water and detergent.

TINFOIL SUBSTITUTE

Apply AL-COTE® Separator liberally to areas of the warm mould [approx. 120°F (49°C)] that will contact the resin. Remove excess AL-COTE film with a dry brush. (AL-COTE Separator does not form a film on teeth.) Allow AL-COTE Separator

to dry thoroughly (approx. 3 minutes following application) before packing the resin.

IMPORTANT

Measure powder and liquid accurately. Do not whip air into mix by excessive spatulation.

MIXING

To assure uniformity, shake powder jar (with a rocking motion) to prevent concentration of small particles at the bottom of jar. Add 1 unit of powder (21 gm or 32 cc) to 10 ml LUCITONE Liquid. Stir sufficiently (15-30 seconds) to assure wetting of all powder particles. Do not over-mix. Cover mixing jar and allow material to reach packing consistency [approx. 9 minutes at room temperature of $73^{\circ} \pm 2^{\circ}\text{F}$ ($23^{\circ} \pm 1^{\circ}\text{C}$)].

PACKING

LUCITONE CLEAR Dental Resin should be packed in a warm flask [approx. 110°F (43°C)]. The correct packing consistency is identified when the resin dough is free from stickiness but not rubbery. Remove resin dough from jar and condense with finger pressure into mould or place into Success Injection Cartridge.

TRIAL PACKING IN COMPRESSION PACK CASES

Place assembled flask in bench press and apply pressure slowly. If trial pack indicates insufficient resin in some areas, remove cellophane and add resin. Replace cellophane and repeat trial packing. Make certain there is a slight excess of resin. **CLOSE FLASK SLOWLY.** Open flask, remove cellophane, and trim off excess

resin. Apply AL-COTE film after trial packing has been completed. Be sure AL-COTE film is thoroughly dry prior to final closure. Then place assembled flask in bench press for final closure. The closed flask should remain under pressure for approximately 30 minutes before the heat-cure procedure is started.

SUCCESS® INJECTION PROCESS

Follow directions for the Success Injection System. Inject for 1 minute.

WORK TIME

The work time of LUCITONE CLEAR Dental Resin is approximately 8 minutes. Never pack or inject after the work time has elapsed or the material has become rubbery.

CURING

Submerge closed flask (locked by compress, spring clamps or Success Flask Brackets) in water at $163^{\circ} \pm 2^{\circ}\text{F}$ ($73^{\circ} \pm 1^{\circ}\text{C}$) for 1 ½ hours. Follow by ½ hour in boiling water. A periodic check of water bath temperature with an accurate thermometer is recommended.

IMPORTANT: DO NOT USE THE LONG CURE CYCLE [9 HOURS AT 163°F (73°C)] FOR PROCESSING LUCITONE CLEAR DENTAL RESIN. THE LONG CURE CYCLE WILL NOT PRODUCE OPTIMUM CLARITY OF LUCITONE CLEAR DENTAL RESIN.

COOLING

Bench cool the flask at room temperature for 30 minutes. Then immerse in cool water [$60^{\circ}\text{-}80^{\circ}\text{F}$ ($16^{\circ}\text{-}27^{\circ}\text{C}$)] for 15 minutes before deflasking.

FINISHING AND POLISHING

LUCITONE CLEAR Dental Resin should be finished and polished in the same manner as denture base resins.

SHADE RECOMMENDATIONS	
Dentsply Sirona REPAIR MATERIAL	DENTURE BASE RESINS
Pink	LUCITONE® (pink) and HY-PRO LUCITONE® Plain
Pink Fibered	HY-PRO® LUCITONE® Fibered Characterized LUCITONE® Lt. Red, Pink, Red, and others
Light Fibered	Characterized LUCITONE® Light HY-PRO LUCITONE® Fibered Light and LUCITONE 199®

Marca Registrada - Chile, Guatemala, Mexico. The Trademark under which this product is sold is registered in various countries.

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DENTSPLY International Inc.
York, PA 17401 U.S.A.
EC REP **DeguDent GmbH**