VITALLIUM®

V.R. INVESTMENT SYSTEM
DIRECTIONS FOR USE

The V.R. Investment System has been especially designed for use when casting Vitallium® Alloy and Vitallium® 2 Alloy Removable Partial Denture frameworks and Vitallium Special Alloy Subperiosteal Implants.

The V.R. Investment System consists of the following.
V.R. Investment: (N022065) High purity silica based refractory investment powder.
V.R. Binder: (N023001) Ethyl silicate binder
V.R. Binder-Thinner: (N024001) Binder system catalyst
Perflex® Duplicating Material: (N021055) Extremely accurate thermally reversible hydrocolloid duplicating material.
Silflex™ III Duplicating Material: (N021501) Silicone rubber duplicating material for precise and accurate model duplication.
All Purpose Protective Coat: (N025925) Refractory powder slurry (wash coat) for improved surface finish on the casting.
Refractory Model Dip: (N010540) Synthetic Resin which is used to harden the refractory model (must be melted before use.)
Cold Model Dip: (N010450) Refractory Model hardener which is used at room temperature.

Technical Data:	
Setting Expansion	0%
Thermal Expansion	1.8%*

*After complete burnout cycle at 2150°F (1177°C)
Compressive strength: 2.3 MPa (335psi)
INSTRUCTIONS FOR V.R. INVESTMENT
Keep container closed when not in use. Store in a cool dry area.

MIXING AND STORAGE INSTRUCTIONS FOR V.R. BINDER-THINNER
Prepare only the amount of liquid which will be used within three days.
Accurately measure the proper proportion of each liquid and add it to a clean and dry glass or plastic container.
NOTE: The container should be capable of being tightly sealed during storage.

Standard Mixing Ratio:**
1 part (750ml) V.R. Binder
1.2 parts (900ml) V.R. Binder-Thinner

****If a too rapid** set time is experienced under conditions of high temperatures and or humidity, the following ratio is recommended:
750ml V.R. Binder
800ml V.R. Binder-Thinner

Thoroughly mix the solution with a magnetic stirrer for 10 minutes. Allow any vapors which are produced during the mixing to escape. As an alternative, the solution can be mixed by shaking it by hand for approximately 15 minutes, until the mixture is a transparent pink (**non-cloudy**) liquid. Before use, store the prepared binder liquid in a refrigerator set to 40°F (4°C) for at least 3 hours.

The shelf-life of the mixed liquid is **72 hours**. Liquid older than 72 hours should not be used. Only clean glass or plastic containers should be used for storing liquid in the refrigerator. DO NOT mix old and new liquids together. AVOID metal contact when storing the liquid. Also, do not use containers that are contaminated with crystals of old liquid.

*For detailed processing instructions consult the Vitallium Alloy Processing Technique Manual.

MODEL DUPLICATION
1. Duplicate the surveyed and relieved master model using Perflex Duplicating Material (ReOrder# N021055). Follow label instructions for proper melt-down, dispensing and chilling procedures.
2. For optimal surface smoothness on the refractory model, use All Purpose Protective Coat (ReOrder# N025925). Shake or mix well before using. Pour a quantity of Protective Coat into the cavity of the Perflex Mold and rotate the mold to uniformly coat all surfaces. Invert the mold and shake or vibrate out any excess to prevent puddling. When the Protective Coat loses its glossy appearance immediately pour the refractory model. Do not allow the Protective Coat to dry to the state where it begins to crack. For additional information, consult the instruction pamphlet which accompanies the Protective Coat.
3. Accurately weigh and measure the appropriate amount of investment and Binder liquid.
V.R. Investment/Liquid ratio for one refractory model:
300 grams of V.R. Investment and 56 ml prepared Binder Liquid.
Material for no more than 4 refractory models should be mixed at one time.
Using a clean dry bowl, thoroughly mix investment powder and investment liquid. Pour the mix into the Perflex mold, being careful not to entrap air bubbles. Transfer case to a vibrating table and allow the mold to vibrate until the investment is set (approximately 15-20 minutes). Remove the case from the vibrating table and allow it to bench set for an additional 20-25 minutes.
4. Remove the refractory model gently from the Perflex mold and trim it on a dry model trimmer. After trimming, the base of the model should be at least 12mm (1/2-inch) thick, and the perimeter of the model should be tapered and be free of any sharp corners.
5. Place the refractory model in a cold oven and heat it to 451°F (232°C) and soak for 45 minutes. Submerge the refractory model in Refractory Model Dip (ReOrder# N010540) that has been heated to a temperature of 325°F (163°C) for 8-10 seconds. Remove the refractory model and allow the excess model dip to drain off. Allow the refractory model to cool to room temperature.
NOTE: REFRACTORY MODEL DIP MUST BE MELTED AND USED UNDER A FUME HOOD.

WAX-UP AND SPRUING
1. Wax-up case on refractory model using Flexseal Patterns.
2. Attach sprues and prepare case for investing by applying a coat of All Purpose Protective Coat according to directions.

INVESTING
1. Wax the prepared refractory model to a base using Amber Wax (N015740).
2. Prepare an investing ring which will allow at least 12mm (1/2-inch) clearance around and above the refractory model. The following rings are recommended. The ring or ring lever should extend at least 35-40mm (1-1/2 inches) above the highest point of the model to allow for the proper tapping and setting of the V.R. Investment.
a) Preformed Foil Lined Investment Rings (ReOrder# N123902 and N125104)
b) Split-Metal Rings (ReOrder# N123210, N123410 and N123610) with Asbestos-free ring liner (ReOrder# N027001).
NOTE: Dip asbestos-free ring liner in amber wax before investing.
c) Disposable rings (ReOrder# N123800).
NOTE: Seal with amber wax before investing.
DO NOT USE SOLID METAL RINGS.
3. Accurately weigh and measure the appropriate amount of Investment and Binder Liquid.
V.R. Investment/Liquid Ratio for Investing: 600 grams V.R. Investment and 112ml Prepared Binder Liquid.
For best results mix only enough investment for two cases at one time.
Mix investment and Binder Liquid thoroughly and pour the mix, carefully filling ring. Transfer the case to the vibrating table and allow the mold to vibrate until the investment begins to set (approximately 15-20 minutes). Then remove the case from the vibrating table and allow it to completely harden by bench setting for an **additional hour**.
4. Grind down the excess investment from the top of the case (rings) using a model trimmer. Leave at least 12mm (1/2-inch) of investment above the highest point of the refractory model.

BURNOUT
1. Place the rings in a cool furnace and heat to 1200°F (649°C) and hold for 30 minutes.
2. Increase the furnace temperature to 2150°F (1177°C) and soak for 30-60 minutes.
3. Cast the cases following the instructions provided with each bag of Vitallium, Vitallium 2 and Vitallium Special Alloys.
4. Allow the cases to bench cool and then proceed with divesting and finishing operations.

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