Prisma® Gloss™ & Prisma® Gloss™ Extra Fine Composite Polishing Pastes

DIRECTIONS FOR USE – ENGLISH

Caution: U.S. federal law restricts this device to sale by or on the order of a dentist.

Prisma® Gloss™ Composite Polishing Pastes are water-soluble, aluminum oxide pastes designed for final polishing of anterior and posterior composite restorations, diastema closures, esthetic incisal lengthenings or veneers. When applied with the Prisma Gloss Foam Polishing Cup, Prisma Gloss Polishing Paste improves the smoothness of the finished surface of microfill composites to a final, high luster surface, or pre-polishes by smoothing finished surfaces of hybrid and micro-hybrid composites to low luster. For hybrid and micro-hybrid composites, following with Prisma Gloss Extra Fine Polishing Paste results in a final, high luster polished surface.

COMPOSITION

Glycerin; Dispersed Alumina; Dispersed Silicon Dioxide

INDICATION

Final polishing of direct or indirect anterior and posterior composite resin restorative materials.

CONTRAINDICATIONS

None known.

Safety Alert Symbol

This is the safety alert symbol. It is used to alert you to potential injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

WARNINGS

1. Prisma Gloss Polishing Paste contains aluminum oxide and fumed silica particulates in glycerin base. Avoid prolonged or repeated contact with skin and eyes. Do not take internally.
Eye Contact: Prisma Gloss Polishing Paste contains aluminum oxide and fumed silica particulates in glycerin base which may be temporarily irritating to eyes. Before using this product wear protective glasses as well as covering the patient's eyes to protect from splashing material. In case of contact with eyes, rinse immediately with plenty of water.
Skin Contact: Prisma Gloss Polishing Paste contains aluminum oxide and fumed silica particulates in glycerin base which can cause transient skin sensitization (allergic contact dermatitis) in susceptible individuals. If contact with skin occurs wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention.
Oral mucosa contact: Avoid contact with oral soft tissues. If accidental contact occurs, flush mucosa with plenty of water and expectorate water. If sensitization of mucosa persists, seek medical attention immediately.
2. Prisma Gloss Polishing Paste should not be used with patients who have a history of severe allergic reaction to any of the components.

PRECAUTIONS

1. This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and is the sole responsibility of the practitioner.
2. Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
3. The Prisma Gloss Polishing Paste syringes should be tightly closed with the original caps immediately after use.
4. Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.

STORAGE

Prisma Gloss Polishing Pastes should be stored in the original syringe only, kept out of direct sunlight, and stored in a well ventilated place at temperatures between 10-24°C/50-75°F. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

ADVERSE REACTIONS

Product may irritate the eyes and skin. Eye contact: transient irritation. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous Membranes: inflammation, edema, sloughing. (See Warnings)

INTERACTIONS

Do not use Prisma Gloss Extra Fine Polishing Paste on microfilled composites.

STEP-BY-STEP INSTRUCTIONS FOR USE

1. Place composite restoration as per standard techniques. Complete gross reduction of excess and general outline form of the restoration using carbide finishing burs or diamond finishing instruments. Additional finishing and polishing is obtained by use of Enhance® Discs, Cups or Points. Alternatively, other standard aluminum oxide disc series may be used (see complete Directions for Use of finishing instruments). Prior to proceeding to polishing with Prisma Gloss

Composite Polishing Pastes, the restoration should be finished to a final contour. Surface should be smooth and defect free.

2. Attach an Prisma Gloss Foam Polishing Cup to the provided mandrel by inserting the bayonet end of the mandrel into the opening on the narrow end of the cup. Rotating the mandrel 1/4 turn makes insertion easier. Make sure the mandrel inserts fully into the cup. Insert the mandrel into a conventional speed contra-angled handpiece. **NOTE: PRISMA GLOSS POLISHING CUPS MUST NOT BE USED IN HIGH SPEED HANDPIECES.**
3. Dispense a portion of Prisma Gloss Polishing Paste onto a clean paper pad or dampen dish. Recap syringe promptly.
4. Apply a small amount of Prisma Gloss Polishing Paste to the surface of the Prisma Gloss Cup. Without running the handpiece, apply material to the surface to be polished. Running slowly at first, using overlapping, continuous circular motions, work the surface of the restoration, dry initially, up to moderate speed and pressure. Use flat end and corner edge of cup.
5. To increase surface luster, add water in small amounts (i.e. dropwise) to dilute the paste using a light circular buffing action. Repeat as needed to produce a smooth surface.
6. Rinse Prisma Gloss Polishing Paste from the tooth surface and Prisma Gloss Cup. **This is the final step for microfilled composites.** Continue with the next step for hybrid and micro-hybrid composite materials.
7. Dispense and apply Prisma Gloss Extra Fine Paste to the Prisma Gloss Cup as described in Steps 3 and 4 and polish the surfaces dry at first, then adding increasing amounts of water for 15-30 seconds at higher speed, with a lighter buffing pressure for final luster. The Prisma Gloss Cup should be discarded after use.
8. Clean excess material from oral surfaces with air/water spray and/or instruct patient to rinse and expectorate.

Technique Tips:

- As with any rotary instrument, heat will build up when using the Prisma Gloss Foam Polishing Cup with prolonged contact. Use with intermittent pressure.
- Use the flat end and corner edge of the Prisma Gloss Polishing Cup. Excessive lateral pressure may dislodge cup from mandrel.

DISINFECTION/STERILIZATION

Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

The mandrel for the Prisma Gloss Foam Polishing Cup exposed to splatter or spray of body fluids or that may have been touched by contaminated hands, or oral tissues, should be steam autoclaved following autoclave manufacturer's instructions. Use of other sterilization procedures/methods is at the discretion and sole responsibility of the practitioner.

The Prisma Gloss Polishing Paste syringes, if exposed spatter or spray of body fluids or that may have been touched by contaminated hands, or oral tissues, should be disinfected with a hospital-level disinfectant. Acceptable disinfectants are EPA-registered as tuberculocidal. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternaries are approved disinfectants. Disinfect the syringes by spraying the syringes with any recommended hospital-level surface disinfectant for the contact time recommended by the disinfectant manufacturer for optimum results. Spraying with glutaraldehyde is not recommended. Some phenolic-based agents and iodophor-based products may cause surface staining. Agents containing organic solvents, such as alcohol, may tend to dissolve the plastic. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

NOTE: Wipe syringe GENTLY, vigorous wiping may destroy the label.

LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO Standard used: "YYYY/MM"
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number on the syringe
 - Expiration date

Manufactured by
DENTSPLY
38 West Clarke Avenue
Milford, DE 19963 USA
www.dentsply.com
Tel.: 1-302-422-4511

DENTSPLY
De-Trey-Str.1
78467 Konstanz
Germany
www.dentsply.de
Tel.: 49-7531-583-0

Distributed by
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada

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