For better dentistry

---

**DENTSPLY**

# Dyract®

**THE Evidence-Based Restorative**

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**Caution:** For dental use only

**USA:** Rx only

Content	Page
1 <b>Product description</b> .....	2
2 <b>Safety notes</b> .....	3
3 <b>Step-by-step instructions</b> .....	4
4 <b>Hygiene</b> .....	6
5 <b>Lot number and expiration date</b> .....	7

### 1. Product description

**Dyract®** restorative is a light-curing restorative material for all cavity classes in front and posterior teeth.

**Dyract®** restorative combines the fluoride release of glass-ionomer materials with the strength and aesthetics of a light-curing composite.

#### Evidence-Based

With more than 220 million Compules® Tips sold, the Dyract restorative family<sup>1</sup> belongs to the most frequently used restorative materials worldwide. 17 years of experience and clinical success are documented in 45 clinical investigations and more than 500 scientific publications.

#### 1.1 Delivery forms

- Predosed Compules® Tips

#### 1.2 Composition

- Urethane dimethacrylate (UDMA)
- Carboxylic acid modified dimethacrylate (TCB resin)
- Camphorquinone
- Ethyl-4(dimethylamino)benzoate
- Butylated hydroxy toluene (BHT)
- UV stabilizer
- Strontium-alumino-sodium-fluoro-phosphor-silicate glass
- Highly dispersed silicon dioxide
- Strontium fluoride
- Iron oxide pigments and titanium oxide pigments

#### 1.3 Indications

- Direct restorations of all cavity classes in anterior and posterior teeth including small Class I and II restorations of permanent teeth where the restoration has no direct occlusal contact and the width is less than 1/3 of the intercuspal distance.

<sup>1</sup> This includes Dyract®, Dyract® AP, Dyract® Posterior and Dyract® eXtra.

## 1.4 Contraindications

- Use with patients who have a history of severe allergic reaction to methacrylate resins.
- Use in Class I and II cavities the width of which exceeds 1/3 of the intercuspal distance.
- Use as Core build-up.

## 1.5 Compatible adhesives

The material is chemically compatible with methacrylate-based dentin/enamel adhesives including DENTSPLY adhesives designed for use with visible light cured composite restoratives.

Xeno<sup>®</sup>, self-etching adhesives and Prime&Bond<sup>®</sup> NT, etch and rinse adhesive are specially recommended.

## 2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

### Safety alert symbol

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

## 2.1 Warnings

The material contains methacrylates and polymerizable monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton (and alcohol) and wash thoroughly with water and soap. In case of skin sensitisation or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

## 2.2 Precautions

This product is intended to be used only as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Use of Compules<sup>®</sup> Tips with the Compules<sup>®</sup> Tips gun is recommended.
- Compules<sup>®</sup> Tips may be used for direct intraoral application of restorative material into a cavity or for indirect application by first placing the restorative material on a pad. After intraoral use, discard Compules<sup>®</sup> Tips and do not reuse in other patients in order to prevent cross-contamination.

- The material should extrude easily. **DO NOT USE EXCESSIVE FORCE.** Excessive pressure may cause the Compules® Tip to rupture or to eject from the Compules® Tips gun.
- **Interactions:**
  - Do not use eugenol- and hydrogen peroxide-containing materials in conjunction with this product since they may interfere with hardening and cause softening of the product.

### 2.3 Adverse reactions

- **Eye contact:** Irritation and possible corneal damage.
- **Skin contact:** Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- **Contact with Mucous Membranes:** Inflammation. (See Warnings)

### 2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store in sealed packs at temperatures between 10 °C and 24 °C (50 °F and 75 °F).
- Keep out of sunlight and protect from moisture.
- Do not freeze.
- Do not use after expiration date.

Humidity can adversely affect the properties of Compules® Tips. For unsealed Compules® Tips, the following is recommended:

- Store in a dry environment (< 80% relative humidity).
- Use unsealed Compules® Tips within 4 weeks.

## 3. Step-by-step instructions

### 3.1 Shade selection

Shade selection is made whilst the teeth are hydrated and prior to the restorative procedure. Dyract® restorative is available in 10 VITA®<sup>2</sup> shades (A1, A3, A3,5, A4, B1, B3, C2, C3, C4) including one opaque shade (O-B3).

1. Remove any extraneous plaque or surface stain.
2. Use the VITAPAN®<sup>2</sup> classical shade guide. Use the central part of the respective VITA® teeth for shade selection.  
Alternatively use the Dyract® shade guide that is available on request. The Dyract® shade corresponds to and contains samples of original Dyract® restorative. The color coding dot on the shade guide matches the colored cap on the Compules® Tip.

### 3.2 Cavity preparation

1. Clean the tooth surface with a fluoride free prophylaxis paste (e.g. Nupro® Prophylaxis Paste).
2. Prepare cavity (unless not necessary, e.g. cervical lesion).  
Use adequate isolation such as rubber dam.
3. Rinse surface with water spray and carefully remove rinsing water.  
Do not desiccate the tooth structure.

---

<sup>2</sup> VITA® and VITAPAN® are a registered trademarks of VITA Zahnfabrik, H. Rauter GmbH & Co.

### 3.2.1 Matrix placement and wedging (Class II restorations)

1. Place a matrix (e.g. Automatrix® matrix system or Palodent® matrix system) and wedge. Burnishing of the matrix band will improve contact and contour. Pre-wedging or BiTine® ring placement is recommended.

### 3.3 Pulp protection, tooth conditioning/dentin pre-treatment, adhesive application

Refer to adhesive manufacturer's Instructions for Use for pulp protection, tooth conditioning and/or adhesive application. The DENTSPLY Xeno® adhesives are applied in the self-etch technique. XP Bond® is applied in the etch-and-rinse technique. Prime&Bond® NT may in combination with the Dyract® restoratives be applied in both techniques. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of the material.

### 3.4 Application

Apply material in increments immediately after the application of the adhesive. The material is designed to be cured in increments up to a 3 mm depth/thickness.

#### Using Compules® Tips

1. Insert Compules® Tip into the notched opening of the Compules® Tips gun. Be certain that the collar on the Compules® Tip is inserted first.
2. Remove the colored cap from the Compules® Tip. The Compules® Tip may be rotated 360° to gain the proper angle of entrance into the cavity or to the mixing pad.  
**DO NOT USE EXCESSIVE FORCE.**  
Dispense the necessary amount of material onto a mixing pad and protect against light. Apply material into the cavity with a plastic instrument.  
- or -  
Dispense the material directly into the cavity preparation using a slow, steady pressure.
3. To remove the used Compules® Tip, be sure that the Compules® Tips gun plunger is pulled back completely by allowing the handle to open to its widest position. Apply a downward motion to the front end of the Compules® Tips and remove.

#### Excessive force

Injury

1. Apply slow and steady pressure on the Compules® Tips gun.
2. Do not use excessive force. Compules® Tip rupture or ejection from Compules® Tips gun may result.

### 3.5 Light curing

1. Light cure each area of the restoration surface with a suitable curing light designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470 nm. Minimum light output must be at least 500 mW/cm<sup>2</sup>. Refer to curing light manufacturer's Instructions for Use for compatibility and curing recommendations.
2. Cure each layer according to the table below.

#### Insufficient curing

Inadequate polymerization

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output.

## Curing Light Output $\geq 800$ mW/cm<sup>2</sup>

Shades	Increment/layer	Curing time
All (A1, A3, A3,5, A4, B1, B3, C2, C3, C4, O-B3)	2 mm	10 sec
All (A1, A3, A3,5, A4, B1, B3, C2, C3, C4, O-B3)	3 mm	40 sec

## Curing Light Output 500-800 mW/cm<sup>2</sup>

Shades	Increment/layer	Curing time
All (A1, A3, A3,5, A4, B1, B3, C2, C3, C4, O-B3)	2 mm	20 sec

### 3.6 Finishing and polishing

1. Contour the restoration using finishing burs or diamonds.
2. Use Enhance<sup>®</sup> or Enhance<sup>®</sup> Multi finishing instruments and interproximal finishing strips. In patients with adequate oral hygiene, the high final luster of the restoration comes with use.

For contouring, finishing, and/or polishing, follow the manufacturer's Instructions for Use.

## 4. Hygiene

### 4.1 Compules<sup>®</sup> Tips gun

#### 4.1.1 Disassembling and cleaning

1. Partially close the Compules<sup>®</sup> Tips gun and place your thumb on the rear part of the hinge.
2. Push upwards and lift hinge. The Compules<sup>®</sup> Tips gun is separated in two parts and the plunger is exposed.
3. Remove restorative material, if any, with a soft paper tissue and alcohol of 70% vol.
4. Clean Compules<sup>®</sup> Tips gun by scrubbing with water and soap.

#### 4.1.2 Disinfection

1. Disinfect the Compules<sup>®</sup> Tips gun with a water-based hospital-level disinfection solution according to national/local regulations. Agents containing organic solvents, such as alcohol, may tend to dissolve the plastic.
2. Follow the disinfectant manufacturer's Instructions for Use.
3. Rinse and dry.

#### 4.1.3 Sterilization (optional)

1. Compules<sup>®</sup> Tips gun can be sterilized by steam autoclave at 134 °C/2 bar = 273 °F/29 psi, unwrapped; sterilization time 3 minutes 30 seconds.

#### 4.1.4 Reassembling

1. After sterilization insert plunger into the Compules<sup>®</sup> Tips gun cylinder.
2. Press components together and snap hinge mechanism in place.

Do not continue to use damaged or soiled Compules<sup>®</sup> Tips guns.

Compules<sup>®</sup> Tips gun can be reprocessed for up to 100 times.

## 4.2 Compules® Tips – after intraoral use

### Cross-contamination

#### Infection

1. After intraoral use do not reuse Compules® Tips.
2. Dispose the intraorally used and/or contaminated Compules® Tips in accordance with local regulations.

## 4.3 Shade guide holder and individual tabs

1. Disinfect the shade guide holder and individual tabs with a water-based hospital-level disinfection solution according to national/local regulations. Use of chlorine- (bleach), glutaraldehyde-, phenolic-, iodophor- or organic solvent (e.g. alcohol) based or containing products may over time affect the color shade.
2. During disinfection, remove individual tabs from the shade guide holder.
3. Follow the disinfectant manufacturer's Instructions for Use.
4. After disinfection, clean the shade guide holder and individual tabs by scrubbing with water and soap.
5. Rinse and dry.
6. Do not autoclave the shade guide holder and individual tabs.

## 5 Lot number ( ) and expiration date ( )

1. Do not use after expiration date. ISO standard is used: "YYYY/MM".
2. The following numbers should be quoted in all correspondence:
  - Reorder Number
  - Lot number
  - Expiration date

© DENTSPLY DeTREY 2010-06-02

[These Instructions for Use are based on Master Version 05]

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**DENTSPLY DeTREY GmbH**