# Esthet•X® HD

## High Definition Micro Matrix Restorative

DENTSPLY

### DIRECTIONS FOR USE – ENGLISH

For dental use only.

USA: Rx only.

#### 1. PRODUCT DESCRIPTION

Esthet•X® HD High Definition Micro Matrix Restorative is a visible light cured, radiopaque, composite restorative material for anterior and posterior restorations of primary and permanent teeth. It is to be used following the application of a suitable dentin/enamel adhesive and is compatible with all DENTSPLY adhesives designed for use with visible light cured composite restoratives (see complete Directions for Use of selected adhesive). Use of other dentin/enamel adhesive systems with Esthet•X® HD Restorative material is at the discretion and sole responsibility of the dental practitioner.

##### 1.1 Delivery forms

This one-component, visible light cured composite restorative is available in:

- Predosed Compules® Tips
- 3g Easy•Twist syringes

The regular body shades are optimized to match the Vita® system. Additional available body, opacious dentin and translucent enamel shades are available for esthetically demanding applications existing beyond the traditional Vita® System. A custom blended shade guide is available to guide the clinician in shade selection and blending technique.

Esthet•X® HD Restorative material may be used with separately available BisGMA-based tints and opaques at the discretion and responsibility of the dental practitioner for individual custom characterization of the ultimate esthetic direct-placement restoration.

##### Available Shades

| Opacious Dentin | Regular Body   | Translucent Enamel  |
|-----------------|----------------|---------------------|
| White-O = W-O   | White=W        | B5/Dark Yellow      |
| A2-O            | Extra Light=XL | C1                  |
| A4-O            | A1             | C2                  |
| B2-O            | A2             | C3                  |
| C1-O            | A3             | C4                  |
| C4-O            | A3.5           | C5/Extra Grey-Brown |
| D3-O            | A4             | D2                  |
|                 | B1             | D3                  |
|                 | B2             | Universal           |
|                 | B3             |                     |

##### 1.2 Composition

The resin matrix of Esthet•X® HD Restorative material consists of a Bis-GMA adduct, a Bis-EMA adduct, and triethylene glycol dimethacrylate, Camphorquinone (CQ), Photoinitiator, Stabilizer, Pigments.

The filler combination consists of barium fluoroborosilicate glass with a mean particle size below 1µm and nanofiller silica (particle size 0.04µm).

##### 1.3 Indications

1. Esthet•X® HD Restorative material is indicated as a direct restorative for all cavity classes in anterior and posterior teeth.
2. Esthet•X® HD Restorative material is indicated as a direct esthetic veneering restorative material and cosmetic reshaping (e.g., diastema closure, incisal lengthening).
3. Esthet•X® HD Restorative material is indicated for the indirect fabrication of inlays and onlays.

##### 1.4 Contraindications

Esthet•X® HD Restorative material is contraindicated for use with patients who have a history of severe allergic reaction to methacrylate resins.

## 2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions for use.

### 2.1 Warnings

1. Esthet•X® HD Restorative material contains polymerizable methacrylates which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
2. Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
3. Avoid skin contact to prevent irritation and possible allergic response. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitisation or rash, discontinue use and seek medical attention.

4. Avoid contact with oral mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitisation of mucosa persists, seek medical attention.
5. When using a Compules® Tips Gun, dental professionals and patients should wear protective eyewear.

### 2.2 Precautions

1. This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and is the sole responsibility of the practitioner.
2. Contact with saliva and blood during composite placement may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
3. Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
4. The Esthet•X® HD Restorative material is a light-cured material. Proceed immediately once material has been placed on pad or protect from ambient light.
5. Material should extrude easily from the Compules® Tip. Use a gentle, even motion when exerting pressure on the Compules® Tips Gun. DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause the Compules® Tip to eject from the Compules® Tips Gun.
6. Compules® Tips are designed for one time use. Discard after use. Do not reuse.
7. The Esthet•X® HD Restorative material Easy•Twist syringe should be tightly closed immediately after use.
8. Use of Compules® Tips with the DENTSPLY Caulk Compules® Tips Gun is recommended. Use with another manufacturer's dispensing device is at the discretion and sole responsibility of the practitioner. Prior to use, consult the respective manufacturer's instructions for compatibility, use and sterilization information.
9. **Interactions:**  
Eugenol-containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.

### 2.3 Storage

Esthet•X® HD Restorative material Compules® Tips are for one time use and should be discarded after use. Keep out of direct sunlight and stored in a well ventilated place at temperatures between 2°-24°C/35°-75°F. Allow material to reach room temperature prior to use. Excessive humidity can adversely affect the properties of Esthet•X® HD Restorative material. Protect from moisture. Do not freeze. Do not use after expiration date. The TruMatch® Restorative material shade guide should be disinfected by following the instructions in the Cleaning and Disinfection of the TruMatch® Shade Guide section. The Compules® Tips Gun should be sterilized by following the instructions in the Sterilization and Maintenance of Compules® Tips Gun section. Compules® Tips Gun has an expected useful life of one year from date of purchase.

### 2.4 Adverse Reactions

Product may irritate the eyes and skin. **Eye contact:** irritation and possible corneal damage. **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin. **Mucous membranes:** inflammation, edema, sloughing. See Warnings.

### 2.5 Technique Tips

1. Esthet•X® HD Restorative material is a "sculptable," non-sticky, material that resists slumping and allows the operator to create precise anatomical form.
  - When using Esthet•X® HD Restorative material, it is best not to "spread" the material, but to "push" or "sculpt" the material.
  - It is recommended to dispense/extrude material from a Compules® Tip with a slow, steady force.
  - Rapid and/or forceful extrusion of the material from a Compules® Tip may cause material to appear whitish and/or separated. Upon placement and manipulation, material will coalesce without compromise in strength or esthetics.
2. Esthet•X® HD Restorative material was developed to provide the highest surface luster in order to replicate natural enamel. To achieve maximum polish and maintenance it is highly recommended to fully complete the polish of Esthet•X® HD Restorative material by finishing the restoration with Enhance® Finishing System and polishing with the PoGo® One Step Diamond Micro-Polisher System and/or Prisma•Gloss™ Composite Polishing Pastes.

### 2.6 Optimizing Shade Selection

1. The opacity of Esthet•X® HD Restorative material Body Shades has been scientifically designed to allow the tooth structure enamel and dentin to blend together with the composite with a life-like translucent result. Additional Opaque and Enamel Shades allow replacement of each component of lost tooth structure. Layered opacities of material visually blend to form a restoration with natural vitality and translucence.

2. The final shade will ultimately result from an interaction of the thickness of the layer(s) of restorative composite and surrounding tooth structure (enamel and/or dentin). The intensity of the shade is related to the thickness or thinness of opacous, body and/or translucent shade layers of the Esthet•X® HD Restorative material used. Restoration shade approximation to the selected shade tab of the TruMatch® shade guide occurs with restorations allowing adequate thicknesses of each recipe component, e.g., Class IV, labial-lingual Class III, etc. For restorations less than 2mm total thickness, e.g., conservative cosmetic veneering, and/or where significant shade lightening is desired, placement of a liquid opaquer may provide a more desirable foundation for subsequent shade layering.
3. In all cases, the use of a mock-up is an excellent method for shade verification. The technique includes the use of the shades and opacities as prescribed by the selected TruMatch® shade tab in the appropriate thicknesses to approximate that of the final Esthet•X® HD Restorative material restoration. The materials are placed on the clean, hydrated unetched tooth in the area of the cavity or restoration on the tooth. Each layer is applied without etch or bonding agent, then cured thoroughly. Layering and/or altering shades in varying thickness will allow trial customization of the final shade. The dentist, patient and ancillary dental support personnel can co-jointly view the shade selection result. After viewing, the composite material can be removed by the use of a dental explorer or scaler.
4. Be aware of room and ambient light effects on shade selection; incandescent versus fluorescent lighting versus natural daylight (ideal light is northern exposure/outdoors/overcast).
5. Room accessories can influence shade selection as with colored walls/wallpaper reflection. Also, the patient dental napkin should be removed. Observe tooth shades for short periods of time, using a blue background to "neutralize" the effect of extended shade viewing comparisons. Viewing a blue background will have a relaxing effect on the viewing operator's photo/color optic discrimination. Eyes should be rested. It is useful to have ancillary corroboration on shade selection by dental personnel as well as by the dental patient.

### 3. STEP-BY-STEP INSTRUCTIONS

#### 3.1 Direct Restoration

1. Shade Selection: Before selecting the shade, teeth should be clean, hydrated and free of extrinsic material or stain. The TruMatch® shade guide included in the Introductory Kit and the Complete Organizer System has been created to simulate the result of layering specific shades. Select the shade tab that simulates best the final desired result. Following the recipe on the shade tabs will create the desired shade. Mixing layered shades and various thicknesses will allow customized shade adjustability. Please refer to the Optimizing Shade Selection section above. Additionally, use of the Technique Guide will detail the step-by-step procedure for accurate shade reproduction.
2. Cavity Preparation
  - 2.1 Anterior Restorations: Use conservative cavity preparations for all Class III, IV and Class V restorations. Refinement of the cavo-surface margin (beveling) for enhancement of acid-etching, enamel bonding, and finishing procedures is recommended. **Technique Tip:** Esthetic masking of underlying intrinsic staining may require deeper preparation, allowing for adequate restorative material blending.
  - 2.2 Posterior Restorations: Cavity design requirements are essentially a conventional preparation with refinement of the cavo-surface margin for enhancement of acid etching. No residual amalgam or other base material should be left in the internal forms of the preparation which would interfere with light transmission and the hardening of the restorative.
3. Placement of Matrix: The use of a Mylar, Bimatrix, sectional (e.g. Palodent® Sectional Matrix System) or thin matrix band (e.g. AutoMatrix® Retainerless Matrix System) and subsequent burnishing of the matrix band will improve final interproximal contact and contour. PRE-WEDGING/Bitine® RING PLACEMENT IS ADVOCATED TO ACHIEVE SLIGHT SEPARATION AND FACILITATE ACCEPTABLE PROXIMAL CONTACT.
4. Pulp protection, Tooth Conditioning/Dentin Pretreatment, Adhesive Application: Refer to adhesive manufacturer's directions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of Esthet•X® HD Restorative material.
5. Placement of Esthet•X® HD Micro Matrix Restorative Material
  - 5.1 Easy•Twist Syringe: Remove the cap from the Easy•Twist syringe by using lateral force. It is easier to remove the cap in a "snap-off" rather than a "pull-off" motion. Dispense the necessary amount of Esthet•X® HD Restorative material from the Easy•Twist syringe onto the mix pad by turning the handle slowly in a clockwise direction. To prevent oozing of the material when dispensing is completed, point the front tip of the Easy•Twist syringe upwards and turn the handle counter-clockwise and recap the Easy•Twist syringe. **Technique Tip:** Protect dispensed Esthet•X® HD Restorative material from premature polymerization via ambient light while current increment is adapted and light cured.
  - 5.2 Pre-dosed Compules® Tips: Load Compules® Tips Gun with predosed Compules® Tip. Insert a Compules® Tip into the notched opening of the Compules® Tips Gun barrel. Be certain that the collar on the Compules® Tip is inserted first. Remove the colored cap from the Compules® Tip. The Compules® Tip may be rotated 360° to gain the proper angle of entrance into the cavity. To dispense the material into a cavity preparation, use a slow, steady pressure. Do not use excessive force. To remove the used Compules® Tip, be sure that the Compules® Tips Gun plunger is pulled back completely by allowing the handle to open to its widest position. Apply a downward motion to the front end of the Compules® Tip and remove.
  - 5.3 Optional Flowable Liner: Placement of a compatible flowable liner such as Esthet•X® Flow Flowable Restorative or Dyract•Flow Flowable Compomer (available separately) prior to placement of Esthet•X® HD Restorative material is optional. Follow manufacturer's directions for use.

- 5.4 Anterior Placement: Dispense Esthet•X® HD Restorative material directly into the cavity preparation/tooth surface from the Compules® Tip using slow, steady pressure. Excessive force is not necessary. Alternatively, material may be expressed onto a clean pad from the Compules® Tip or Easy•Twist syringe and carried to the preparation with suitable placement instrument. Adapt, contour and shape with appropriate composite instruments. Material may be placed and light cured in increments up to 2mm (see Curing, Step 6). **Technique Tip:** For esthetic blending of shades, shade layers may be individually cured to form a foundation for subsequent shade placement.
- 5.5 Posterior Placement: Dispense Esthet•X® HD Restorative material directly into the cavity preparation from the Compules® Tip using slow, steady pressure. Excessive force is not necessary. Alternatively, material may be expressed onto a clean pad from the Compules® Tip or Easy•Twist syringe and carried to the preparation with suitable placement instrument.
  - 5.5.1 Class I and V: Esthet•X® HD Restorative material may be placed and adapted in 2mm increments, light curing each increment. (See Curing, Step 6)
  - 5.5.2 Class II: Beginning in the proximal box, firmly adapt a 2mm increment of Esthet•X® HD Restorative material to cavity walls and stabilize matrix band by stepping a clean non-serrated condenser (plugger) from middle to edges, reducing incorporated air and eliminating voids. Light cure (see Curing, Step 6). The remainder of the preparation may be filled in up to 2mm increments, each followed by light cure (see Curing, Step 6). Prior to light curing, contour and shape the final increment with the operator's choice of clean carving and burnishing instruments. **Technique Tips:** To minimize finishing time, contour marginal ridge first before occlusal anatomy. Carve margins and anatomy to final form. Esthet•X® HD Restorative material resists slumping, allowing carving of the majority of anatomical form prior to visible light curing.
6. Curing: Light cure each area of the restoration surface with a suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470nm. Minimum light output must be at least 550mW/cm² exposure for at least 20 seconds. Some advanced performance curing units (e.g. SmartLite® iQ™ L.E.D. Curing Unit) have been shown to cure 2mm increments of most shades of Esthet•X® HD Restorative material in 10 seconds. Refer to curing light manufacturer's recommendations for compatibility and curing recommendations. The Esthet•X® HD Restorative material should be additionally exposed to the curing unit through the proximal, lingual, and buccal enamel walls following matrix removal for the recommended time.
7. Finishing and Polishing
  - 7.1 Begin finishing immediately after curing. Gross excess may be removed and general outline form established with Prisma® Finishing Burs or other carbide or diamond finishing instruments. Additional finishing is recommended by the use of Enhance® Finishing System. See manufacturer's complete directions for use.
  - 7.2 To achieve a very high luster on Esthet•X® HD Restorative material, it is necessary to complete the polishing. PoGo® One Step Diamond Micro-Polisher System and/or Prisma®•Gloss™ Composite Polishing Pastes are recommended. See manufacturer's complete directions for use. Once polished with PoGo® polisher and/or both Prisma®•Gloss™ pastes, Esthet•X® HD Restorative material will maintain its high surface gloss due to its micro matrix technology.

#### 3.2 Fabrication of Indirect Inlay/Onlay

1. Shade Selection: Final desired shade selection should be accomplished prior to tooth preparation. Before selecting the shade, teeth should be clean, hydrated and free of extrinsic material or stain. The TruMatch® shade guide included in the Introductory Kit and the Complete Organizer System has been created to simulate the result of layering specific shades. Select the shade tab that simulates best the final desired result. Following the recipe on the shade tabs will create the desired shade. Mixing layered shades and various thicknesses will allow customized shade adjustability. Please refer to the Optimizing Shade Selection section.
2. Cavity Preparation: Cavity design requirements are essentially a conventional preparation. Rounding of internal angles and refinement of the cavo-surface margin for enhancement of enamel bonding and finishing procedures is recommended. No residual amalgam or other base material should be left in the internal forms of the preparation which would interfere with light transmission and the hardening of the luting cement. Refer to adhesive and/or luting cement manufacturer's directions for base/liner/pulp protection requirements.
3. Impression, Master Cast Fabrication
  - 3.1 Make an accurate impression of the preparation following manufacturer's directions for impression material chosen. A material that allows two pours is recommended. Alternatively, two impressions may be made.
  - 3.2 Prepare and place provisional restoration on the prepared tooth. Use non-eugenol containing material(s).
  - 3.3 Two casts will be fabricated. One working cast will be indexed and separated, per usual crown and bridge procedures. The second master cast will be used for final contour adjustment. Following impression material manufacturer's instructions for disinfection and impression casting, use of die stone is recommended for cast fabrication. Plating of the impression is not recommended. Additionally, an opposing cast should be prepared and articulated.
4. Restoration Fabrication
  - 4.1 Block out any preparation undercuts present. Apply separating medium and/or die spacer to separated die. In most cases, the restoration may be fabricated in no more than 3 increments, each up to 4mm depth.
  - 4.2 Apply first increment/shade, creating restoration body. Adapt to die, staying short of margins. Recommended light curing is accomplished by placing in Triad® Light Curing Unit (DENTSPLY Trubyte) for 2 minutes. Apply second, body layer, allowing cut-back for occlusal anatomy and final proximal contour. Repeat Triad® Light Curing. Prior to

final "enamel" shade placement, custom staining with appropriate compatible stains may be accomplished per manufacturer's instructions. Apply final layer, slightly overfilling and covering all margins. Replace die into articulated model. Establish all external, proximal and occlusal contacts and anatomy. Slight lubrication of adjacent and opposing model teeth is recommended. Cure briefly (10 seconds) with hand-held visible curing light unit to fix contours. Remove die with restoration, place in Triad Unit for final 2 minute curing.

- 4.3 When removing restoration from die, it may be necessary to scrape die stone away from restoration margins, to prevent accidental chipping of restoration. Clean any residual die stone from restoration. Carefully trim away visible flash beyond preparation margins with acrylic bur.
- 4.4 Gently sandblast the internal surfaces with 50m alumina abrasive. Remove visibly undercut material.
- 4.5 Seat restoration onto preparation of uncut master cast, making adjustments as needed. Check for marginal integrity and overall fit and contour. Additional increments may be added if needed by roughening surface and applying compatible adhesive per manufacturer's instructions, followed by placement and light curing as outlined above.
5. Finishing and Polishing (laboratory): Complete as outlined in the above Finishing and Polishing Section, Step 7, for Direct Restorations.
6. Cementation: Recommended cementation technique is bonding using a dual cure adhesive and esthetic resin cement such as Calibra® Esthetic Resin Cement and Prime&Bond® NT™ Dual Cure Universal Dental Adhesive System. Follow adhesive and cement manufacturer's instructions for prepared tooth and restoration surface pretreatments.
7. Adjusting, Finishing and Polishing (clinically): Following cementation, make all necessary occlusal adjustments and polish any clinically adjusted surfaces as outlined in the above Finishing and Polishing Section, Step 7, for Direct Restorations.

## 4. HYGIENE

### 4.1 Cleaning and Disinfection

To prevent Esthet•X® HD Restorative material Easy•Twist syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid package contamination. Repeated disinfection may damage label. Do not attempt to clean, disinfect or re-use Compules® Tips. Properly dispose of used Compules® Tips.

### 4.2 Cleaning and Disinfection of TruMatch® Shade Guide

The TruMatch® shade guide holder and individual tabs may be cleaned by scrubbing with hot water and soap or detergent.

1. The shade guide holder and individual tabs if exposed spatter or spray of body fluids or that may have been touched by contaminated hands, or oral tissues, should be disinfected with a hospital-level disinfectant. Acceptable disinfectants are EPA-registered as tuberculocidal. Disinfect the shade guide holder and individual tabs by spraying the shade guide holder and individual tabs with any recommended hospital-level dual or synergized quaternary ammoniums or immersing in hospital-level disinfectant such as dual or synergized quaternaries or glutaraldehyde for the contact time recommended by the disinfectant manufacturer for optimum results. Spraying with glutaraldehyde is not recommended. Chlorine based (bleach), some phenolic-based agents and iodophor-based products should not be used as they may cause surface staining and/or color alteration. Agents containing organic solvents, such as alcohol, should be avoided, as they may tend to dissolve the plastic. Water-based disinfectant solutions are preferred. The disinfectant manufacturer's directions should be followed properly for optimum results.
2. Following disinfection, thoroughly rinse and dry the TruMatch® shade guide holder and individual tabs before storage. Autoclaving the TruMatch® shade guide holder and individual tabs is not recommended.

### 4.3 Sterilization and Maintenance of Compules® Tips Gun

1. To clean the Compules® Tips Gun, the following procedure is recommended. Partially close the syringe and place thumb under the rear portion of the hinge.
2. Push upward and lift hinge separating the Compules® Tips Gun, exposing the plunger. If the plunger has excess material build-up clean with an alcohol moistened gauze. The Compules® Tips Gun may be cleaned by scrubbing with hot water and soap or detergent. It is not recommended to submerge Compules® Tips Gun into disinfection solutions.
3. The Compules® Tips Gun may be autoclaved. **NOTE:** As with any plastic instrument, the Compules® Tips Gun may weaken over time.
4. To reassemble, insert plunger into Compules® Tips Gun barrel, press components together and snap hinge mechanism in place. Prior to each use check to make sure that the Compules® Tips Gun is fully engaged and in good working order.

### 5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard uses: "YYYY/MM"
2. Compules® Tips Gun is warranted for one year from date of purchase.
3. The following numbers should be quoted in all correspondence:
  - Reorder Number
  - Lot number on Compules® Tip/Easy•Twist syringe
  - Expiration date on Compules® Tip/Easy•Twist syringe

Vita is not a registered trademark of DENTSPLY International.  
©2008 DENTSPLY International. All Rights Reserved. Printed in U.S.A.

Manufactured by:  
DENTSPLY Caulk  
38 West Clarke Avenue  
Milford, DE 19963 USA  
Tel.: 1-302-422-4511  
www.dentsply.com

Distributed by:  
DENTSPLY Canada  
161 Vinyl Court  
Woodbridge, Ontario  
L4L 4A3 Canada

DENTSPLY DeTrey GmbH  
De-Trey-Str. 1  
78467 Konstanz  
Germany  
Tel.: 49-7531-583-0  
www.dentsply.de

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## Complete System

# Enhance® Finishing & PoGo® Polishing Systems

DENTSPLY  
CAULK

The Enhance® Finishers (Discs, Cups and Points) are pre-mounted single use aluminum oxide impregnated cured urethane dimethacrylate resin finishers designed for preparing composite surfaces for their final polish. PoGo® Polishers (Discs, Cups and Points) are pre-mounted, single-use diamond impregnated cured urethane dimethacrylate resin polishing devices designed for use in the final polishing of all accessible composite resin restorations. When used as directed, the Enhance® Finishing and PoGo® Polishing System will yield a high surface luster on hybrid composites such as Esthet•X® Micro Matrix, TPH® Spectrum®, Prisma® TPH®, and Prisma® AP.H® Composites, SureFil® High Density Restorative Material and microfilled composite restoratives. The Enhance® Finishing & PoGo® Polishing Systems are also effective in final polishing of compomer and flowable composite such as Esthet•X® flow and compomer restoratives such as Dyract® AP and Dyract®flow.

**Cautions:** U.S. Federal law restricts this device to sale by or on the order of a dentist.

### INDICATIONS

1. Enhance® Finishing Discs, Cups and Points are indicated for contouring and finalizing the form of composite and compomer restorations.
2. PoGo® Polishing Discs, Cups and Points are indicated as the final step in polishing all properly contoured and finished composite resin and compomer restorations.

### COMPOSITION

Enhance®:

Polymerized Urethane Dimethacrylate Resin; Aluminum Oxide; Silicon Oxide; Plastic latch-type mandrel

PoGo®:

Polymerized Urethane Dimethacrylate Resin; Fine Diamond Powder; Silicon Oxide; Plastic Latch-type mandrel

### CONTRAINDICATIONS

Do not reuse. Enhance® Finishers & PoGo® Polishers are indicated as a single use only, disposable item. Attempts to sterilize and/or reuse these items can adversely affect their performance. Enhance® FINISHERS & PoGo® POLISHER SHOULD NOT BE USED IN HIGH-SPEED HANDPIECES.

### WARNINGS

1. Enhance® Finishing Discs, Cups, Points and PoGo® Polishing Disc, Cups and Points are designed for use without water. As with any rotary instrument, heat will build up with prolonged contact. Use with light intermittent pressure. Application of an air cooling stream during finishing and polishing assists with residue removal.
2. Use of the Enhance® Finishers and PoGo® Polishers with excessive lateral pressure may dislodge finisher and/or polisher from the mandrel or cause quick deterioration and decreased effectiveness. This type of excess pressure may cause gouging of the restoration surfaces.
3. Care should be taken to have good control on any rotary instrument to protect patients from injury. Contact with soft tissues (gingiva, mucosa) may result in injury to the tissue. If contact occurs, thoroughly wash the affected area with water and seek appropriate medical attention.

### PRECAUTIONS

1. This product is intended to be used as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Use in a well ventilated area. Avoid inhaling airborne particulates.
3. Protect adjacent restorations from accidental contact.
4. Storage & Shelf Life: "Dry Storage is Recommended". Store the Enhance® Finishing and the PoGo® Polishing System between 65°F-80°F (18°C-25°C) and in relative humidity of 50% ± 20%. Avoid direct sunlight. Do not use after expiry date.

### ADVERSE REACTIONS

Allergic contact dermatitis and other allergic reactions may occur in susceptible individuals.

### STEP BY STEP INSTRUCTION FOR USE

STEP 1: Place restoration or composite veneer as per standard techniques. Complete gross reduction of excess and general outline form of the restoration using carbide finishing burs (i.e. Prisma® Finishing Burs), diamond finishing burs, etc.

STEP 2: Insert an Enhance® Finisher (Disc, Cup or Point) into a conventional speed contra-angled handpiece and continue finishing. NOTE: ENHANCE® DISCS, CUPS AND POINTS MUST NOT BE USED IN HIGH SPEED HANDPIECES. The aggressiveness of the Enhance® Finishers is controlled by the pressure applied to the surface of the composite. The greater the pressure, the more material is removed; lighter pressure leaves a smooth surface without removing bulk.

#### Technique Tip:

1. Prior to proceeding to the PoGo® Polishing System restoration surface should be finished to a final contour. Surface should be smooth and defect free.
2. Enhance® Finishing Points and PoGo® Polishing Points are ideal for posterior occlusal and concave lingual surfaces.
3. Enhance® Finishing Cups and PoGo® Polishing Cups are well suited for cervical and proximal line angle areas.

STEP 3: Insert a PoGo® Polishing Disc, Cup or Point into a conventional speed contra-angled handpiece. NOTE: PoGo® POLISHERS MUST NOT BE USED IN HIGH SPEED HANDPIECES. Begin polishing by applying light intermittent pressure at moderate speed to the restoration surface. When using the disc, the flat broad surface of the PoGo® Polisher should be used in contact with the restoration to optimize its effectiveness. To increase surface luster, decrease the pressure applied to the restoration surface using a light buffing motion

**Technique Tip:** The flexibility of PoGo®, when used lightly, allows polishing of established line angles without alteration of form.

### STERILIZATION AND MAINTENANCE

Do not attempt to sterilize and/or reuse the Enhance® Finishers and/or PoGo® Polishers (disc, cup or point.) Properly dispose after use.

### LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard used: "YYYY/MM"
2. The following numbers should be quoted in all correspondences.
  - Reorder number
  - Lot number on package
  - Expiration date on package

Manufactured by:  
DENTSPLY Caulk  
38 West Clarke Avenue  
Milford, DE 19963 USA  
Tel.: 1-302-422-4511  
www.dentsply.com

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## Light Cured Dental Adhesive

### DIRECTIONS FOR USE - ENGLISH

Prime&Bond® NT™ Adhesive is a light cure self-priming dental adhesive designed to bond light cured composite and compomer materials to enamel and dentin as well as to metals and ceramic. Prime&Bond® NT™ Adhesive combines primer and adhesive in a single container. The reduction of components and treatment steps simplifies use, maintaining superior bond strengths and protection against microleakage. Prime&Bond® NT™ Adhesive is available in a 4.5 ml tri-laminated bottle and in pre-dosed, 0.125 ml disposable containers for single patient use. For routine dual cure procedures, please see the Prime&Bond® NT™ Dual Cure Universal Dental Adhesive System instructions that include the Prime&Bond® NT™ Adhesive and Self Cure Activator components.

Caulk® 34% Tooth Conditioner Gel is included in both the Prime&Bond® NT™ Adhesive Economy Kit.

**Caution:** U.S. federal law restricts this device to sale by or on the order of a dentist.

### COMPOSITION

#### Prime&Bond® NT™ Adhesive:

Di- and Trimethacrylate resins; PENTA (dipentaerythritol penta acrylate monophosphate); Photoinitiators; Stabilizers; Nanofillers - Amorphous Silicone Dioxide Cetylamine hydrofluoride; Acetone

**Caulk® 34% Tooth Conditioner Gel** (See complete directions for use enclosed):  
Water; 34% Phosphoric acid; Silicon Dioxide; Surfactants; Blue Colorant

### INDICATIONS

1. Direct, light cured composite and compomer restorations.
2. Indirect Restorations; light cured, resin cemented veneers.
3. Composite, ceramic and amalgam repairs.
4. Cavity varnish for use with fresh amalgam.

### CONTRAINDICATIONS

1. Prime&Bond® NT™ Adhesive is contraindicated for use with patients who have a history of severe allergic reaction to methacrylate resins or any of the components
2. Prime&Bond® NT™ Adhesive is contraindicated for direct application to dental pulp tissue (direct pulp capping).
3. Caulk® 34% Tooth Conditioner is contraindicated for use on dental pulp tissue.

### WARNINGS

1. Caulk® 34% Tooth Conditioning Gel contains phosphoric acid, which may cause burns. Avoid contact with oral tissues, eyes and skin. Do not take internally.  
**Eye Contact:** Tooth Conditioning Gel contains phosphoric acid, which may cause burns or may be irritating to eyes. Before using this product wear protective glasses as well as covering the patient's eyes to protect from squirting the material inadvertently. In case of accidental contact with eyes, rinse eyes immediately with plenty of water and seek medical attention.  
**Skin Contact:** Tooth Conditioning Gel contains phosphoric acid, which can cause burns or skin sensitization in susceptible individuals. If contact with skin occurs immediately wipe off and flush with generous amounts of water, then wash well with soap and water after contact. If skin rash, irritation, sensitization or other allergic reaction occurs, seek medical attention immediately.  
**Oral Mucosa Contact:** Avoid contact with oral soft tissues. If accidental contact occurs, flush mucosa with plenty of water and expectorate water. If sensitization of mucosa persists, seek medical attention immediately.  
**Ingestion:** Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. If nausea or illness develop, seek medical attention immediately. Contact regional Poison Control Center immediately.
2. Prime&Bond® NT™ Adhesive contains polymerizable methacrylate monomers. Avoid prolonged or repeated contact with skin (allergic contact dermatitis), oral soft tissues, and eyes. Avoid prolonged inhalation. Do not take internally.  
**Eye contact:** Prime&Bond® NT™ Adhesive contains methacrylates which may be irritating to eyes. Before using this product wear protective glasses as well as covering the patient's eyes to protect from splashing material. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.  
**Skin contact:** Prime&Bond® NT™ Adhesive contains polymerizable monomers which can cause skin sensitization (allergic contact dermatitis) in susceptible individuals. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water after contact. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention.  
**Oral mucosa contact:** Avoid contact with oral soft tissues. If accidental contact occurs,

flush mucosa with plenty of water and expectorate water. If sensitization of mucosa persists, seek medical attention immediately.

**Ingestion:** Do not swallow or take internally. If accidental swallowing occurs, drink lots of water. If nausea or illness develop, seek medical attention immediately. Contact regional Poison Control Center if necessary.

3. Prime&Bond® NT™ Adhesive contains acetone. Do not breathe vapor. **Inhalation:** Give oxygen or artificial respiration if necessary.
4. Caulk® 34% Tooth conditioning Gel in syringes should extrude easily: DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material.
5. Caulk® 34% Tooth Conditioning Gel should not be used with patients who have a history of severe allergic reaction to any of the components.

### PRECAUTIONS

1. This product is intended to be used only as specifically outlined in the *Directions for Use*. Any use of this product inconsistent with the *Directions for Use* is at the discretion and sole responsibility of the practitioner.
2. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
3. Eugenol-containing dental materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
4. Insufficient data exist to support the use of desensitizing agents and/or cavity cleansing agents with Prime&Bond® NT™ Adhesive. Recommended pre-treatment is pumice or prophylaxis paste with a rubber cup. (See *Step-by-Step Instructions*)
5. Contact with saliva, blood and/or some astringent solutions during adhesive procedures may cause failure of the restoration. Use of a rubber dam or adequate isolation is recommended.
6. Exercise caution if mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures. In vitro studies have suggested compromised adhesive bond strength to enamel and dentin contaminated with hemostatic agents. Marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.
7. In vitro studies have shown a relationship between adhesive performance and air-thinning/drying/solvent evaporation techniques. Improper over-thinning or under-evaporating may compromise bond strength, leading to microleakage, post-operative sensitivity and/or restoration failure. (See *Step-by-Step Instructions*)
8. Replace original cap of Caulk® 34% Tooth Conditioning Gel tightly after each use to avoid evaporation. Discard dispensing tip after use, as dispensing tips may clog if gel is allowed to dry inside.
9. DENTSPLY Caulk supplies the appropriate dispensing tip for the Tooth Conditioning Gel. This tip is the only tip that should be utilized for placement of the gel.
10. The Prime&Bond® NT™ Adhesive bottle should be tightly closed immediately after use.
11. Prime&Bond® NT™ Adhesive is light-cured material. Proceed immediately once material has been placed in mixing well or protect from ambient light. The components are cured by visible light.
12. Use only in well ventilated areas.
13. **Flammable:** Prime&Bond® NT™ Adhesive contains acetone. Keep away from sources of ignition.
14. Single dose containers are intended for single use only, and should be discarded after use. Do not re-seal or re-use.
15. The efficacy of Caulk® 34% Tooth Conditioning Gel as an etchant/cleanser of indirect restorations has not been demonstrated.
16. Variable in-vitro data exist regarding use of light-cured-only adhesives such as Prime&Bond® NT™ Adhesive without Self Cure Activator in conjunction with dual-cured or self-cured resins such as Calibra® Esthetic Resin Cement and FluoroCore®2 Fluoride-Releasing Dual Cure Core Build-Up Material in limited or no light curing applications. Chemical/Product incompatibility may adversely affect product efficacy, leading to premature restoration failure.
17. **Storage and Shelf Life:** Prime&Bond® NT™ Adhesive should be kept out of direct sunlight and stored in a well ventilated place at room temperature not exceeding 25°C/77°F. Refrigerated storage is not required, but is acceptable when not in use. Allow material to reach room temperature prior to use.  
Caulk® 34% Tooth Conditioning Gel: Store syringe with original cap only. Do not store syringe with dispensing tip in place, as this may harden the material making extrusion difficult (see Warnings). Not to be stored at temperature exceeding 25°C/77°F. Protect from moisture. Do not freeze. Do not use after expiration date.

## ADVERSE REACTIONS

1. Product may irritate the eyes and skin. **Eye contact:** irritation and possible corneal damage. **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin. **Mucous membranes:** inflammation, edema, sloughing. (See *Warnings*)
2. Product may cause serious health effects if ingested. (See *Warnings*)
3. Caulk® 34% Tooth Conditioning Gel may cause pulpal effects. (See *Contraindications and Warnings*)
4. Inhalation of vapors may cause varying degrees of damage to the affected tissue and also increased susceptibility to respiratory illness. (See *Precautions*)
5. The following medical conditions are generally aggravated by exposure: Individuals with pre-existing skin disorders, eye problems or impaired liver, kidney and respiratory functions.

## STEP-BY-STEP INSTRUCTIONS FOR USE (LIGHT CURED)

1. **DIRECT RESTORATION (Light Cured Composite Resin and Compomers)**
  - 1.1 **Cleaning:** Clean uninstrumented enamel and dentin with a rubber cup and pumice or a cleaning paste such as Nupro® Prophylaxis Paste. Wash thoroughly with water spray and air dry. Clean freshly instrumented enamel and dentin with water spray and then air dry.
  - 1.2 **Pulp Protection:** For direct and indirect pulp capping, cover the dentin close to the pulp (less than 1mm) with a hard setting calcium hydroxide liner (Dycal® Liner) leaving the rest of the cavity surface free for bonding with Prime&Bond® NT™ Light Cured Dental Adhesive.
  - 1.3 **Tooth Conditioning/Dentin Pretreatment:** When used as a bonding agent for composite materials, it is recommended to follow the total etch technique described in 1.3.1: When used as a bonding agent for compomer materials in non-stress bearing situations, acid etching is optional.
    - 1.3.1 **Application of Caulk® 34% Tooth Conditioner Gel:** After application of rubber dam or other suitable isolation technique, apply Caulk® 34% Tooth Conditioner Gel. Attach disposable needle to end of syringe. Needle tip may be bent for easy access. Gently extrude Caulk® 34% Tooth Conditioner Gel (34% phosphoric acid) to the cavity surfaces starting at the enamel margins. For best results, condition enamel for at least 15 seconds and dentin for 15 seconds or less. Alternatively the conventional Enamel Etch Technique can be followed. In this case, the enamel margins only are treated with Caulk® 34% Tooth Conditioner Gel for at least 15 seconds.
    - 1.3.2 **Rinsing and Blot Drying:** Remove gel with aspirator tube and/or vigorous water spray and rinse conditioned areas thoroughly for at least 10 seconds. Blot dry conditioned areas with a moist cotton pellet. Proceed immediately to application of Prime&Bond® NT™ Light Cured Dental Adhesive. **NOTE:** The goal of blot drying is to provide the correct amount of "wetness" on the tooth surface by removing all excess moisture and yet avoiding desiccating the dentin. Do not rub the tooth surface when blot drying. Dentin should be blotted until there is no pooling of water, leaving a moist, glistening surface. Once the surfaces have been properly treated, **they must be kept uncontaminated.** If salivary contamination occurs, repeat procedure beginning at Step 1.1.
  - 1.4 **Application of Prime&Bond® NT™ Light Cured Dental Adhesive**
    - 1.4.1 **Unit Dose Container:** Grasp container at each end, or insert into holder, placing thumb along center score. Firmly apply pressure until container separates. Holder may be placed on the tabletop or held between fingers for convenience. Insert disposable applicator into opening to saturate applicator tip. **Conventional Bottle:** Dispense Prime&Bond® NT™ Adhesive directly onto a clean disposable brush provided, making sure that the bottle does not come in direct contact with the brush, or place 2-3 drops of Prime&Bond® NT™ Adhesive into a clean well. Replace cap promptly.
    - 1.4.2 Using the disposable brush or applicator tip supplied, immediately apply generous amounts of Prime&Bond® NT™ Adhesive to thoroughly wet all the tooth surfaces. This surface should remain fully wet for 20 seconds and may necessitate additional applications of Prime&Bond® NT™ Adhesive.
    - 1.4.3 Remove excess solvent by gently drying with clean, dry air from a dental syringe for at least 5 seconds. Recommended technique is to begin drying/evaporation with the syringe approximately 4-6 inches (10-15cm) away from the surface, gradually bringing the source to within 0.5 inches (10mm) of the surface over 5-10 seconds. Surface should have a uniform glossy appearance. If not, repeat application and air dry. Surface should not show areas of excessive adhesive thickness or pooling. Repeat air drying/evaporation outlined above if necessary.
      - 1.4.4 Cure Prime&Bond® NT™ Adhesive for 10 seconds<sup>1</sup> using a curing light.
    - 1.5 **Completion:** Place restorative material over the cured Prime&Bond® NT™ Adhesive as per restorative material manufacturer's directions.
2. **VENEER CEMENTATION (Light Cured Resin Cement)**
  - 2.1 **Cleaning:** See Section 1.1
  - 2.2 **Tooth Conditioning/Dentin pretreatment:** See Section 1.3
  - 2.3 **Application of Prime&Bond® NT™ Light Cured Dental Adhesive:** Apply and light cure Prime&Bond® NT™ Adhesive as described for direct restorations, see Section 1.4.
  - 2.4 **Preparation of restoration**
    - 2.4.1 Treat bonding surface of restoration according to manufacturer's or dental laboratory's instructions, i.e. etching or mechanical roughening.
    - 2.4.2 Apply Calibra® Silane Coupling Agent (available separately) to ceramic or porcelain restoration following manufacturer's instructions.
    - 2.4.3 Apply a single coat of Prime&Bond® NT™ Adhesive to the internal bonding surface of the restoration. Immediately air dry for 5 seconds.
    - 2.4.4 Cure Prime&Bond® NT™ Adhesive for 10 seconds<sup>1</sup> using a curing light.

- 2.5 **Cementation:** Prepare and apply light-cured resin cement according to manufacturer's instructions.

3. **COMPOSITE, CERAMIC AND AMALGAM REPAIRS (Light Cured Composite Resin and Compomers)**
  - 3.1 **Preparation, Cleaning:** Roughen the fractured surface as much as possible with a diamond bur. Create where possible. For best results, sandblast the area to be repaired with an intraoral microetcher (50µ alumina). Rubber dam is recommended with high speed evacuation. Rinse microetched areas for 15-20 seconds with water. Air dry. **NOTE:** Bevel porcelain margins before etching.
  - 3.2 **Tooth Conditioning/Dentin Pretreatment:** Etch tooth with Caulk® 34% Tooth Conditioner Gel as outlined in Section 1.3. Etch ceramic restoration repair area with hydrofluoric acid according to manufacturer's instructions.
  - 3.3 **Treatment of the Restoration:** Rinse with water for 10 seconds. Air dry. **NOTE:** Apply Calibra® Silane Coupling Agent to porcelain surfaces to be repaired following manufacturer's instructions.
  - 3.4 **Application of Prime&Bond® NT™ Light Cured Dental Adhesive:** Apply and light cure Prime&Bond® NT™ Adhesive as described for direct restorations, see Section 1.4.
  - 3.5 **Completion:** Complete repair with placement and curing of desired shade(s) of light cured composite restorative following manufacturer's directions.
4. **CAVITY VARNISH FOR USE WITH FRESH AMALGAM**

When used as a cavity varnish, Prime&Bond® NT™ Light Cured Adhesive is not an amalgam adhesive.

  - 4.1 **Preparation:** Finish preparation by removing existing restorations and/or caries.
  - 4.2 **Cleaning:** Clean preparation and place Dycal® Liner if needed, following Sections 1.1 and 1.2.
  - 4.3 **Tooth Conditioning/Dentin Pretreatment:** Rinse and carefully air-dry cavity preparation, but do not desiccate exposed dentin. **NOTE:** Acid etching of prepared cavity (enamel or enamel/dentin) is optional prior to placement of Prime&Bond® NT™ Adhesive. If desired, follow Section 1.3.
  - 4.4 **Application of Prime&Bond® NT™ Light Cured Dental Adhesive:** Apply and light cure Prime&Bond® NT™ Adhesive as described for direct restorations, see Section 1.4.
  - 4.5 **Completion:** Place and condense amalgam (e.g. Dispersalloy® Dispersed Phase Alloy) as per manufacturer's directions.

## CLEANING AND DISINFECTION

To prevent Prime&Bond® NT™ Adhesive bottle from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid package contamination. Repeated disinfection may damage label.

Do not attempt to clean, disinfect or re-use applicator brush or unit dose container. Properly dispose used brushes and containers. The re-usable brush handle may be cleaned by scrubbing with hot water and soap or detergent. Do not autoclave brush handle or unit dose holder. Disinfect as outlined below.

## DISINFECTION OF RE-USABLE BRUSH HANDLE AND UNIT DOSE HOLDER

The brush handle or holder, if exposed to spatter or spray of body fluids or that may have been touched by contaminated hands, or oral tissues, should be disinfected with a hospital-level disinfectant. Acceptable disinfectants are EPA-registered as tuberculocidal. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternaries are approved disinfectants. Disinfect the brush handle by spraying with or immersing in any recommended hospital-level disinfectant except neutral glutaraldehyde for the contact time recommended by the disinfectant manufacturer for optimum results. Some phenolic-based agents and iodophor-based products may cause surface staining. Agents containing organic solvents, such as alcohol, may tend to dissolve the plastic. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

## LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard is used: "YYYY/MM"
2. The following numbers should be quoted in all correspondence:
  - Reorder Number
  - Lot number on syringe/bottle
  - Expiration date on syringe/bottle

<sup>1</sup> Check curing light for minimum curing output of at least 550 mw/cm2 and a spectral output including 470nm (peak absorption of the CQ photoinitiator).

©2006 DENTSPLY International Inc. All Rights Reserved. Printed in U.S.A.

0120

Manufactured in the U.S.A by:  
DENTSPLY Caulk  
38 West Clarke Avenue  
Milford, DE 19963  
Tel.: 1-302-422-4511

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©2006 DENTSPLY International Inc. Tous droits réservés. Imprimé aux U.S.A.