Jeltrate® Chroma

DENTSPLY
CAULK

Dustless Color-Changing Alginate Impression Material – Fast Set

DIRECTIONS FOR USE – ENGLISH

Caution: For dental use only.
USA: Rx only.

1. PRODUCT DESCRIPTION

Jeltrate® Chroma Dustless Color-Changing Alginate Impression Material is a conventional alginate dental impression material intended to be placed on a preformed impression tray and used to reproduce the structure of a patient's teeth. Jeltrate® Chroma Alginate is also intended to provide models for study and for the production of restorative prosthetic devices, such as dental inlays and dentures.

Jeltrate® Chroma Alginate is formulated to change colors for each preparation phase: fuchsia during mixing, lavender during tray loading and seating to light blue during setting in the mouth. Jeltrate® Chroma Alginate is a dustless formulation, to reduce airborne particles when dispensing and mixing, for a healthier environment and easier cleanup.

1.1 Delivery forms

Jeltrate® Chroma Dustless Color-Changing Alginate Impression Material is available in:

- 454g pouches
Contents: 26 Two-Scoop Average Impressions
17 Three-Scoop Large Impressions

1.2 Composition

Diatomaceous Earth; Sodium Alginate; Sodium Carboxymethylcellulose; Calcium Sulfate; Trisodium phosphate; Magnesium hydroxide; Phenolphthalein; Blue Organic Pigment; Spearmint Oil; Mineral Oil

1.3 Indications for use

Jeltrate® Chroma Dustless Color-Changing Alginate Impression Material is indicated for making dental impressions to provide models for:

- study
- the production of restorative prosthetic devices

1.4 Contraindications

None known.

1.5 Technical specifications

ISO Specification 1563-1990

FAST SET

Mixing Time 0'30"-0'45"

Total working time 1'00"

Initial setting time 1'45"

Setting time*. 2'00"

*Setting time includes mixing time and the minimum time for handling based upon mixing powder with distilled or deionized water at 23°C. The temperature and hardness (tap water) of the mixing water may vary the color change timing and total working and setting time.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. (The following statement is required by the State of California). This product contains phenolphthalein which is known by the State of California to cause cancer. Prior to use, read the Safety Data Sheet (SDS) information and product instructions.
2. Jeltrate® Chroma Alginate contains crystalline silica. Avoid prolonged exposure to eyes. Avoid prolonged inhalation. Do not take internally. Mix water and powder carefully to minimize dust inhalation risk and decrease airborne particles to exposed individuals.
 - **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - **Avoid inhalation** to prevent irritation. In case of inhalation, move subject to fresh air. Drink water to clear throat and blow nose to evacuate material. If irritation or soreness persists, seek medical attention.
 - **Avoid ingestion/swallowing** of material to prevent irritation or obstruction. If accidental swallowing occurs, drink lots of water. This material is not hazardous when small quantities are ingested.
 - **Avoid multiple impressions at one sitting** (more than three of the same arch/area) to prevent irritation to the oral mucosa. If symptoms occur discontinue use immediately. If irritation does not subside consult a physician.
3. Jeltrate® Chroma Alginate contains coloring agents and should not be used with patients who have a known hypersensitivity to them. In the event of allergic reaction, seek immediate medical attention.

2.2 Precautions

- This product is intended to be used as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
- Wear suitable protective eyewear, clothing, mask and gloves. Protective eyewear is recommended for patients.
- Use in a well-ventilated area. Avoid inhaling airborne particulates. (See Warnings)
- Replace original canister lid tightly immediately after dispensing.
- The temperature and hardness of the water used for mixing will affect the color change timing and the work time and set time. Recommended water temperature is 23°C/73°F. (See Step-by-Step instructions.)
- Jeltrate® Chroma Alginate is not recommended for transfer impressions or final impressions for fixed restorations.
- Impressions should be cast immediately (upon completion of disinfection) in a suitable gypsum material. Do not store impressions submerged in water. Impression will remain dimensionally stable for 5 days if wrapped in a moist paper towel and stored in a tightly sealed plastic bag.
- Prior to each use, fluff powder before dispensing by gently tumbling tightly closed canister.

Interactions:

- Do not expose material to, or use in conjunction with hydrofluoric acid. Exposure of the silica contained within Jeltrate® Chroma Alginate to hydrofluoric acid renders the product unstable.

2.3 Adverse reactions

- **Inhalation:** Inhalation of crystalline silica has been classified as carcinogenic. Inhalation of powder may cause varying degrees of damage to the affected tissue and also increased susceptibility to respiratory illness.
- **Eye contact:** irritation and possible corneal damage.
- **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin.
- **Bowel obstruction** or other digestive distress may result from ingestion of mixed impression material. (See Warnings)

2.4 Storage

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store at temperatures between 10°C/50°F – 24°C/75°F.
- Use the product at room temperature.
- Transfer pouch contents to the canister container immediately after opening.
- Store canister with the lid tightly closed.
- Keep out of direct sunlight.
- Do not freeze.
- Protect from moisture.
- Do not use after expiration date.

Store clean, dry, properly disinfected plastic water and powder measures at room temperature. Do not store scoop and/or cup in the canister containing Jeltrate® Chroma Alginate.

3. STEP-BY-STEP INSTRUCTIONS

1. Initial Use (454g pouch). Cut pouch top and transfer contents into a suitable container (clean, dry canister with tight fitting lid is recommended). Keep empty pouch or record Lot number and Expiration date for future reference. Thereafter, follow Steps 2-12.
2. Select and/or prepare a suitable impression tray. Rigid trays (with sides) are recommended for most applications. Trays with retentive features such as Caulk Rim-Lock® Trays or perforations are recommended. Commercially available alginate tray adhesive is not necessary with these trays, but may be used if desired following manufacturer's instructions.
3. **Before Dispensing Material:** Tumble closed container gently to fluff powder. (See Warnings, Precautions).
4. Dip supplied scoop lightly into powder until the scoop is full. Do not pack scoop against container wall.
5. Tap scoop against rim to insure full measure without voids. Scrape excess off with a spatula to achieve level scoop. **Technique tip:** Tap spatula lightly on the filled scoop until the powder is fully condensed. Then scrape excess if necessary.
6. Empty the scoop into a clean, dry mixing bowl. Repeat steps 4 and 5 until the desired amount of powder is in the mixing bowl.
7. For each one (1) scoop (8g) of powder add 1/3 measure of water (19mL). For each two (2) scoops of powder add 2/3 measure of water. For each three (3) scoops of powder add one (1) whole measure of water. For best results, use distilled water at approximately 73°F (23°C). **Technique tip:** Colder water retards setting and elongates color change; warmer water hastens setting and shortens color change.
8. Mix water and powder carefully at first, until all powder is wetted. Once all powder is incorporated, continue to spatulate in a vigorous fashion. Spatulate against the bowl sides to reduce air incorporation. Do not whip or stir material. **NOTE COLOR CHANGE:** Blue powder immediately turns fuchsia upon water addition. Mix until fuchsia color is uniform and streak-free, approximately 0'30".
9. Once the material appears to be creamy but with consistency, apply the Jeltrate® Chroma Alginate onto the tray. **NOTE COLOR CHANGE:** Mixed material turns from fuchsia to lavender to indicate tray loading, approximately 0'45". **Technique Tip:** To minimize air incorporation, wipe the loaded spatula against the tray rim, allowing material to flow from the spatula into the tray.

10. Insert filled tray into the mouth. **NOTE COLOR CHANGE:** Mixed material turns from lavender to light blue to indicate mouth insertion, approximately 0'60". Press tray gently into position. Do not overseat tray. Immediately muscle trim. Hold seated tray immovable for 1 minute after stickiness vanishes. There is no further color change. When mixed with distilled water at 23°C, the set time of Jeltrate® Chroma Alginate is approximately 2'00" from start of mix.

Timing Guide

0'00"	0'45"	1'00"	2'00"
Fuchsia	Fuchsia to Lavender	Lavender to Light Blue	
Begin mixing	Complete mixing, load tray	Seat loaded tray	Remove set impression

11. Once material sets, break seal, snap and remove impression from mouth. Have patient rinse with warm water. Check for and remove residual material from mouth and/or between teeth.
12. Rinse impression thoroughly under running water. Immediately disinfect following disinfection procedures below before pouring. Impression should be poured immediately after disinfection, following casting procedures below. If later pouring is necessary, wrap impression in wet towel and place into a plastic bag, sealing the bag tightly. Impression will remain dimensionally stable for up to 5 days. **DO NOT STORE IMPRESSION SUBMERGED IN WATER.**

4. HYGIENE

Cross-contamination

Reprocess reusable products as described below.

4.1 Disinfection of the impression

The impression should be disinfected with a hospital-level tuberculocidal disinfectant according to national and local regulations.

- Sodium hypochlorite (0.525%), phenols, and dual or synergized quaternary ammoniums are approved disinfectants.
- Thoroughly soak by spraying the impression in any recommended hospital level disinfectant for the contact time recommended by the disinfectant manufacturer for optimum results.
- Water-based disinfectant solutions are preferred.

4.2 Cleaning the water and powder measures, canister

The water (cup) measure, powder (scoop) measure and storage canister may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.

4.3 Disinfecting the water and powder measures, canister

- Disinfect water and powder measures and storage canister with a hospital-level tuberculocidal disinfectant according to national and local regulations. Use impregnated wipes not sprays.
- **NOTE:** Vigorous wiping can destroy the label. Wipe canister gently.
- Some phenolic based agents and iodophor based products may cause surface staining.
- Agents containing organic solvents, such as alcohol may tend to dissolve the plastic. The disinfectant manufacturer's directions should be followed properly for optimum results.
- Following disinfection, thoroughly rinse and dry the cup and/or scoop before storage.
- Autoclaving the cup or scoop is not recommended.

4.4 Impression casting

- The impression should be removed from the disinfectant, rinsed with water and lightly air dried.
- Pour impression, following gypsum manufacturer's directions, within one minute upon removal from the disinfectant and rinsing. If later pouring is necessary, wrap impression in wet towel and store in a sealed plastic bag.
- Do not delay pour more than 5 days. Impression distortion is proportional to the delay in pouring and amount of water evaporation.
- If the impression is to be shipped, use suitable packaging to preclude distortion.

The material is compatible with the following die stones:

ALPHA MODIFIED GYPSUM (Type IV): Glastone®

ALPHA MODIFIED GYPSUM: Alpha Dur 700*

Alpha Gypsum (Type III): Castone

Alpha Gypsum: Calstone, E.P.1*

5. LOT NUMBER AND EXPIRATION DATE

1. Located on pouch. Save empty pouch or record numbers before discarding pouch.
2. Do not use after expiration date. ISO standard used: "YYYY/MM"
3. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date

Manufactured by:
DENTSPLY Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
Tel: 1-302-422-4511
www.dentsply.com

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Distributed by:
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada

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DENTSPLY DeTrey GmbH