1. **Veneer cementation – light-cured**
 Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin/Pretreatment/Adhesive Application. Adhesives compatible with light-cured, CO initiated methacrylates are recommended. If veneer restorations exceed 10mm thickness, are heavily opaquel, or otherwise impede the transmission of light, operator is urged to follow section 2.0. Dual-Cure cementation technique. Apply adhesive to the internal bonding surface of the restoration if directed by adheiv manufacturer's directions for use.
 1) Cementation technique

- 1.1.1 Dispense the desired shade of Calibra® Esthetic Resin Cement base paste from the syringe directly onto the veneer. Protect cement from exposure to light.
- 1.1.2 Seal the loaded syringe in place. Remove excess from the gingival margin with a blunt instrument. Briefly (10 seconds or less) light cure the gingival portion only to seal restoration in place. Remove any excess from proximal and lingual margins. **Technique Tip:** Placement of mylar strips between preparation and adjacent teeth prior to seating veneer aids in isolation and an excess cement clean up. After removing gingival excess and tacking veneer in place, remove interproximal excess cement by pulling mylar strip towards the facial, i.e., from tooth to restoration surface.

- Inadequate polymerization due to insufficient curing**
 - Check compatibility of curing light
 - Verify adequate curing time
 - Check curing output before each procedure

- 1.1.3 Using a visible light, curing unit designed to cure CO initiated methacrylates (spectral output including 470nm), with 11.5 watt/cm² of visible (400nm) light cure 20 seconds each surface of the restoration including proximal and lingual aspects.

- 1.1.4 Following the light cure check and adjust occlusion as necessary. Proceed to Finishing and Polishing section.

2. Crown/bridge/inlay/onlay cementation dual cured

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment/Adhesive Application. Adhesives (plus activators, if applicable) compatible with both light-cured, CO initiated methacrylates and peroxide-amine initiated self-cured methacrylates are recommended. Apply adhesive/activator (if applicable) to the internal bonding surface of the restoration if directed by adhesive manufacturer's directions for use.

2.1 Dispense the desired shade of Calibra® Esthetic Resin Cement base paste from the syringe onto a clean mixing pad. Dispense an equal amount of the desired viscosity of catalyst paste (Regular or High Viscosity). Mix the cement for 20-30 seconds. **Technique Tip:** As with any dual cure resin cement system, the use of a dual cure adhesive system can shorten working time (See Interactions). Calibra® Esthetic Resin Cement base may be mixed with High Viscosity catalyst yielding a High Viscosity "stay-put" thixotropic dual cured cement suitable for inlays, shallow onlays and veneers. Mixing Calibra® Esthetic Resin Cement base with Regular Viscosity catalyst yields an easier flowing dual cured cement suitable for large, multi-surface and full coverage restorations with higher hydraulic strength/pressure.

- Danger of injury due to excessive force**
 - Apply slow and steady pressure on the syringe
 - Do not use excessive force - dual barrel syringe rupture may result

Dual Barreled syringe dispensed

- 1. Remove syringe cap. Dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
- 2. Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mixing tip to with the v-shaped notch on the syringe flange. Turn colored mixing cap 90 degrees in a clockwise direction to lock in place on syringe.
- 3. Gently depress syringe plungers to begin the flow of material. **DO NOT USE EXCESSIVE FORCE.** If force is encountered, remove syringe from operating field; remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.
- 2.1.2 Apply a thin layer of cement on the entire internal surface of the restoration. For inlays/onlays, it may be helpful to apply a 1/4 mil layer of cement to the internal portions of the tooth preparation to allow any porosity or voids. **Technique Tip:** At room temperature, Calibra® Esthetic Resin Cement, in ambient light, offers a minimum work up time of 2 min. 30 sec. on the mixing pad when protected from ambient light. If loaded into restoration lined with clear (or transparent) restorative material, working time will be less (See Precautions).
- 2.1.3 Seal the restoration with gradual pressure. A gentle rocking or vibratory motion may be helpful to insure optimal seating.
- 2.1.4 Remove gross excess from marginal areas. Use an instrument such as a bluntended explorer, periodontal probe or a curet with a broad tip. Restoration should not be removed or torqued during any gross excess cement. Special attention should be paid to interproximal areas using floss to remove excess cement. **Technique Tip:** A 10 second light "pre-cure" of excess cement at the margins will cause cement to "gel", allowing easy cleanup. **Technique Tip:** Apply moderate and consistent pressure to the restoration throughout the self-cure set time of approximately 6 minutes from the beginning of mixing.
- 2.1.5 Once stabilized, light cure, using a visible light, curing unit designed to cure CO initiated methacrylates (spectral output, including 470nm), with a minimum output of 550mW/cm². Cure all marginal areas of the restoration for 20 seconds from each direction – buccal, lingual and occlusal.
- 2.1.6 Following the self-cure set check and adjust occlusion as necessary. Proceed to Finishing and Polishing, section.

3. Endodontic posts

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment/Adhesive Application. Adhesives (plus activators, if applicable) compatible with both light-cured methacrylates and peroxide-amine initiated self-cured methacrylates are recommended. Apply adhesive/activator (if applicable) to the bonding surface of the post if directed by adhesive manufacturer's directions for use.

3) Cementation technique

3.1.1 Dispense the desired shade of Calibra® Esthetic Resin Cement base paste from the syringe onto a clean mixing pad. Dispense an equal base shade e.g. Translucit will allow maximum light transmission. Dispense an equal amount of Regular Viscosity catalyst paste. Mix the cement for 20-30 seconds. **Technique Tip:** As with any dual cure resin cement system, the use of a dual cure adhesive system can shorten working time. (See Interactions)

- Danger of injury due to excessive force**
 - Apply slow and steady pressure on the syringe
 - Do not use excessive force - dual barrel syringe rupture may result

Dual Barreled syringe dispensed

- 1. Remove syringe cap. Dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
- 2. Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mixing tip to with the v-shaped notch on the syringe flange. Turn colored mixing cap 90 degrees in a clockwise direction to lock in place on syringe.
- 3. Gently depress syringe plungers to begin the flow of material. **DO NOT USE EXCESSIVE FORCE.** If force is encountered, remove syringe from operating field; remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.
- 3.1.2 Squeeze mix Calibra® Esthetic Resin Cement components on surface of post and/or into the post preparation with a syringe tip, Lento/Spiral or metal file.
- 3.1.3 Seal post immediately. Clean up excess with appropriate instruments. A 10 second light exposure "pre-cure" of excess cement at the margins will cause cement to "gel", allowing easy cleanup.
- 3.1.4 Following the self-cure set time of approximately 6 minutes from the beginning of mixing, check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.
- 3.1.5 Once stabilized, light cure, using a visible light, curing unit designed to cure CO initiated methacrylates (spectral output, including 470nm), with a minimum output of 550mW/cm². Cure all accessible areas of the post for 20 seconds from each direction – buccal, lingual and/or preparation.
- 3.1.6 Proceed with core build-up and/or preparation.

FINISHING AND POLISHING

- 1. Removal of resin flash is best accomplished with the Enhance® Finishing System of points, cups and discs. The Enhance® System will remove flash and finish restoration margins without removal of trauma to the enamel.
- 2. Polish final restoration using Poso® Polishing System or Prisma®-Gloss® Polishing Paste and Prisma®-Gloss® Extra Fine Polishing Pastes. (See complete Directions for Use supplied with polishing product chosen)

5. **HYGIENE**
 Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

- 4.1 Cleaning and disinfection**
 - For dual-barreled syringes, remove used mixing tip and discard appropriately. Replace original syringe cap prior to storing.
 - The dual-barreled syringe may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.
 - To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid package contamination.
 - Disinfect contaminated syringes with a water-based hospital-level disinfection solution according to national/local regulations.
 - Repeated disinfection may damage label. **NOTE:** Destruction of the label by vigorous wiping. Wipe syringe gently.
 - The supplied mixing spatula may be cleaned by scrubbing with hot water and soap or detergent. Do not autoclave spatula. Disinfect as outlined below.

4.2 Disinfection of re-usable mixing spatula

Disinfect spatula with a hospital-level, tuberculocidal disinfectant solution according to national/local regulations, lodophors, sodium hypochlorite (5.25%), chlorine dioxide, and/or quaternary ammonium disinfectants.
 Some phenolic-based agents and iodochlor-ben products may cause surface staining. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

5. LOT NUMBER AND EXPIRATION DATE

- 1. Do not use after expiration date. ISO standard is use: "YYYY/MM"
- 2. The following numbers should be quoted in all correspondence:
 - Reorder Number
 - Lot Number
 - Expiration date

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