ProPexTM
APEX LOCATOR
User Manual

CONTENTS


GB	3
F	12
D	21
E	30
I	39
P	Este manual de utilização está também disponível, quando requisitado, nas seguintes línguas: português, holandês, dinamarquês, sueco, finlandês e grego.	
NL	Deze gebruiksaanwijzing is, op aanvraag, eveneens verkrijgbaar in devolgende talen: Portugees, Nedelands, Deens, Zweeds, Fins en Grieks.	
FIN	Käyttöohje on saatavana myös seuraavilla kielillä: portugali, hollanti, tanska, ruotsi, suomi ja kreikka.	
S	Denna bruksanvisning finns även att tillgå på följande språk:portugisiska, tyska, danska, svenska, finska och grekiska.	
DK	Denne brugsanvisning kan også rekvireres på følgende sprog: portugisisk, hollandsk, dansk, svensk, finsk og græsk.	
GR		

► Visit our website: www.mallefer.ch

Introduction

GB

Congratulations on the purchase of Propex, an apex locator which is precise, ergonomic, and able to determine the exact working length in seconds.

The measurement of the apex is made at the cemento-dental junction (Foramen Apex), by the Propex apex locator.

Indications

Propex is an electronic device used for apex location determination.

Contraindications

In case where a patient has been fitted with an implanted heart pacemaker (or other electrical equipment) and has been cautioned against the use of small electrical appliances (such as electric shavers, hair dryers etc.) it is recommended that the Propex not be used.

Warnings

See above.

Precautions

In cases where a patient has been fitted with an implanted heart pacemaker (or other electrical equipment) and has been cautioned against use of small electrical appliances (such as electrical shavers, hair dryers, etc.) it is recommended that the Propex not be used. To ensure that short circuits do not impair the measurements, be particularly careful with patients fitted with metal crowns or bridges.

- Make sure that the canal is wet enough to ensure reliability of the measurement.
- Ensure that the file does not touch another instrument.
- Do not expose Propex to any liquid.
- Propex must be stored in normal temperature (< 70°C) and humidity conditions.
- Do not use in the presence of flammable products.
- Teeth with open apices will give imprecise results.

Adverse Reactions: none known

A Contents

Before use, please check the exact contents of the new equipment:

- one Propex apex locator
- one charger
- one measurement cable including 2 lip clips, 2 connection hooks and 2 connection forks.

B Recharging the Battery

The Propex is delivered with a rechargeable battery.

The Propex screen has an indicator showing the level of charge of the battery. When this is flashing, the battery requires recharging. However, it is still functional for 3 or 4 treatments before the battery dies.

Procedure for recharging the battery:

- 1) Disconnect the measurement cable
- 2) Connect the charger cable to the equipment (connection on the right side of the Propex)
- 3) Connect the charger to the mains

An orange light indicates that the battery is charging.
A green light indicates that charging is complete.

Duration of charging: 12 hours (24 hours after long periods without use).

Note: Propex cannot be used while charging.

C Getting Started

- 1) Disconnect the charger from the device if connected
- 2) Connect the measurement cable and switch on the device
- 3) Link the lip clip to the patient
- 4) Insert the file into the canal
- 5) Connect the hook to the file

A green light (coupled with a 2-beep signal) indicates that the Propex is ready for use
A red light (and no sound signal) indicates a faulty connection. Check your connections, clean your hook connection and lip clip, moisten the canal if necessary and start again.

No other specific adjustment is necessary before commencing measurement. Propex works with all types of canal instruments.

One tone audio signal is emitted when the unit is turned on.

D Sound Adjustment

Propex is equipped with a sonic indicator which enables monitoring of the progression of the file within the canal.

This function, in conjunction with the digital display of progression, enables working "blind" while still monitoring the progression of the file.

The volume can be adjusted to one of four levels: mute, low, normal and high, by successive presses on the . volume key.


Visual of different volumes to be included.

E Search for the Apex

1) Median Zone

Slowly introduce the file into the canal.


In the pre-apical zone, Propex indicates the progression of the file via two vertical arrows. A downward pointing arrow appears when the file descends into the canal, and an upward pointing arrow when the instrument is withdrawn.

2) Pre-apical Zone (approx. 2 to 3 mm before the apex)

When the file reaches an area 2 to 3 mm from the apex, Propex emits an audible tone (2 tones).

On the screen, two blinking horizontal arrows appear.


3) Apical Zone


The apical zone is divided into 10 segments graduated arbitrarily from 0.9 to 0.0 (apex) as visual information of file progression. In conjunction with the visual information, Propex emits audible information of file progression via a series of progressive beeps.

When the apex is reached, Propex emits a solid tone.


Apical Zone start

Middle of apical zone


Apex reached

4) Beyond the Apex

A red light (in conjunction with a warning signal) indicates that you have passed the Apex.

F Automatic Shutdown

Propex automatically shuts down after 5 minutes without use. It is advisable, however, to manually switch off equipment after measurement by simply pressing the 0/I key.

G Precautions for Use

- We do not advise the use of Propex on patients fitted with pacemakers (or other electrical equipment) for safety reasons.
- To ensure that short circuits do not impair the measurements, be particularly careful with patients fitted with metal crowns or bridges.
- Avoid excessive liquids inside the tooth cavity to prevent overflow and incorrect measurements.
- Ensure that the file does not touch another instrument.
- Do not expose Propex to any liquid.
- Propex must be stored in normal temperature (< 70°C) and humidity conditions.
- Do not use in the presence of inflammable products.
- Measurements in canals with open apex may provide results with reduced precision.

H Maintenance of your Propex.

The device doesn't contain user serviceable parts. The service and repair should be provided by trained service personnel only.

The equipment may be cleaned with a cloth impregnated with soap and water. Use of chemical agents may cause damage to the equipment.

The accessories (lip clip, connection hook) must be sterilized between treatments (by autoclave at 134 °C).

The measuring cable can't be autoclaved.

I Guarantee

Propex is guaranteed for 24 months from the date of purchase. The accessories (cables, battery etc.) are guaranteed for 6 months.

The guarantee is valid for normal usage conditions. Any modification or accidental damage will render the guarantee void.

J Technical Specifications

Propex complies to IEC60601-1 safety standard and the requirement of CE Marking of Conformity.

CE
0459

Propex electronic apex locator belongs to following category of medical devices:

- Internally powered equipment (3.6 V NiMH rechargeable battery)
- Type LF (Low Frequency) applied parts
- Not suitable for use in the presence of flammable anaesthetic mixtures with air or with oxygen or nitrous oxide
- Continuous operation
- Ingress of liquids – not protected
- Environmental conditions during transportation: temperature: -20°C to +60°C (0 to 140°F); relative humidity: 10% to 90%, non-condensing

Technical specifications

Type of screen: liquid crystal

Supply: NiMH 3.6 V battery

External charger: 120 V/50 Hz or 220V /50 Hz input– 6V DC output

Warning: The charger should comply with the requirements of IEC 60950 standard

K Standard symbols

On the device label appear standard symbols as follows:

Class II equipment

Type BF applied part

Attention, consult accompanying documents

Direct current

Disinfection, cleaning and sterilization procedure for dental instruments and their supports

Foreword

For hygiene and sanitary safety purposes, all medical devices commercialised by Dentsply Maillefer owe to be disinfected, cleaned or yet sterilized before using them into the mouth. Some devices are considered as one usage instruments , anyhow they must follow the same procedure.

Area of application

Cutting instruments, manual or mechanical such as:

Root canal files and broaches, manual or mechanical;

Diamond burs, made of tungsten carbide, stainless steel or carbon steel;

Reamers, drills;

Condensation and filling instruments.

Supports, kits and systems of instrumental organization.

Hand instruments.

Exclusion

Paper points

Standard and Regulation References

ISO 9001, EN 46001, Quality System

European Medical device directive for CE marking, 93/42 EEC of 14th June 1993

ISO 1942-3 Dental vocabulary

ISO/DIS 17664, Sterilization of medical devices.

Procedure

	Operation	Operating mode	Warning	First usage	Second usage and following ones
1	Disinfection	Soak all instruments in a disinfectant solution and then discard broken and bent instruments.	Follow instructions and observe concentrations given by the manufacturer (an excessive concentration may cause cracks on instruments). Observe time of immersion. Do not use products containing Phenol. <u>Instruments made of NiTi:</u> Do not soak more than 5 minutes in a solution of NaOCl at more than 5%.		X
2a	Manual cleaning	Manually brush the instruments in order to eliminate any dental debris. Place them in a support or in a kit.	In the case any dentine tissues are still present on the instrument, use an Air-Flow device.		X
2b.	Manual cleaning	Automatic Cleaning Set the instruments in the kit, and put them in: A) an ultrasonic washing device or B) a washing machine.	Follow instructions and observe concentrations given by the manufacturer. After 5 cycles, degradation of marking may appear. Do not clean the instruments in bulk with an ultrasonic device. Instruments in aluminium : Do not use caustic soda nor solutions with mercury salt nor acid (pH < 6) or alkaline (pH > 8) solutions. If dentine tissues are not totally removed, use an Air-Flow device.	X	X

	Operation	Operating mode	Warning	First usage	Second usage and following ones
3	Rinsing	Rinse the instruments thoroughly with clean, demineralised or distilled water and then dry them.	<u>Instruments made of carbon steel:</u> Must be protected by an anticorrosion agent	X	X
4	Inspection	Sort out and control the instruments, discard those which show any deformation, damage or mark affecting the resistance or performance of the instrument.			X
5	Sterilization	Sterilize the instruments with their support, protected from the environment (in bags).	The instruments and the plastic supports may be sterilized according to the packaging labelling. After 5 cycles, degradation on marking may appear.	X	X
6	with autoclave	124°C during 15 min at 2 bars or 134°C during 5 min at 3 bars	Avoid plastic devices which are not resistant to a temperature of 135°C. <u>Instruments made of carbon steel:</u> Protect them with an anticorrosion agent.	X	X
7	with Chemiclave	132°C during 20 min at 2 bars, soak beforehand instruments on the same solution of Chemiclave in order to avoid any marks	Recommended technique for carbon steel instruments. Avoid plastic devices which are not resistant to a temperature of 135°C	X	X
8	with hot air	180°C during 120 min at atmospheric pressure	Avoid plastic devices which are not resistant to a continuous temperature of 180°C	X	X

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Dentsply France SAS