# Neo Spectra™ ST flow

## Flowable Composite Restorative

Composite fluide de restauration

Composite restaurador fluido

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CAUTION: For dental use only.

USA: Rx only.

Content	Page
1 Product description.....	2
2 Safety notes .....	3
3 Step-by-step instructions .....	5
4 Hygiene.....	9
5 Lot number, expiration date and disposal.....	10

### 1 Product description

**Neo Spectra™ ST flow** composite restorative is a visible light cured, radiopaque composite resin restorative material with flow characteristics that make it ideal for use in anterior and posterior teeth. **Neo Spectra™ ST flow** composite restorative adapts to the cavity wall without the use of hand instruments. The unique "stackable" flow characteristics allow intimate adaptation without slumping.

**Neo Spectra™ ST flow** is a hybrid composite possessing the high strength and wear resistance required by anterior and certain posterior restorations, combined with an easily achieved high surface luster and smoothness.

#### 1.1 Indications

Neo Spectra™ ST flow composite restorative is indicated for the following procedures:

- Direct restoration of cavities (except occlusal stress bearing Class II cavities), lesions.
- Filling of defects and undercuts in crowns, inlays, and onlays.
- Blockouts.
- Repair of defects.
- Pit and fissure sealants.
- Cementation of light transmissible indirect restorations.

#### 1.2 Contraindications

- Neo Spectra™ ST flow composite restorative is contraindicated for use with patients who have a known hypersensitivity to methacrylate resins.

#### 1.3 Delivery forms

Neo Spectra™ ST flow composite restorative is available in:

- Predosed Compula® Tips
- Syringes

## 1.4 Composition

- The resin matrix of Neo Spectra™ ST flow composite restorative consists of:
- Urethane modified BisGMA-adduct
  - BisEMA and diluents
  - Camphorquinone (CQ) Photoinitiator
  - Stabilizers
  - Pigments

### Filler system based on SphereTEC®.

The filler matrix of Neo Spectra™ ST flow composite restorative consists of barium-aluminum-borosilicate glass, ytterbium fluoride and iron oxide pigments and titanium oxide pigments according to shade.

The filler system based on SphereTEC® fillers ranges from 0.1 to 3.0 µm.

Including catalysts, additives, stabilizers and pigments (~1 weight-%, according to shade), the total filler is:

62.5 weight-% or 40.0 volume-%.

## 1.5 Compatible adhesives

Neo Spectra™ ST flow composite restorative is used following application of a suitable dentin/enamel adhesive and is chemically compatible with conventional (meth)acrylate-based dentin/enamel adhesives including Dentsply Sirona adhesives designed for use with visible light cured composite restoratives (see complete Directions for Use of selected adhesive).

## 2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Directions for Use.

### Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

## 2.1 Warnings

The material contains methacrylates and polymerizable monomers which may be irritating to skin, eyes and oral mucosa and may cause sensitization by skin contact and allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes rinse with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, red-dish rashes may be seen on the skin. If contact with skin occurs, remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

## **2.2 Precautions**

- This product is intended to be used as specifically outlined in these Directions for Use. Any use of this product inconsistent with these Directions for Use is at the discretion and sole responsibility of the practitioner.
- If refrigerated, allow material to reach room temperature prior to use.
  - Contact with saliva and blood during composite placement may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
  - Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
  - Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
  - The syringes cannot be reprocessed. To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringes are handled with clean/disinfected gloves. Discard syringes if contaminated.
  - As additional precautionary measure, syringes may be protected from gross debris but not from all contamination by applying a protective barrier.
  - Use of Compula® Tips with the Dentsply Sirona Compules® Tips Gun is recommended.
  - The Compules® Tips Gun is only intended to be used with Dentsply Sirona Compules® or Compula® Tips.
  - For further information please refer to the Dentsply Sirona Compules® Tips Gun Directions for Use.
  - The shade guide and the individual shade guide tabs cannot be reprocessed. To prevent the shade guide and the individual shade guide tabs from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the shade guide and the individual shade guide tabs are handled with clean/disinfected gloves.
  - The material should extrude easily. DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause the Compula® Tip to rupture or to eject from the Dentsply Sirona Compules® Tips Gun.
  - Tightly close syringes immediately after use. Do not allow dispensing tip to remain as a syringe cap. Replace original cap tightly after each use. Discard and properly dispose dispensing tips after use, as the tips may clog if material is allowed to dry or set inside.
  - Dentsply Sirona supplies the appropriate dispensing tip for the Neo Spectra™ ST flow composite restorative material syringe. This tip is the only tip that should be utilized for placement of the material.
  - Insufficient data exist to support the use in Class I cavities with isthmus widths beyond central groove areas, or when replacing a functional cusp. Excessive wear or restoration failure may result.
  - Insufficient data exist to support the use in Root canal fillings. Restoration failure may result.
  - Interactions:
    - Do not use eugenol- or hydrogen peroxide-containing materials in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
    - The Neo Spectra™ ST flow composite restorative is light cured material. Therefore, it should be protected from ambient light. Proceed immediately once material has been placed.
    - If mineral-impregnated (e.g., ferric compounds) retraction cords and/or hemostatic solutions are used in conjunction with adhesive procedures, marginal seal may be adversely affected, allowing microleakage, subsurface staining and/or restoration failure. If gingival retraction is necessary, use of plain, non-impregnated cord is recommended.

## **2.3 Adverse reactions**

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see Warnings).

## 2.4 Storage conditions

- Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.
- Store in a well ventilated place at temperatures between 2 °C and 28 °C (36 °F and 82 °F).
  - Keep out of direct sunlight and protect from moisture.
  - Do not freeze.
  - Do not use after expiration date.

## 3 Step-by-step instructions

### 3.1 Direct Application – Cavities, Lesions, Repair

#### 3.1.1 Shade selection

1. Prior to shade selection, clean the tooth surface with prophylaxis paste to remove plaque or surface stain.
2. Select the shade while teeth are hydrated.

Neo Spectra™ ST flow composite restorative has been scientifically designed to allow the tooth structure enamel and dentin to blend together with the composite with a life-like translucent result. The universal CLOUD shades of Neo Spectra™ ST flow composite restorative adequately reproduce the traditional 16 VITA<sup>®</sup> shades.

- A1 (= VITA<sup>®</sup> shades A1, B1, C1)
- A2 (= VITA<sup>®</sup> shades A2, B2, D2)
- A3 (= VITA<sup>®</sup> shades A3, C2, D3, D4)
- A3,5 (= VITA<sup>®</sup> shades A3,5, B3, B4, C3)
- A4 (= VITA<sup>®</sup> shades A4, C4)

A Bleach White shade (BW), as well as opaque dentin shades D1 and D3 and translucent enamel shade E1 are also available.

For shade selection, there are two tools that may be used:

- The VITA<sup>®</sup> classical shade guide in combination with the i-shade label. The Spectra ST i-shade label is applicable for the Neo Spectra™ ST flow composite restorative and facilitates shade selection by assigning one of five universal CLOUD shades to each of the 16 VITA<sup>®</sup> shades (see table below).

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- The Spectra ST composite restorative shade guide consists of 9 shade tabs, which are manufactured from composite with optical properties close to those of Neo Spectra™ ST flow composite restorative.

#### Avoid contamination of the shade guide.

To prevent the shade guide from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the shade guide is handled offsite the dental unit with clean/disinfected gloves. For shade selection, remove individual tabs from shade guide holder.

<sup>1</sup> Not a registered trademark of Dentsply Sirona, Inc.

### **3.1.2 Cavity preparation**

1. Prepare the cavity so that no residual amalgam or restorative is left (unless repairing an otherwise sound restoration).
2. Use adequate isolation such as rubber dam.
3. Rinse surface with water spray and carefully remove rinsing water. Do not desiccate the tooth structure.

### **3.1.3 Matrix placement and wedging**

1. Place a matrix (e.g., Mylar, AutoMatrix® matrix system or Palodent® matrix system) and wedge. Burnishing of the matrix band will improve contact and contour. Pre-wedging or BiTine® ring placement is recommended.

### **3.1.4 Pulp protection, tooth conditioning/dentin pre-treatment, adhesive application**

Refer to adhesive manufacturer's directions for pulp protection, tooth conditioning and/or adhesive application. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately to placement of restorative material.

### **3.1.5 Application – Delivery**

#### **Using syringes**

##### **Avoid contamination of the syringe.**

To prevent the syringe from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringe is handled offside the dental unit with clean/disinfected gloves.

1. Remove cap from the end of the syringe. To assure free flow of material from syringe, express a small amount onto pad, away from the patient field.
2. Attach disposable dispensing tip to end of the syringe. Turn tip clockwise 1/4 to 1/2 turn to assure that it is fully seated. Tug on tip to be sure that it is locked into the collar of the syringe.
3. Neo Spectra™ ST flow composite restorative should flow freely with gentle pressure. DO NOT USE EXCESSIVE FORCE. If more than gentle pressure is required, remove from patient field and check for obstruction.
4. DISCARD AND PROPERLY DISPOSE OF DISPENSING TIP IMMEDIATELY AFTER USE. REPLACE ORIGINAL CAP. Do not store syringe with gray dispensing tip in place. STORE ONLY WITH ORIGINAL CAP.

**Note:** It is recommended to pull back slightly on the syringe plunger after use to prevent excessive flow of material.

#### **Using Compula® Tips**

##### **Excessive force.**

Injury.

1. Apply slow and steady pressure on the Compules® Tips Gun.
2. Do not use excessive force. Compula® Tip rupture or ejection from Compules® Tips Gun may result.

1. Insert Compula® Tip into the notched opening of the Compules® Tips Gun. Be certain that the collar on the Compula® Tip is inserted first.
2. Remove the colored cap from the Compula® Tip. The Compula® Tip may be rotated 360° to gain the proper angle of entrance into the cavity.
3. Dispense the material into the cavity preparation using a slow, steady pressure.  
DO NOT USE EXCESSIVE FORCE.
4. To remove the used Compula® Tip, be sure that the Compules® Tips Gun plunger is pulled back completely by allowing the handle to open to its widest position. Apply a downward motion to the front end of the Compula® Tip and remove.

### **3.1.5.1 Applying restorative**

1. Dispense Neo Spectra™ ST flow composite restorative directly into cavity preparation/tooth surface from the Compula® Tip or syringe tip using slow, steady pressure. Avoid lifting the tip while dispensing to minimize air entrapment. Any visible air bubbles should be pierced with a clean, sharp explorer prior to curing. Material may be placed and light cured in increments up to 2 mm. In deeper preparations it is recommended to double the curing time to compensate for lower irradiance (see Curing, Step 3.1.6). For maximum light transmission and to minimize finishing, avoid overfilling cavity.

### **3.1.6 Light curing**

1. Light cure each area of the restoration surface with a suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. peak spectral output 440–480 nm. Cure universal CLOUD shades, BW, and E1 increments for 10 seconds, and D1 and D3 shades for 30 seconds when using a curing light that has a minimum irradiance of 800 mW/cm<sup>2</sup>. Cure universal CLOUD shades, and BW for 20 seconds, E1 for 10 seconds, and D1 and D3 shades for 40 seconds when minimum irradiance is between 550 mW/cm<sup>2</sup> and 800 mW/cm<sup>2</sup>. Refer to curing light manufacturer's Directions for Use for compatibility and curing recommendations.

#### **Curing time for 2 mm increments**

<b>Neo Spectra™ ST flow Curing Recommendations</b>		
<b>Shades</b>	<b>Irradiance mW/cm<sup>2</sup></b>	<b>Curing Time</b>
<b>A1-A4 and BW</b>	≥ 550	20 seconds
	≥ 800	10 seconds
<b>D1, D3</b>	≥ 550	40 seconds
	≥ 800	30 seconds
<b>E1</b>	≥ 550	10 seconds

#### **Insufficient curing.**

Inadequate polymerization.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check minimum irradiance.
4. Cure each area of each increment for the recommended curing time.
5. Check distance to surface to be cured.

### **3.1.7 Finishing and polishing**

1. Contour the restoration using finishing burs or diamonds.
2. Use Enhance® finishing devices for additional finishing.
3. Prior to polishing, check that surface is smooth and free of defects.
4. For achieving a very high luster of the restoration, we recommend Enhance® and Enhance® PoGo®, finishing and polishing instruments, Prisma® Gloss™ Regular and Extrafine polishing paste.

For contouring, finishing, and/or polishing, follow the manufacturer's Directions for Use.

All shades of Neo Spectra™ ST flow composite restorative are radio-opaque, with a 1 mm radio-opacity equivalent to 1.8 mm radio-opacity of aluminum. Aluminum has a radio-opacity equivalent to that of dentin. Thus 1 mm of material having a radio-opacity equivalent to 1 mm of aluminum has a radio-opacity equivalent to that of dentin.

### **3.2 Cementation of light transmissible indirect restorations**

- I. Following removal of the temporary restoration and any remaining temporary cement, clean enamel and dentin as directed by the adhesive manufacturer's directions.
- II. Check the fit and esthetics of the restoration. Technique Tip: Occlusal adjustment of veneers and onlays is best accomplished after final cementation.

#### **3.2.1 Treatment of the restoration**

Ceramic/Composite Restorations

Follow the dental laboratory or restoration manufacturer's instructions for pre-treatment, if required. Restorations designed to be silanated or if the internal silanated surface has been disturbed during try-in, apply Calibra® Silane Coupling Agent (available separately) according to the manufacturer's instructions.

#### **3.2.2 Tooth conditioning/dentin pre-treatment/adhesive application**

1. Proceed with adhesive application as directed in the complete Directions for Use.

##### **Delayed or reduced adhesion.**

Adhesive incompatibility.

1. Select compatible light cure adhesive.
2. Follow selected adhesive Directions for Use.

2. Always refer to and follow adhesive manufacturer's complete Directions for Use for selection and application of the appropriate light cure adhesive. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately following adhesive application to cementation with Neo Spectra™ ST flow material.

#### **3.2.3 Cementation (light cured)**

Before proceeding refer to above sections for Treatment of the restoration and Tooth conditioning/dentin pre-treatment/adhesive application.

- Adhesives compatible with light cured, CQ initiated methacrylates are recommended.
- If restorations exceed 1.0 mm thickness, are heavily opaqued, or otherwise impede the transmission of light, operator is urged to select Calibra® Esthetic Resin Cement and follow the Dual-Cured cementation technique, mixing base with catalyst (available separately).
- Apply adhesive to the internal bonding surface of the restoration if directed by adhesive manufacturer's Directions for Use.

1. Dispense the desired shade of Neo Spectra™ ST flow material from the syringe directly onto the restoration. Protect cement from exposure to light.
2. Seat the loaded restoration in place. Remove excess from the margins with a blunt instrument. Briefly (10 seconds or less) light cure the gingival portion only to tack restoration in place. Remove any remaining excess from margins. Technique Tip: Placement of mylar strips between preparation and adjacent teeth prior to seating veneers aids in isolation and in excess cement clean up. After removing gingival excess and tacking veneer in place, remove interproximal excess cement by pulling mylar strip towards the facial, i.e., from tooth to restoration surface.

### 3.2.4 Curing

#### Insufficient curing.

Inadequate polymerization.

- 1. Check compatibility of curing light.
- 2. Check curing cycle.
- 3. Check minimum irradiance.
- 4. Cure each area of each increment for the recommended curing time.
- 5. Check distance to surface to be cured.

1. Light cure with a suitable visible light curing unit designed to cure materials containing camphorquinone (CQ) initiator, i.e. peak spectral output 440-480 nm. Cure each area of the restoration surface (buccal, lingual, and interproximal) for at least 20 seconds at an irradiance of minimum 800 mW/cm<sup>2</sup>. To ensure proper irradiance follow curing light instructions (e.g., increase exposure time to compensate lower irradiance due to increased distance between light guide/tip and restorative surface or curing lights with lower radiant emittance).
2. Following the light cure check and adjust occlusion as necessary. Proceed to finishing and polishing as outlined in section 3.1.7 above.

## 4 Hygiene

#### Cross-contamination.

Infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Syringes cannot be reprocessed. Dispose of contaminated syringes in accordance with local regulations.
- Reprocess reusable products according to the instructions.

### 4.1 Syringes – cross-contamination

#### Cross-contamination.

Infection.

- The syringes cannot be reprocessed.
- To prevent the syringes from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the syringes are handled with clean/disinfected gloves. Do not reuse syringes if contaminated.
- Dispose of contaminated syringes in accordance with local regulations.

To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended. The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination.

Incidental contact of the syringe with water, soap or a water-based hospital-level disinfection solution will not damage syringe body. Do not allow any solution contact with contained material. Discard composite material that has been in contact with any fluid or non-sterile instrument.

Repeated liquid contact may damage label. Dry the syringe with a lint-free single use cloth.

**NOTE: Vigorous wiping can destroy the label.** Wipe syringe gently.

#### 4.2 Compules® Tips Gun

For reprocessing instructions please refer to the Directions for Use of the Compules® Tips Gun, which is available on our webpage at [www.dentsplysirona.com](http://www.dentsplysirona.com) and [www.dentsply.eu/IFU](http://www.dentsply.eu/IFU). If requested, we will send you a free printed copy of the Directions for Use in the language you require within 7 days. In the USA, call 1-800-532-2855. Outside North America, use the order form provided on [www.dentsply.eu/IFU](http://www.dentsply.eu/IFU) for this purpose.

#### 4.3 Compula® Tips – cross-contamination

##### Cross-contamination.

###### Infection.

- Do not reuse Compula® Tips. Dispose of Compula® Tips in accordance with local regulations.

#### 4.4 Shade guide and individual shade guide tabs – cross-contamination

##### Cross-contamination.

###### Infection.

- The shade guide and the individual shade guide tabs cannot be reprocessed.
- To prevent the shade guide and the individual shade guide tabs from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the shade guide and the individual shade guide tabs are handled with clean/disinfected gloves. Do not reuse shade guide and the individual shade guide tabs if contaminated.
- Dispose of contaminated shade guide and the individual shade guide tabs in accordance with local regulations.

### 5 Lot number (), expiration date () and disposal

1. Do not use after expiration date.

ISO standard is used: "YYYY-MM" or "YYYY-MM-DD".

2. The following numbers should be quoted in all correspondence:

- Reorder number
- Lot number
- Expiration date

3. Dispose in accordance with local regulations.

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© Dentsply Sirona 2018-12-10
Dentsply DeTrey GmbH