# cercon<sup>®</sup> ht

CE 0124 RX<sup>only</sup>

Gebrauchsanweisung  
Instructions for use - multilingual

 Dentsply  
Sirona

<b>AE</b>	Cercon® ht تَعْلِيمَاتُ الْإِسْتِعْمَال .....	1
<b>CZ</b>	Cercon® ht Návod k použití .....	11
<b>DE</b>	Gebrauchsanweisung Cercon® ht .....	21
<b>DK</b>	Cercon® ht Brugsanvisning .....	31
<b>ES</b>	Instrucciones de uso de Cercon® ht .....	41
<b>FR</b>	Notice d'utilisation Cercon® ht .....	51
<b>GB</b>	Cercon® ht Instructions for Use .....	61
<b>GR</b>	Cercon® ht Οδηγίες Χρήσης .....	71
<b>HR</b>	Cercon® ht Upute za upotrebu .....	81
<b>HU</b>	Cercon® ht Használati útmutató .....	91
<b>IT</b>	Cercon® ht - Istruzioni per l'uso .....	101
<b>JP</b>	Cercon® ht 取扱説明書 .....	111
<b>LT</b>	„Cercon® ht“ naudojimo instrukcija .....	121

<b>LV</b>	Cercon® ht Lietošanas instrukcija .....	131
<b>NL</b>	Cercon® ht gebruiksaanwijzing .....	141
<b>NO</b>	Cercon® ht Bruksanvisning .....	151
<b>PL</b>	Instrukcja postępowania - Cercon® ht .....	161
<b>PT</b>	Cercon® ht Instruções de uso .....	171
<b>RC</b>	泽康高透氧化锆使用说明 .....	181
<b>RO</b>	Instructiuni de utilizare Cercon® ht .....	191
<b>RU</b>	Cercon® ht Инструкция по применению .....	201
<b>SE</b>	Cercon® ht användningsinstruktioner .....	211
<b>SI</b>	Navodila za uporabo Cercon® ht .....	221
<b>SK</b>	Cercon® ht Návod na používanie .....	231
<b>TR</b>	Cercon® ht Kullanım için talimatlar .....	241

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GB

## Product description:

Cercon ht blanks are made of yttrium oxide- (yttria-) stabilized zirconium oxide (zirconia) (Y-TZP). They are used in fabricating frameworks for fixed prosthetic restorations.

The material is an oxide ceramic material characterized by its particularly high strength. Depending on the framework design, Cercon ht frameworks can be ceramically veneered or delivered as fully contoured restorations. Which blank is selected will depend on the tooth shade to be reproduced and the space available for the veneer. With fully contoured restorations, no space is required for the ceramic veneer, which may allow the dentist to preserve more of the tooth substance during preparation.

Framework material	Zirconium oxide (Y-TZP)
Temporary cementation	Possible (for fully contoured restorations)
Definitive cementation	<ul style="list-style-type: none"><li>Adhesive cementing</li><li>Conventional cementing</li></ul>

*Objects are individually fabricated to your digital design specifications such as anatomic contour, framework and wall thickness, connector diameter and cementing gap.*

## Technical specifications:

- Type II, class 5 (pursuant to DIN EN ISO 6872:2015)
- CTE: 10.5 µm/m · K (25–500 °C)
- Modulus of elasticity: 210 GPa
- Flexural strength: approx. 1200 MPa (three-point flexural testing)

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- 2-piece abutments\*
- Conical\* and telescope primary crowns
- Crowns
- Multi-unit bridges (with no more than two pontics between abutment crowns)\*\*
- Bridges with cantilevered pontics, with min. 2 abutment teeth (up to the second premolar)

\* Not valid for the U.S.

\*\* For Canada: Restricted to 6 units or less

**Contraindications:**

- This medical product may not be used in patient hypersensitivity to zirconia (Y-TZP) or one of the other ingredients
- Bruxism or recalcitrant parafunctional habits (for ceramically veneered frameworks)
- Insufficient available space
- Endodontic posts
- Endosseous implants
- Inlay bridges

**Warnings:**

Possible cross-reactions or interactions of this medical product with other medical products or material already present in the oral environment must be taken into consideration by the dentist when selecting this medical product.

**Precautions:**

Please note:

- Keep product dust away from eyes.
- Avoid any contact with mucosa.
- After use, wash your hands and apply a hand cream.
- Do not smoke, eat or drink while handling the product.
- Do not swallow the product.
- Do not inhale dust particles during grinding.
- Use local vacuum suction and suitable mouth/face protection during manual machining at the workplace.

The safety and warning notes listed here describe how to use our medical product in a safe and risk-free manner. Notify the dentist in charge of all factors described above if you use this medical product for a custom design and make sure to comply with the pertinent Material Safety Data Sheets (MSDS).

**Adverse effects:**

If properly processed and used, adverse effects of this medical product are highly unlikely. However, reactions of the immune system (such as allergies) to substances contained in the material or localized paraesthesia (such as taste disturbances or irritation of the oral mucosa) cannot be completely ruled out as a matter of principle. Should you hear or be informed of any adverse effects – even when doubtful – we would like to request notification.

## Special notes on framework design

### **Framework design for fully contoured restorations:**

Indicated especially where limited occlusal space is available, for frameworks that are not to be veneered or for frameworks with body stains.

The surface of fully contoured frameworks can carefully be optimized with rotary tools such as fine cutters prior to sintering.

Make sure that the occlusal surface is not changed by subsequent deepening of fissures, as notching action may compromise the strength of the material. Please note that flat occlusal reliefs can extend the live expectancy of fully contoured restorations. When making manual adjustments, make sure never to separate the interdental spaces of the frameworks with cutting discs or other rotary instruments. Doing so might damage the framework and compromise the strength of the material!

### **Important note:**

Please ensure that the minimum wall thickness of the framework in the area of the occlusal surface is respected even after occlusal adjustments.

### **Framework design for veneered restorations:**

Frameworks to be ceramically veneered are designed to reduced anatomical contour to provide maximum support for the veneer.

**The frameworks can be veneered using the press-on technique or the build-up technique.**

#### **Wall and border thickness:**

Wall thickness, single crowns*	0.4 mm
Margin thickness, single crowns	0.2 mm
Wall thickness, bridges*	0.5 mm
Margin thickness, bridges	0.2 mm

\* For monolithic crowns and bridges treatment add in each case 0.1 mm reserve for grind in.

#### **Additional dimensional requirements for the anterior region:**

Number of pontics	2
Connector cross-section	6 mm <sup>2</sup>

#### **Additional dimensional requirements for the posterior region:**

Number of pontics	2
Connector cross-section	9 mm <sup>2</sup>
Cantilever pontic at tooth position (only one pontic, up to one premolar width)	Up to the second premolar
Connector cross-section for this cantilever pontic	12 mm <sup>2</sup>

## Nesting

---

### Special nesting notes for large-span bridges (9 units or more)

For distortion-free sintering of long-span bridges (9 units or more) when nesting the object in a Cercon ht 98 disk, make sure to add a sintering reinforcement in the form of a "tongue" to create.

Always place objects (single copings, three-unit frameworks) within the sintering reinforcement ("tongue") in order to achieve a homogeneous shrinkage of large-span objects during sintering.

### Entering the shrinkage factor for the milling unit:

In CAM software allowing the entry of 3 dimensions, enter X, Y and Z values.

In CAM software allowing the entry of 2 dimensions, enter X, Y or Z values.

In CAM software allowing the entry of 1 dimension, enter the X value.

## Finishing

---

### Special notes on finishing:

Please read the respective Instructions for Use for your device regarding further processing.

We recommend the following milling strategies for secure processing:

Procedure	Tool, spherical	Feed rate Z mm/min	Feed rate F mm/min	Speed rpm	$a_p$ mm	$a_e$ mm	Strategy	Dimension
Contour rough-machining, occlusal side (OS)	HM Ø2	800	1800	22000	0.8	1	Contour-parallel rough-machining	0.3
Contour rough machining, cavity side (CS)	HM Ø2	800	1800	22000	0.8	1	Contour-parallel rough-machining	0.3
Pre-finishing (CS)	HM Ø2	400	1200	22000		0.2	Residue material rough-machining	0.1
Pre-finishing (OS)	HM Ø2	400	1200	22000		0.2	Residue material rough-machining	0.1
Finishing (OS)	HM Ø 1	400	1200	27000		0.15	3D offset	0
Finishing without cavity (CS)	HM Ø 1	400	1200	27000		0.15	3D offset from inside	0
Finishing cavities	HM Ø 1	250	1000	27000		0.1	3D offset	0
Fissures	HM Ø 0,5-0,6	250	1000	35000		0.1	Partial finishing, 3D	0

*The milling strategies are recommendations. If necessary, carry out trial millings and adapt the parameters as needed.*

## Separating

---

### Notes on separating the objects:

Separate the objects from the disk by sandblasting with aluminium oxide (50 µm, max. 1.5 bar). The sandblasting aid facilitates the removal of objects from the blank. It helps prevent framework fracture or other damage to the objects. For large-span bridges (9 units or more), separate only the labial and buccal sprues of the objects and the "tongue" connector, because the objects must be sintered together with that "tongue". Any protruding ridges at the bottom of the "tongue" must be removed in order to assure that the objects will stay firmly put on the sintering tray. The smaller objects nested inside the "tongue" are detached completely and sintered separately.

Sandblasting aid

Milling and sandblasting the object

## Sintering

---

### Sintering in the Cercon heat plus P8:

- 1500 °C in the Cercon heat plus P8
  - Program #4 for bridges up to 8 units,  $T_{max} = 1500 \text{ }^{\circ}\text{C}$
  - Program #5 for bridges of 9 units or more,  $T_{max} = 1500 \text{ }^{\circ}\text{C}$

### Special sintering notes for large-span bridges (9 units or more)

Two wide-span bridges (9 or more units) can be sintered in the Cercon heat plus P8 at one time. Place the objects on the sintering block with due regard to the internal vertical clearance of the Cercon heat plus P8 (130 mm) and the need to facilitate contraction without any mechanical constraints. The objects must not touch the sintering block during the sintering.

Sintertray with sintering block

Correct position on the sintering tray

## **Sintering in the heat DUO/Multimat2Sinter:**

- Program #6: Speed program for bridges up to 6 units,  $T_{max} = 1540^{\circ}\text{C}$
- Program #7: Standard sintering for bridge frameworks up to 8 units,  $T_{max} = 1520^{\circ}\text{C}$
- Program #8: Sintering program for bridges of 9 units or more,  $T_{max} = 1520^{\circ}\text{C}$

## **Special sintering notes for large-span bridges (9 units or more)**

Two wide-span bridges (9 or more units) can be sintered in the heat DUO or Multimat2Sinter at one time. Place the objects on the special sintering bar for large-span bridges and keep in mind that the objects must not touch the sintering block during the sintering.

Note the maximum dimensions for milling:

Height: 65 mm

Width: 90 mm

Placement of bridges with  
sintering bar (9 or more units)

## **Sintering in third-party furnaces:**

Sintering results may be adversely affected by e.g.:

- Incorrect sintering temperatures
- Insufficient heating power
- Incorrect temperature curves
- Incorrect object placement
- Insufficient heat-storing capacity of the furnace over the time of the sintering cycle
- Manufacturer-related or age-related variations in heating performance
- Object contamination by oxidation products emitted by non-enclosed heating elements

Any of these factors by itself or in combination may reduce the maximum strength of our above-mentioned zirconium dioxide materials and compromise the life expectancy of the frameworks.

For these reasons, we cannot grant general approval for the use of third-party furnaces for sintering two-piece abutments (mesostructures) and crown and bridge frameworks made of Cercon ht. We will, however, open the system technically for the use of furnaces of third-party furnaces only on condition that the following requirement are met:

**IMPORTANT!**

The sintering of **Cercon ht** in third-party sintering furnaces is done at the user's own responsibility and risk. DeguDent can accept no responsibility for any kind of damage caused by the sintering of zirconia materials in third-party furnaces, including but not limited to damage to or damage caused by the sintered objects such as abutments, crowns or bridge frameworks.

**Additional notes:**

The programming you use for your sintering furnace should be analogous to the DeguDent sintering programs.

**General Sinter programs Cercon ht all shades**

Material:	Start-Temp.	Ramp time	Temp. 1	Holding time	Ramp time	Temp. 2	Holding time	Cooling
	°C	min	°C	min	min	°C	min	
Cercon ht, all shades								
Sinter program for bridge frameworks up to 8 units	RT <sup>1)</sup>	40	900	0	55	1500	145	With closed furnace cooling down to 200 °C
Sinter program for bridge frameworks for 9 or more units	RT <sup>1)</sup>	120	860	0	320	1500	120	With closed furnace cooling down to 200 °C
Speed-Sinter program für bridge frameworks up to 6 units	RT <sup>1)</sup>	90	1540 <sup>2)</sup>	35	20	1150	0	Gradual opening of the furnace within 35 Min down to 200 °C

1) Room temperature

2) valid for closed sinter bowls, otherwise 1520°C

Material:	Start-Temp.	Ramp time	Temp. 1	Holding time	Ramp time	Temp. 2	Holding time	Cooling
	°C	°C/min	°C	min	°C/min	°C	min	
Cercon ht, all shades								
Sinter program for bridge frameworks up to 8 units	RT <sup>1)</sup>	22	900	0	11	1500	145	With closed furnace cooling down to 200 °C
Sinter program for bridge frameworks for 9 or more units	RT <sup>1)</sup>	7	860	0	2	1500	120	With closed furnace cooling down to 200 °C
Speed-Sinter program für bridge frameworks up to 6 units	RT <sup>1)</sup>	17	1540 <sup>2)</sup>	35	18	1150	0	Gradual opening of the furnace within 35 Min down to 200 °C

1) Room temperature

2) valid for closed sinter bowls, otherwise 1520°C

Material:	Start-Temp.	Ramp time	Temp. 1	Holding time	Ramp time	Temp. 2	Holding time	Cooling
	°C	°C/h	°C	h:min	°C/h	°C	h:min	
Cercon ht, all shades								
Sinter program for bridge frameworks up to 8 units	RT <sup>1)</sup>	1320	900	00:00	660	1500	02:15	With closed furnace cooling down to 200 °C
Sinter program for bridge frameworks for 9 or more units	RT <sup>1)</sup>	420	860	00:00	120	1500	02:00	With closed furnace cooling down to 200 °C
Speed-Sinter program für bridge frameworks up to 6 units	RT <sup>1)</sup>	1020	1540 <sup>2)</sup>	00:35	1080	1150	00:00	Gradual opening of the furnace within 35 Min down to 200 °C

1) Room temperature

2) valid for closed sinter bowls, otherwise 1520°C

Sintering programs, Multimat2Sinter / heat DUO / Sirona HTC-speed sintering furnace:

Cercon ht speed sintering of bridge frameworks with up to 6 units:

Step	Heating rate (°C/min)	Temperature (°C)	Holding time (min)
S4	70	750	0
S3	70	1100	0
S2	70	1540	35
S1	70	0	0

Cercon ht standard sintering:

Step	Heating rate (°C/min)	Temperature (°C)	Holding time (min)
S4	31	300	0
S3	30	750	0
S2	11	1520	130
S1	22	880	0

Cercon ht sintering of bridge frameworks with 9 or more units:

Step	Heating rate (°C/min)	Temperature (°C)	Holding time (min)
S4	11	170	0
S3	11	1100	0
S2	2	1520	125
S1	7	820	0

*Sintering temperatures are recommendations. If necessary, carry out a trial sintering cycle and adapt the sintering temperatures or times as needed.*

## **Separating the sintering reinforcement in the case of large-span bridges:**

Objects are separated from the "tongue" after sintering using irrigated rotary diamond cutters.

### **Manual finishing after sintering:**

- Sandblast the inside and outside of the framework with aluminium oxide (110–125 µm, max. 2–3 bar, 45° angle).
- Eliminate premature contacts one by one until the framework has reached its final position on the die(s).
- During the try-on and fitting of the framework, keep the dies on the cast and try on the framework as a whole.
- Once the try-on and fitting has been completed, do not perform any additional adjustments such as finishing the entire framework.

**Note:** Zirconia crown or bridge abutments should fit passively, without friction. The rationale for requiring a passive fit if zirconia frameworks lies in the physical properties of the material itself: Ceramics will tolerate compressive stress but not tensile stress. In crowns that exhibit a friction fit, this friction is generated by certain "bearing" parts of the framework, as the crown itself never rests completely on the die due to processing-related surface roughness (the same actually holds for cast-metal frameworks). So it is only the peaks of the surface roughness (the named "bearing" parts) that are in contact with the die surface. This converts the compressive forces generated by masticatory pressure into tensile forces, which may damage the crown.

Adapt the framework margins exactly to the preparation margins.

**Note:** Sintered zirconia should be finished using diamond instruments under proper irrigation only. Keep the pressure on the framework material to a minimum and work only in one direction.

- Sandblast the finished areas once again with aluminium oxide (110–125 µm, max. 2–3 bar, 45° angle).
- Finally, clean the framework using a steam cleaner.

### **Veneering/staining:**

Adjust your build-up technique to the respective Cercon ht shade.

Note that the light transmission of Cercon ht is very high due to its translucency.

### **General**

The reproduction of the individual patient's tooth colour can be significantly influenced by:

- Shade of the die
- Shade of the material used for temporary or definitive cementing
- Framework wall thickness
- Liner, if applied

### **Staining technique**

For tooth-coloured staining of fully contoured Cercon ht restorations we recommend Cercon ceram body stains, Cercon TCT Stains and Celtra Stains.

### **Layering technique**

We recommend Cercon ceram Kiss/Cercon ceram press/Cercon ceram love ceramic veneering materials for veneering zirconia frameworks. Make sure to follow the applicable Instructions for Use.

### **Annealing:**

Based on the results of our scientific testing of zirconia frameworks, we consider a separate annealing ("heating") step both unnecessary and inappropriate.

### **Polishing in the laboratory:**

Unveneered Cercon ht frameworks should be highly polished or glazed with a high-gloss ceramic glaze. This also makes the restoration more amenable to proper oral hygiene procedures.

### **Polishing in the dental office:**

Extensive studies have shown that the abrasive action of Cercon ht on antagonists is less than that of conventional ceramic veneers and not more than standard lithium disilicate ceramics even after finishing and polishing.

Important: After making minor final occlusal adjustments at chairside, we recommend prophylactic polishing of the adjusted sites to a high lustre or adding glaze with a smooth surface prior to temporary or definitive cementing to protect the antagonists from possible abrasion. This also makes the restoration more amenable to proper oral hygiene procedures.

Tooth surfaces whose shade was obtained by staining may present brighter spots where occlusal adjustments were made.

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DeguDent GmbH