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1. Description

1.1 Field of application

Production of residue-free-combustible cast objects for precision casting technology by means of a DLP process.

1.2 Suitable for Dentsply Sirona Primeprint Solution

Primeprint Cast is suitable for the Primeprint 385 nm DLP printer and the Primeprint PPU finishing unit by Dentsply Sirona.

1.3 Note on processing

- ☒ Phosphate-bonded investment materials can be used to manufacture the mold. A preheating temperature of approx. 800 °C, independent of the final casting temperature, must be maintained for approx. 30–45 min, depending on the size of the mold.
- ☒ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for optimal properties. Therefore, it must be ensured that the exposure unit is in proper condition, and the objects are completely cured (refer to chapter “Processing [2]”).
- ☒ After storage, the material in the bottle should be shaken intensively and homogenized before use.
- ☒ Before inserting the cartridge into the material unit, check the same for integrity.
- ☒ Remove the seal label before inserting the cartridge into the material unit.
- ☒ Processing temperature 23 °C ± 2 °C.

1.4 Safety instructions

- ☒ Only for the specified field of use by trained specialists.
- ☒ Product contains (meth)acrylates and phosphine oxides. The ingredients of Primeprint Cast may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ☒ Avoid direct contact with the liquid material and the components before post-curing. Irritating to eyes and skin (sensitization possible).
- ☒ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ☒ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ☒ In case of contact with skin, wash immediately with plenty of soap and water.
- ☒ The processing instructions and precautionary measures must be observed.
- ☒ Refer to the relevant safety data sheet for hazard and safety information.

1.5 Notes

- ☒ DETAX shall not be held liable for any damage caused by misuse.
- ☒ Product may cause allergic reactions.
- ☒ Read and understand the safety data sheet!

1.6 Storage

Store Primeprint Cast in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.

1.7 Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

2. Processing

2.1 Manufacturing process

Refer to the Primeprint operating instructions (REF 6745173) by Dentsply Sirona for the insertion, setup and correct handling of the Primeprint cartridge.

Data processing and creation of the support structure according to the information of the Dentsply Sirona software.

2.2 Building process

Creation of a print job in compliance with the machine and material parameters.

2.3 Finishing process in the Primeprint PPU

If possible, the finishing must take place immediately after the building process in the Primeprint PPU. To do this, remove the Primeprint box (transport container with build platform) from the Primeprint including build job and transfer it to the Primeprint PPU. In the Primeprint PPU, the pre-cleaning and main cleaning steps are performed in two separate vessels with isopropanol and fresh isopropanol (purity >98%). The build job is then automatically fed to the drying and post-exposure process. These steps are performed in the system according to the material. No user interaction is required.

2.4 Component cleaning

Following the steps that have been automatically completed in the Primeprint PPU, the Primeprint box is removed from the Primeprint PPU. The objects are now removed from the build platform, and the support structures are detached.

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Sirona Dental Systems GmbH
DETAX GmbH