

X-Smart® Apex Locator FactFile



X-Smart® Apex Locator

Precise root canal apex localization and correct working length determination are key factors for successful endodontic treatment. The European Society of Endodontists (ESE) and American Association of Endodontists (AAE) recommend using an apex locator in conjunction with verifying radiographs [1, 2]. Dentists report to use this advanced technology device in their daily practice for more than two decades [3-5]. Modern apex locators are valuable devices not only for primary endodontics, but also serve as a useful adjunct during retreatment and have been shown to be highly successful in the determination of root perforations. These devices produce accurate results in different environments, such as presence of blood, pulp tissue or sodium hypochlorite in the root canal space regardless of the file size, the alloy or the file type used [6-9].





Figure 1: X-Smart® Apex Locator; A: Front-side view; B: During apex localization

The X-Smart® Apex Locator model (Figure 1) is based on a unique patented technology with more than 20 years of established clinical history, which uses electrical properties of different tissues inside the root canal system and provides correct results in canals with various conditions - dry and wet [10]. X-Smart® Apex Locator is enhanced with the cutting-edge technology for improved ergonomics and user interface. It has a modern design including cable management. The screen with capacitive touch panel presents larger viewing area with advanced graphical image of the tooth and root canal. The high-resolution display provides all viewing angles [11]. Visual and audio feedback guides the dentist during root canal measurements. X-Smart® Apex Locator is equipped with modern lithium-ion rechargeable battery, which can be charged 4 times faster and lasts longer compared to nickel-metal hydride battery.

“Dr’s Choice”  is a convenient optional function. This feature allows the user to mark a certain predefined position inside the root canal and to obtain additional visual and audio indication when the endodontic file tip reaches this position (Figure 2).




Figure 2: X-Smart® Apex Locator Dr's Choice.

***In vitro* [12] and *ex vivo* [13] data**

Accuracy, repeatability and stability of X-Smart® Apex Locator was tested on plastic blocks (root canal) filled with 0.9% Sodium Chloride (NaCl) and compared to the measurement of a gold-standard apex locator (Dentaport ZX, Morita Co, Japan) using K-files ISO size 10 and 40. The measurements of X-Smart® and the gold-standard apex locator were equivalent, demonstrating precise performance of the X-Smart® Apex Locator.

Accuracy and efficiency to detect the apical foramen position of X-Smart® Apex Locator was compared to the established competitor apex locators in 10 human extracted teeth [13]. Root canals were filled with 0.9% NaCl, different K-files of ISO size 8, 10, 15, 20 and 25 were used. X-Smart® Apex Locator presented 100% identical readings when compared to the reference devices [13].

Conclusion

X-Smart® Apex Locator could be your apex locator of choice, which has gold-standard equivalent performance, while providing improved ergonomics and user interface [12, 13].

References:

1. European Society of Endodontology, Quality guidelines for endodontic treatment: consensus report of the European Society of Endodontology. *Int Endod J*, 2006. 39(12): p. 921-30.
2. American Association of Endodontists. Treatment Standards: Executive Summary. 2019 3 Jan 2023].
3. Christensen, H.Q., et al., Influence of dental education on adoption and integration of technological aids in the delivery of endodontic care by dental practitioners: a survey. *Acta Odontol Scand*, 2022. 80(8): p. 611-618.
4. Savani, G.M., et al., Current trends in endodontic treatment by general dental practitioners: report of a United States national survey. *J Endod*, 2014. 40(5): p. 618-24.
5. Chandler, N.P. and S. Koshy, Radiographic practices of dentists undertaking endodontics in New Zealand. *Dentomaxillofac Radiol*, 2002. 31(5): p. 317-21.
6. Tufenkci, P. and A. Kalayci, Evaluation of the accuracy of different apex locators in determining the working length during root canal retreatment. *J Dent Res Dent Clin Dent Prospects*, 2020. 14(2): p. 125-129.
7. Koç, S., A. Kuş tarıcı, and K. Er, Accuracy of different electronic apex locators in determination of minimum Root perforation diameter. *Aust Endod J*, 2022.
8. Çınar, F. and Y. Üstün, Ex Vivo Evaluation of the Accuracy of 3 Electronic Apex Locators in Different Environments: A Micro-Computed Tomography Study. *Eur Endod J*, 2020. 5(3): p. 226-230.
9. Ferreira, I., A.C. Braga, and I. Pina-Vaz, The Precision of Propex Pixi with Different Instruments and Coronal Preflaring Procedures. *Eur Endod J*, 2019. 4(2): p. 75-79.
10. Internal data: Clinical Evaluation Report for ALC 2021-2022 Rev. 01.
11. Internal data: Product Claims and Comparison.
12. Internal data: In-Vitro Measurements Test Report 5-ORX-005.TR.
13. Internal data: ORYX Ex-Vivo performance verification 5-ORX-000.CLI.