MAILLEFER

ProTaper Next®

Internet

Flexible performance



THE ProTaper PHILOSOPHY REMAINS

A SINGLE FILE SEQUENCE FOR MOST CLINICAL CASES

VARIABLE TAPERS FOR OPTIMIZED CROWN-DOWN TECHNIQUE

APPROVED APICAL FINISHING DIAMETERS

The next generation of the endo reference

ProTaper Next[®] is the successor to the ProTaper[®] Universal system, which has been the gold standard in endodontics for many years. It is an efficient^{*} solution for endo practitioners seeking a versatile, flexible system that will handle the vast majority of root canal treatments.





* Elnaghy, A. M., & Elsaka, S. E. (2014). Assessment of the mechanical properties of ProTaper Next nickel-titanium rotary files. Journal of endodontics, 40(11), 1830-1834.

The next generation...

...covers more difficult clinical cases*

The unique swaggering movement and the greater flexibility of the ProTaper Next® files make it possible to shape more severely curved narrow canals than what was possible before with most NiTi systems on the market.

...brings improved safety*


The risk of file breakage is significantly decreased, at the same time, respect of the original root canal anatomy is greatly improved. Instrument performance and patient safety are- essential to good endodontics.

...shortens the shaping time*

A shorter clinical sequence means that less time is spent changing instruments. The high cutting efficiency also reduces the shaping time. For the dentist, this is valuable time that can be used for other procedures, such as irrigation.

...improves debris removal

The unique off-centred rectangular cross section gives the file a snake-like swaggering movement and provides enlarged space for the debris.



0 1 2 3 4 5 available volume for debris (mm³)

Source: Internal data

Next generation technology

Swaggering Effect

ProTaper Next®'s innovative off-centred rectangular cross section gives the file a snake-like "swaggering" movement as it moves through the root canal. The off-centred cross section and the unique design of the file generate enlarged space for debris hauling.

M-WIRE[®]

The M-Wire^{*} NiTi material improves file flexibility, while still retaining cutting efficiency. This is achieved through an advanced Dentsply Sirona proprietary thermal treatment process.





ProGlider* ProTaper Next®

Motor Settings: 300 rpm / 2 - 5.2 Ncm

ProTaper Next[®] XA (Accessory)

Recommended for preflaring the orifice, removing triangles of dentin, relocating the coronal aspect of a canal away from an external root concavity, and creating more shape as desired.



The next clinical outcome



PRE-OP

Image of maxillary left first premolar. Nonvital pulp had resulted in radiographic lesion of endodontic origin.

POST-OP



Obturation downpack image demonstrating two canals and five visible portals of exit sealed. ProTaper Next* provides predictable shaping preparations for confident and consistent obturation regardless of obturation technique.

Courtesy of Dr. John West





Entering the world of increasing medical standards

Dentsply Sirona Endodontics advocates medical standard care by packaging ProTaper Next® files in pre-sterilised blister packs for single patient use. Using two to three ProTaper Next[®] shaping files for each clinical procedure, instead of four to five previously, makes it affordable for the practitioner to implement a single use policy for endodontic instruments with numerous advantages, namely:



There is no need for reprocessing and re-organising ProTaper Next® NiTi files.



Endodontic files are subject to cyclic fatigue, which is the main cause of file breakage during endodontic treatment. This increases dramatically with multiple use. Single use therefore reduces the risk of file breakage and thus increases patient safety.



Dr. Peet van der Vyver



Dr. Cliff Ruddle

Solution of the second second

Using endodontic files for multiple cases actually reduces their cutting efficiency, affecting the performance and satisfaction for both dentist and patient. With ProTaper Next®, the cutting efficiency stays the same as a new instrument is used for each patient.



No matter how well an instrument is reprocessed, research shows that certain bacteria and tissue remnants may remain on the instrument. Additional uses of instruments* therefore increase the potential of cross-contamination to every subsequent patient.



Ready to go with a complete solution

ProTaper Next[®] is a complete solution, covering all the steps of the endodontic procedure:



Glide path file

ProGlider' is suitable for most root canals, including the highly curved ones* whilst respecting the root canal anatomy. It's therefore an excellent choice combined with ProTaper Next".**

Shaping files

The ProTaper Next^{*} shaping files make it possible to treat the vast majority of endodontic cases.

Endo motor

The X-Smart IQ[™] is a continuous and reciprocating motion cordless motor controlled by an Apple iOS^{*} application developed by Dentsply Sirona Endodontics.

It is a complete solution designed to:

- Support you throughout the endodontic treatment
- Prepare your practice for the future



ProTaper Next[®] Conform Fit[®] gutta-percha

have a greater heat extension and an increased accuracy to match precisely the shape created by the ProTaper Next^{*} files allowing you to achieve a superior apical fit and an accurate tug back response that you can feel.

GuttaCore[®] for ProTaper Next[®] obturators

feature a crosslinked gutta-percha core that sends the surrounding warm, flowable gutta-percha three-dimensionally to follow the natural canal anatomy allowing you to obturate even challenging canals with confidence.

ThermaPrep[®] 2 Oven

is for use with the new GuttaCore[®] for ProTaper Next[®]. In addition to the GuttaCore[®] obturators, the ThermaPrep[®] 2 Oven works with Thermafil[®] and ProTaper Universal obturators.



In the event of retrea

In the event of retreatment, ProTaper[®] Universal retreatment files are the recommended instruments for the obturation unfilling procedure.

* Alovisi M et al. 2016 et Kirchhoff AL et al., 2015

** Pasqualini D et al., 2015, Alovisi M et al., 2016



MAILLEFER

and the second second

ProTaper Next®

A complete solution

LITI	11111	NiTi files				
F	- 931994					
ProTaper Next®	X1		N haaroor	ProTaper Next®	assortment X1 - X	(3 •
A080022100103	File Ster. 21mm	017 .04	6x	A08002219A003	Ass. File Ster.21	
A080022500103	File Ster. 25mm	017 .04	6x	A08002259A003	Ass. File Ster.25	
A080023100103	File Ster. 31mm	017 .04	6x	A08002319A003	Ass. File Ster.31	
A080322100103	File Ster. 21mm	017 .04	Зx	A08032219A003	Ass. File Ster.21	
A080322500103	File Ster. 25mm	017 .04	3x	A08032259A003	Ass. File Ster.25	
A080323100103	File Ster. 31mm	017 .04	3x	A08032319A003	Ass. File Ster.31	
ProTaper Next®	X2					
A080022100203	File Ster. 21mm	025 .06	6x	ProTaper Next®		
A080022500203	File Ster, 25mm	025.06	6x	A080321900A03	File Ster. 19mm	
A080023100203	File Ster. 31mm	025 .06	6x	A080021900A03	File Ster. 19mm	
A080322100203	File Ster. 21mm	025 .06	Зx			
A080322500203	File Ster. 25mm	025 .06	3x			
A080323100203	File Ster. 31mm	025.06	3x	ProGlider ®		
				A092622100103	File Ster. 21mm	
Due Teis en Merut®	V7		The second second	A092622500103	File Ster. 25mm	
ProTaper Next®	X3	-		A092623100103	File Ster. 31mm	
A080022100303	File Ster. 21mm	030 .07	6×	A092322100103	File Ster. 21mm	
A080022500303	File Ster. 25mm	030 .07	6×	A092322500103	File Ster. 25mm	
A080023100303	File Ster. 31mm	030 .07	6x	A092323100103	File Ster. 31mm	
A080322100303	File Ster. 21mm	030 .07	3×			
A080322500303	File Ster. 25mm	030 .07	3×			
A080323100303	File Ster. 31mm	030 .07	3×			
ProTaper Next®	XA		Statistics.	ProTaper® Univer retreatment D1 -		
-		040.06	Ex	A141021600103	File Ster.D1	
A080022100403 A080022500403	File Ster. 21mm File Ster. 25mm	040 .06 040 .06	6x	A141021800103	File Ster.D2	
A080022500403 A080023100403	File Ster. 25mm File Ster. 31mm	040.06	6x 6x	A141022200103	File Ster.D3	
A080023100403 A080322100403	File Ster. 21mm	040.06	3x	A141220090003	Ass. File Ster.	
A080322100403 A080322500403	File Ster. 21mm File Ster. 25mm	040.06	3x 3x			
A080322300403 A080323100403	File Ster. 31mm	040.06	3x 3x			
A080323100403	File Ster. Simili	040.00	37			
ProTaper Next®	X5	N	and the second second	NiTi Manual ada	Assorted Handle Ad	dapto
A080022100503	File Ster. 21mm	050 .06	6x			
A080022500503	File Ster. 25mm	050 .06	6x			
A080023100503	File Ster. 31mm	050 .06	6x			
A080322100503	File Ster. 21mm	050 .06	3×			
	File Ster. 25mm	050.06	Зx			



Endodontic Handpiece Set

Kit contains: X-Smart IQ™ Motor Handpiece X-Smart IQ[™] Contra-angle

			A
ProTaper Next®	assortment X1 - X3		
A08002219A003	Ass. File Ster.21	X1 - X3	6x
A08002259A003	Ass. File Ster.25	X1 - X3	6x
A08002319A003	Ass. File Ster.31	X1 - X3	6x
A08032219A003	Ass. File Ster.21	X1 - X3	3x
A08032259A003	Ass. File Ster.25	X1 - X3	3x
A08032319A003	Ass. File Ster.31	X1 - X3	Зx

ProTaper Next® XA					
A080321900A03	File Ster. 19mm	019 .035	Зx		
A080021900A03	File Ster. 19mm	019 .035	6x		

ProGlider®			
A092622100103	File Ster. 21mm	016 .02	6x
A092622500103	File Ster. 25mm	016 .02	6x
A092623100103	File Ster. 31mm	016 .02	6×
A092322100103	File Ster. 21mm	016 .02	3×
A092322500103	File Ster. 25mm	016 .02	3x
A092323100103	File Ster. 31mm	016 .02	3×

ProTaper® Univer retreatment D1			
A141021600103	File Ster.D1	030.09	6x
A141021800103	File Ster.D2	025 .08	6x
A141022200103	File Ster.D3	020 .07	6x
A141220090003	Ass. File Ster.	D1- D3	6x



ter Kit



A105420000000

Kit contains: X-Smart IQ™ Motor Handpiece X-Smart IQ™ Contra-angle Endo IQ™ iPad Mini™ Protective Case 10 blisters ProTaper Next* Assortment X1-X3, 25mm (30 files) 1 blister ProGlider[®], 25mm (3 files)



Dentsply Sirpno

-

ProTaper Next*



ProTaper Next* Conform Fit* gutta-percha

8- -	Pro B0
8	Pro B0
8- -	Pro B0
@_ _	Pro

Obturators Blister

unispiry uning	1	1
Contraction of the local division of the loc	1111	2
and the second	11111	





Maillefer Instruments Holding Sàrl ******* Chemin du Verger 3 - 1338 Ballaigues Switzerland **CE** 2797 Phone: +41 (0)21 843 92 92 endo@dentsplysirona.com dentsplysirona.com files, gutta-percha points





Oven

Thermaprep®2 Oven for GuttaCore® A 0176 2300 0200 A 0176 2300 0233 A 0176 1150 0200

EU (230V) UK (230V) (115V)

ProTaper Next® Paper Point X2 A 022Q 00000 200

ProTaper Next® Paper Point X3 A 022Q 00000 300

ProTaper Next® Paper Point X2 & X3 A 022Q 00090 000

ProTaper Next® Paper Point X4 & X5 A 022Q 00090 100

180x sterile paper points

180x sterile paper points

90x sterile paper points X2 90x sterile paper points X3

90x sterile paper points X4 90x sterile paper points X5

oTaper Next[®] Conform Fit[®] gutta-percha X2 DOPNGPF000X2 60x gutta percha points oTaper Next® Conform Fit® gutta-percha X3

DOPNGPF000X3 60x gutta percha points oTaper Next® Conform Fit® gutta-percha X2 & X3

OOPNGPF00X23 36x gutta percha points X2 24x gutta percha points X3

roTaper Next® Conform Fit® gutta-percha X4 & X5 BOOPNGPF00X45 36x gutta percha points X4 24x gutta percha points X5

utta Core® Blister for PTN X2 1803B 004 25 00 6×	
utta Core® Blister for PTN X3 1803B 004 30 00 6×	r L
utta Core® Blister for PTN X4 1803B 004 40 00 6x	Į.
utta Core® Blister for PTN X5 1803B 004 50 00 6x	

i0x (5 blisters)
i0x (5 blisters)
iOx (5 blisters)



GuttaCore", paper points



Maillefer Instruments Holding Sàrl Chemin du Verger 3 - 1338 Ballaigues Switzerland Phone: +41 (0)21 843 92 92

Phone: +41 (0)21 843 92 92 endo@dentsplysirona.com dentsplysirona.com



DENTSPLY Tulsa Dental Specialties DENTSPLY International, Inc. 608 Rolling Hills Drive Johnson City, TN 37604, USA Phone: 1-800-662-1202 dentsplysirona.com

GuttaCore", paper points

