

Diversey™ Oxivir® 1 Disinfectant Cleaners

Diversey™ Oxivir® 1 Disinfectant Cleaners in ready-to-use (RTU) liquid and wipes are formulated with the patented AHP® (Accelerated Hydrogen Peroxide™) technology which contains a high-performance surfactant system that delivers superior cleaning results versus many other disinfectant technologies. Oxivir is used by many leading hospitals and health systems across the US and Canada. It is chosen for its speed, efficacy, and safety profile.

Features & Benefits

- Meets OSHA Bloodborne Pathogen Standards for decontaminating blood and body fluid spills.
- In all six EPA toxicity categories, Oxivir 1 Disinfectant Cleaners, fall into Category IV, the lowest level of hazard (practically non-toxic, not an irritant) and require no safety warnings.
- Diversey Oxivir 1 Disinfectant Cleaners are on the EPA List N: Disinfectants for Coronavirus (COVID-19). EPA expects all products on List N to kill the coronavirus SARS-CoV-2 (COVID-19) when used according to the label directions.¹
- Broad spectrum, multi-surface cleaner/disinfectant and sanitizer. Harnesses the cleaning power of hydrogen peroxide. Specially designed to penetrate the toughest soils.
- Fast, Effective, Responsible and Sustainable formulation. Environmentally responsible with no Volatile Organic Compounds (VOC's), no added fragrance, no Nonylphenol Ethoxylates (NPE's). Fast-acting, convenient formula, suitable for daily use, and contact time of ONE minute or less.
- Formulated with patented AHP® technology containing no phosphates, quats or phenols. The active ingredient breaks down to water and oxygen after use.



100850923

160 count wipes
6 in. x 7 in.
15.24 cm x 17.78 cm
Canister

Diversey™ Oxivir® 1 Wipes Ready-to-Use Hospital Disinfectant Cleaner Wipe

Next generation ready-to-use, pre-wetted, one-step disinfectant cleaner wipe, based on proprietary AHP® technology. Oxivir® 1 Wipes are formulated to deliver broad spectrum disinfection in ONE minute or less, with 87 total claims.

- Short contact times of 1 minute or less Virucidal, 1 minute Bactericidal, Tuberculocidal and Fungicidal and 10-second for non-food contact sanitizing.

9MSC005-12
Spray
Nozzle



100850916
32 ozs.
Bottle

Diversey™ Oxivir® 1 RTU (ready-to-use) General Virucide, Bactericide, Tuberculocide, Fungicide, Sanitizer

Next generation ready-to-use disinfectant cleaner based on proprietary AHP® technology. Oxivir® 1 is formulated to deliver broad spectrum disinfection in ONE minute or less, with 81 total claims.

- Short contact times of 30-second Virucidal (1 minute for select animal viruses), 1 minute Bactericidal, Tuberculocidal and Fungicidal, 10-second Non-Food Contact Sanitizing and 30-second Soft Surface Sanitizing.

PRODUCT NAME	Oxivir® 1 Wipes	Oxivir® 1 RTU
Company	Diversey	Diversey
Active Ingredients	.5% Hydrogen Peroxide	.5% Hydrogen Peroxide
Signal Word (Toxicity Concern)	NONE	NONE
Required Surface Contact Time	1 minute	1 minute
NON-ENVELOPED VIRUSES		
Adenovirus Type 2/6	30 seconds	30 seconds
Adenovirus Type 7	30 seconds	30 seconds
Adenovirus Type 8	30 seconds	30 seconds
Enterovirus Type 68 (VR-561)	1 minute	1 minute
Norovirus	30 seconds	30 seconds
Poliovirus Type 1	30 seconds	30 seconds
Rhinovirus (Common Cold)	30 seconds	30 seconds
Rotavirus	30 seconds	30 seconds
ENVELOPED VIRUSES		
Avian Influenza A Virus (H5N1)	30 seconds	30 seconds
Cytomegalovirus	30 seconds	30 seconds
Influenza A, Hong Kong Strain	30 seconds	30 seconds
Influenza B Virus	30 seconds	30 seconds
Hepatitis B Virus (HBV)	30 seconds	30 seconds
Hepatitis C Virus (HCV)	30 seconds	30 seconds
Herpes Simplex Virus Type 1	30 seconds	30 seconds
Herpes Simplex Virus Type 2	30 seconds	30 seconds
HIV-1	30 seconds	30 seconds
Human Coronavirus (SARS virus type)	30 seconds	30 seconds
Measles Virus	30 seconds	30 seconds
MERS-CoV	30 seconds	30 seconds
Vaccinia Virus	30 seconds	30 seconds
Pandemic Influenza A Virus	30 seconds	30 seconds
Parainfluenza Virus Type 3	30 seconds	30 seconds
Respiratory Syncytial Virus (RSV)	30 seconds	30 seconds
GRAM NEGATIVE BACTERIA		
Acinetobacter baumannii, Multi-Drug Resistant (MDR)	1 minute	1 minute
Bordetella bronchiseptica	1 minute	1 minute
Bordetella pertussis	1 minute	1 minute
Burkholderia cepacia	1 minute	1 minute
Campylobacter jejuni	1 minute	1 minute
Citrobacter freundii	1 minute	1 minute
Enterobacter aerogenes	1 minute	1 minute
Enterobacter cloacae, New Delhi metallo β lactamase - 1 (NDM 1)	1 minute	1 minute
Escherichia coli	1 minute	1 minute
Escherichia coli, Carbapenem resistant (CRE)	1 minute	1 minute
Escherichia coli, ESBL producing	1 minute	1 minute
Escherichia coli, New Delhi metallo β lactamase - 1 (NDM 1)	1 minute	1 minute
Haemophilus influenzae	1 minute	1 minute
Escherichia coli O157:H7	1 minute	1 minute
Klebsiella aeogenes (MDR)	1 minute	1 minute
Klebsiella oxytoca	1 minute	1 minute
Klebsiella oxytoca, Multi-Drug Resistant Metallo-beta Lactamase (MBL) Positive	1 minute	1 minute
Klebsiella pneumoniae	1 minute	1 minute
Klebsiella pneumoniae, Carbapenem Resistant (CRE)	1 minute	1 minute
Klebsiella pneumoniae, ESBL producing	1 minute	1 minute
Klebsiella pneumoniae, Multi-Drug Resistant	1 minute	1 minute
Klebsiella pneumoniae, New Delhi metallo-beta-lactamase (NDM-1)	1 minute	1 minute
Legionella pneumophila	1 minute	1 minute
Pasteurella multocida	1 minute	1 minute
Pseudomonas aeruginosa	1 minute	1 minute
Pseudomonas aeruginosa, Multi-Drug Resistant Metallo-beta Lactamase (MBL) Positive	1 minute	1 minute
Salmonella enterica (choleraesuis)	1 minute	1 minute
Serratia marcescens	1 minute	1 minute
Shigella dysenteriae	1 minute	1 minute
Stenotrophomonas (Pseudomonas maltophilia)	1 minute	1 minute
Vibrio cholerae	1 minute	1 minute
GRAM POSITIVE BACTERIA		
Enterococcus faecalis, Vancomycin-resistant (VRE)	1 minute	1 minute
Enterococcus faecium	1 minute	1 minute
Enterococcus faecium,MDR	1 minute	1 minute
Listeria monocytogenes	1 minute	1 minute
Staphylococcus aureus, CA-MRSA	1 minute	1 minute
Staphylococcus aureus	1 minute	1 minute
Staphylococcus aureus, (ATCC 33592), MRSA	1 minute	1 minute
Staphylococcus aureus, (ATCC 33591), MRSA	1 minute	1 minute
Staphylococcus aureus, Vancomycin-intermediate resistant (VISA)	1 minute	1 minute
Staphylococcus aureus, Vancomycin-resistant (VRSA)	1 minute	1 minute
Staphylococcus epidermidis, Methicillin-resistant (MRSE)	1 minute	1 minute
Streptococcus pneumoniae	1 minute	1 minute
Streptococcus pneumoniae, Penicillin-resistant (PRSP)	1 minute	1 minute
Staphylococcus pseudintermedius	1 minute	1 minute
Streptococcus pyogenes	1 minute	1 minute
TUBERCULOCIDAL		
Mycobacterium bovis (TB)	1 minute	1 minute
FUNGICIDAL		
Candida albicans	1 minute	1 minute
Candida auris	1 minute	1 minute
Candida glabrata	1 minute	1 minute
Microsporium canis	1 minute	1 minute
Trichophyton mentagrophytes	30 seconds	30 seconds
NON-FOOD CONTACT		
Acinetobacter baumannii	10 seconds	10 seconds
Campylobacter jejuni	10 seconds	10 seconds
Enterobacter aerogenes	10 seconds	10 seconds
Enterococcus faecalis, Vancomycin-resistant (VRE)	10 seconds	10 seconds
Escherichia coli O157:H7	10 seconds	10 seconds
Listeria monocytogenes	10 seconds	10 seconds
Klebsiella pneumoniae	10 seconds	10 seconds
Pseudomonas aeruginosa	10 seconds	10 seconds
Salmonella enterica (choleraesuis)	10 seconds	10 seconds
Staphylococcus aureus	10 seconds	10 seconds
Staphylococcus aureus, Methicillin-resistant (MRSA)	10 seconds	10 seconds
SOFT SURFACE		
Enterobacter aerogenes		30 seconds
Enterococcus faecalis, Vancomycin-resistant (VRE)		30 seconds
Staphylococcus aureus		30 seconds
Staphylococcus aureus, Methicillin-resistant (MRSA)		30 seconds