# Primeprint Guide for Primeprint Solution

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## 1. Description

1.1 INTENDED USE
1.2 INDICATION
1.3 CONTRAINDICATION
1.4 ADVERSE EFFECTS
1.5 PATIENT TARGET GROUP
1.6 INTENDED USERS
1.7 SUITABLE FOR DENTSPLY-SIRONA-PRIMEPRINT-SOLUTION
1.8 SPECIAL MANUFACTURING REQUIREMENTS
1.9 STERILIZATION OF DRILLING TEMPLATES
1.10 TECHNICAL DATA FOR THE STERILISATION PROCESS
1.11 SAFETY INFORMATION
1.12 STORAGE
1.13 MANUFACTURING PROCESS
1.14 FINISHING
1.15 SURFACE PROCESSING

**Processing:**  
at 23 °C ± 2 °C

**Storage:**

☐ Consult Instructions for Use

☐ Catalogue Number

☐ Batch Code

☐ Use by Date

Rx only ☐ Caution: Federal U.S. law restricts this device to sale by or on the order of a dentist (or trained specialist personnel).

☐ Temperature limitation

☐ Keep away from sunlight

☐ Manufacturer

☐ Medical Device

☐ Warning

☐ Warning

**1.8.2 DLP printer and its operation software**

Printer Manufacture	SDS Primeprint
Light source	385 nm
Pixel size x, y	±70 µm
Light intensity	6,2 mW/cm <sup>2</sup>
Operation Software	inLab CAM SW 22

## 1.8.3 Printing Parameters

Parameter data set*	Primeprint Guide
Curing time	2 s
Height base plate	0,5 mm
Support point size	0,7 mm (tip)
Max. support size	1,5 mm
Support density	65 %
Connection between supports	diagonal struts
Slicing	100 µm
Build speed	45 mm/hr
Build path	surface layer drawing path
Recommended orientation angle (degree)	Printing mode optimized quality

\*The set of parameter include all relevant material, printer and post-processing specific information

## 1.8.4 Note on processing

- ☒ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is very important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the chapter "Manufacturing process [1.11]").
- ☒ After storage, intensively shake and homogenize the material in the cartridge before use.
- ☒ Before inserting the cartridge into the material unit, check the same for integrity.
- ☒ Remove the seal label before inserting the cartridge into the material unit.
- ☒ Minimum material thickness for the design is 1.5 mm.
- ☒ Maximum curing depth\* at direct post-exposure: 6 mm \*In case of large objects and exposure on both sides, the material thickness can be up to 12 mm.
- ☒ Polish the surface mechanically if necessary.

## 1.8.5 Environmental Conditions

Processing temperature 23 °C ± 2 °C.

## 1.8.6 Cleaning Kit

Rinse bath tubs, flush cutter, paper towel, squeeze bottle for isopropyl alcohol, scraper.

## 1.8.7 Curing light equipment

SDS Primeprint PPU.

## 1.9 Sterilization of drilling Templates

- ☒ Locally applicable regulations and hygiene guidelines are to be observed.
- ☒ Primeprint Guide is validated for the „W&H Lisa 522, Program UNIVERSAL 121“ sterilization process (see „technical data on sterilization process [1.10]“).
- ☒ Steam sterilization may only be performed with equipment that complies with the standards EN 13060 and EN 285.
- ☒ Responsibility for sterility lies with the user.

## 1.10 Technical data for the sterilisation process

- ☒ A fractionated pre-vacuum with a plateau phase at 121 °C/ 15 minutes is to be selected, whereby this refers to the requirements for a B-cycle small steriliser in terms of its venting performance.

## 1.11 Safety Information

- ☒ Only for the specified intended use by trained specialists.
- ☒ Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- ☒ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ☒ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ☒ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ☒ After contact with skin wash immediately with water and soap.
- ☒ Biocompatibility is only guaranteed with complete polymerisation.
- ☒ MMA-free.
- ☒ Refer to the relevant safety data sheet for hazard and safety information.

## 1.12 Storage

Primeprint Splint is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.

## 1.13 Manufacturing process

Refer to the Primeprint operating instructions (REF 6745173) by Dentsply Sirona for the insertion, setup and correct handling of the Primeprint cartridge. Data processing and creation of the support structure according to the information of the Dentsply Sirona software.

## Mixing

Before inserting the cartridge into the material unit the material must be shaken. Inadequate mixing could cause deviations of colour and print failures.

## Prepare the printer

Inserting the material unit and the Primeprint-Box into the printer. Different batches of the same 3D printing resin may not be mixed. The 3D printing resin can remain in the vat for a maximum of five month.

## Construction process

Visual inspection of print resin material: in case print resin material is visually contaminated, it may not be used and must be discarded.

Generate a print job complying with printer and material parameters.

For setting up the printer follow the instructions for use of the printer.

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer (inLab CAM).

## Post-processing

If possible, the finishing must take place immediately after the building process in the Primeprint PPU. To do this, remove the Primeprint box (transport container with build platform) from the Primeprint including build job and transfer it to the Primeprint PPU. In the Primeprint PPU, the pre-cleaning and main cleaning steps are performed in two separate vessels with isopropanol and fresh isopropanol (purity >98%). The build job is then automatically fed to the drying and post-exposure process. These steps are performed in the system according to the material. No user interaction is required.

**CAUTION:** The use of Primeprint is only approved when applied with the compatible devices mentioned in this instructions for use. Any unauthorized changes to the process equipment, parameters, or software may result in a device that is out of specification.

## 1.14 Finishing

Following the steps that have been automatically completed in the Primeprint PPU, the Primeprint box is removed from the Primeprint PPU. The objects are now removed from the build platform, and the support structures are detached.

## 1.15 Surface Processing

Polish the surface mechanically if necessary.

## 1.16 Cleaning of completed objects in daily use

### By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

### By the patient

For daily cleaning, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards.

### Cleaning

tablets can be used following the manufacturer's instructions.

### Note

- ☒ Do not use toothpaste; the abrasive particles are coarse and scratch the surface.
- ☒ Do not use any acids or solvents.

## 1.17 Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

## 1.18 Notes

- ☒ Always keep container tightly sealed, immediately close the container carefully after each use.
- ☒ Please follow the instructions on the safety data sheet!

### For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to [incident@detax.de](mailto:incident@detax.de) and to the competent authority in which the user and/or patient is established.

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Sirona Dental Systems GmbH
DETAX GmbH