

Care & Cleaning Products

Surface Cleaning and Disinfection

Neutral Detergent (Approved for Use on all Dentsply Sirona Equipment)
Alpro IC-100



B4658 | Alpro IC-100 Spray Bottle 500 mL

Neutral Detergent
For cleaning surfaces and medical devices including upholstery
Spray on instruments and allow the solution to take effect for 5 minutes, then rinse and sterilise

B4659 | Alpro IC-100 Bottle 5 L

Neutral Detergent
Can be used in conjunction with a cloth or pour 1.5 L in to the MAXIWIPES L to create wet wipes

CleanWipes



4153 | CleanWipes Canister with 70 wet wipes

Surface Neutral Detergent
Wet Wipes
Can be used on most surfaces including handpieces, instruments, intra-oral cameras and the treatment centre including upholstery
Solution is Alpro IC-100

4154 | CleanWipes Refill bags - 70 wet wipes per bag

Surface Neutral Detergent
Wet Wipes - Refills
For refilling Clean Wipes canister
Solution is Alpro IC-100

Disinfectants (Approved for use on all Dentsply Sirona Equipment)
PlastiSept eco



4645 | Alpro Plastisept Eco 500mls

Surface Disinfectant Solution
Can be used in conjunction with a cloth or poured into MAXIWIPES to create wet wipes
Alternative is Durr FD366

3063E | Alpro Plastisept Eco 5 L

Surface Disinfectant Solution
Can be used in conjunction with a cloth or poured into MAXIWIPES to create wet wipes
Alternative is Durr FD366

4625 | PlastiSept eco - 2x Canisters of 70 wipes

Surface Disinfection
Wet Wipes
Can be used on most surfaces including handpieces, instruments, intra-oral cameras and the treatment centre including upholstery
Alternative is FD366

4622 | PlastiSept eco - 70 wipe refill bag

Surface Disinfection
Wet Wipes - Refills
For refilling PlastiSept eco canister
Alternative is FD366

Durr FD 366 Sensitive



DURRFD366 | Durr FD 366 Sensitive Bottle 2.5 L

Surface Disinfection solution
Can be used in conjunction with a cloth or pour 2.5 L in to the DURRFDMULTIWIPES to create wet wipes
Alternative is PlastiSept Eco Wipes



DURR366WIPESTUB | FD 366 Wet Wipes Container and 100 wet wipes

Wet Wipes
Presoaked with FD366 solution
Alternative is PlastiSept Eco

Durr FD 366



DURRFD366WIPES | Durr FD 366 4 x 100 Refill packets of Wet Wipes

Surface Disinfection

Wet Wipes - Refills

For refilling the DURRFDMULTIWIPES Container and

DURR366WIPESTUB

Alternative is PlastiSept Eco Wipe

Upholstery Cleaner

Durr FD 360 & Sponges



DURRFD360 | Durr FD 360 500 mL of Solution and 10 x sponges

Upholstery Cleaner

Spray FD 360 on to sponge to clean upholstery on the stools and chairs only.

Don't use sponge on lounge upholstery.

Alternative is PlastiSept eco Wipes or IC-100

Spray Bottle

Durr Spray Bottle



DURRSPRAYBOTTLE | Durr Spray Bottle

Spray Bottle Empty

Can be filled with any approved Cleaning Solution

Dry Wipes and Containers

Dry wipes contain no cleaning or disinfecting agents, these agents can be purchased separately and added to a container to create convenient to use wipes for any purpose. Cleaning and disinfecting agents only need to gently soak the wipe to be effective.

Durr Multi Wipes



DURRFDMULTIWIPES | Durr Multi Wipes Container and 180 Dry Wipes

Dry wipes

Add any disinfectant or neutral detergent solution to create wet wipes

Alternative is MaxiWipes L B4616

Alpro MaxiWipes L



4612 | MaxiWipes L Container + 250 Dry Wipes

Dry Wipes

Add any disinfectant or neutral detergent solution to create wet wipes

Keep container and refill with 4613



4613 | MaxiWipes L 4 x 250 Refill packets of Dry Wipes

Dry Wipes - Refills

Refilling the B4616 Container

Waterline Maintenance and Treatment

Waterline Maintenance Solution

Suitable for use with water bottles and integrated systems for all Dentsply Sirona Treatment Centres

Dentosept S



6610310 | Dentosept S 1 L

Disinfection Solution of Water Lines for internal use in treatment centres. Alternative product Alpron.

Alpron



3183 | Alpron 1 L

Disinfection Solution of Water Lines for internal use in treatment centres. Alternative product Dentosept.

Waterline Testing

Biotest



3079 | Alpro Dip slide Biotest pack of 10x

Dip slide for semi quantitative monitoring of microbial load of water and surfaces. The microbial count is measured in Colony-Forming Units (CFU) per mL.

DS Single Test



5853775 | DS Single Total Count Test

Dip slide for semi quantitative monitoring of microbial load of water and surfaces. The microbial count is measured in Colony-Forming Units (CFU) per mL.

Waterline Treatment

For Treatment Centres experiencing a CFU count of 200 or higher, chairs should be tested and only treated when required.

Bilpron



3181 | Bilpron 1L

For the decontamination and prevention of biofilm formation in Treatment Centre water lines. Minimum contact time for treatment is 12 hours. Can be utilised as a follow up to initial shock treatment in contaminated lines.

BRS Kit



3179 | BRS Treatment Kit

Full set for shock treatment of Treatment Centres experiencing a CFU count of 200 or higher. Includes BRS Remove, precleaner and Bilpron.



B3179WB | BRS Treatment Kit

Full set for shock treatment of Treatment Centres when using with Dentosept S. Includes BRS PreCleaner and Remove.

Water Bottle Cleaner



3199 | BC-San 100 Disinfection Concentrate for Water Bottles

Concentrate for regular disinfection of water bottles.

Suction Line Maintenance

Daily Suction Line Maintenance Solution

Alpro-Jet D



BHY0001 | Alpro-Jet D 5 L

Suction System Disinfectant

Use daily to clean the suction hoses and spittoon bowl.

Use 50 mL to 1 L of water

Can also be placed in the antimicrobial tank - if a tank is fitted in the Treatment Centre

Alternative is Orotol Plus Bottle 2.5 L or 30 L

Durr Orotol Plus



DURROROTOLPLUS | Durr Orotol Plus Bottle 2.5 L

Suction System Disinfectant

Use daily to clean the suction hoses and spittoon bowl. Use 40 mL to 2 L of water

Can also be placed in an antimicrobial tank - if a tank is fitted in the Treatment Centre.

Alternative is Alpro-Jet D



DURROROTOL30LT | Durr Orotol Plus Bottle 30 L

Suction System Disinfectant

Use daily to clean the suction hoses and spittoon bowl. Use 40 mL to 2 L of water

Can also be placed in the antimicrobial tank - if a tank is fitted in the Treatment Centre

Alternative is Alpro-Jet D

Weekly Suction Line Maintenance Solution

Alpro-Jet W



BHY0006 | Alpro-Jet W 1 L

Suction System Cleaner

Use Weekly to clean the suction hoses and spittoon bowl.

Use 50 mL to 1 L of water

Alternative is MD555



BHY0004 | Alpro-Jet W 5 L

Suction System Cleaner

Use Weekly to clean the suction hoses and spittoon bowl.

Use mixture 50 mL to 1 L of water

Alternative is MD555

Durr MD 555



CC555C6155 | Durr MD 555 Bottle 2.5 L

Suction System Cleaner

Use Weekly to clean the suction hoses and spittoon bowl.

Leave in the suction hoses for 30 minutes, then suck up some additional water to flush the hoses

Mixture is 5%

Use 100 mL to 2 L of water

Alternative is Alpro-Jet W

Suction Bottle with Adaptors (for all Suction Solution Types)

Orotol Suction Bottle



DURROROCUP | Orotol Suction Cup

2 L Bottle with suction adaptors
Use with Orotol / MD 555 / Alpro-Jet D / Alpro-Jet W
Attach suction hoses and suck the solution through the hoses
Pour remaining solution in to the spittoon bowl
Use 1 L per Treatment Centre

0780-300-12 | OroCup Insert 6mm

For small suction hoses

0780-300-22 | OroCup Insert 11mm

For standard suction hoses

0780-300-23 | OroCup Insert 16mm

For big suction hoses



CED100A9900 | Pump for Orotol Bottle

For dosing from Orotol Bottle directly

Sterile Filter For Taps

i3 Sterile Filter System

A tap which will produce sterile water from any standard tap line (not a water filtration system only for eliminating bioburden and creating fresh sterile water)



i3ONE | Biodegree i3ONE Tap Filter



i3ONE | M24 Biodegree M24 coupling for i3ONE



i3ONE | M22 Biodegree M22 coupling for i3ONE