


How CEREC Works Sample Mockup



One Hour. One Visit.™
Crowns, inlays, veneers & more

Hard to fit dental visits into your schedule? Here's how we make it easier

Dentistry Powered by CEREC helps us make custom crowns and more in one visit, in about an hour. Here's how it works:

1. We capture an image of your tooth with a digital scanner.
2. A 3D model of your teeth and the new restoration is created on the computer.
3. The restoration is then fabricated from a ceramic block on our milling machine.
4. The dentist places the custom-fit restoration a short time later.

In just one visit using this advanced system, we can prepare your tooth, then make the dental restoration right in our office in about one hour.

That's it. No messy impressions, no temporary crowns, and no long wait for your permanent restoration.

Dentistry Powered by CEREC: Call and ask us about it today.