



Atlantis® CustomBase solution

Why go standard?

 **Dentsply
Sirona**
Implants

Placement of material junction of abutment and crown, considering biological principles

Patient-specific design of core and height for optimized retention and support


Customized emergence profile for enhanced tissue management

30°

Optimal positioning of the screw access channel for improved esthetics, function and installation procedures¹

New!

Angled screw access channel allows the screw access hole to be optimally positioned.¹




Stock titanium-base component provides minimal retention and support for the crown.

Atlantis CustomBase solution designed with emergence profile and patient-specific design of the abutment core and height, for optimized retention and support for the crown.

1. Currently available for Astra Tech Implant System, OsseoSpeed EV only.

Odds are in your favour – beyond titanium bases

Atlantis CustomBase solution for patient-specific, screw-retained single-tooth restorations provides a solid foundation for predictable outcomes. After all, one size does not fit all. The solution is now available with an angled screw access channel, so you can optimally position the screw access channel to help improve esthetics, function and installation procedures – so why not ensure the odds are in your favour?



Atlantis solutions are backed up by a comprehensive warranty. Dentsply Sirona covers the abutment, the crown and the implant.²

51 

51 published clinical studies and case reports documenting success using Atlantis abutments on different implant platforms.



2. Subject to the terms and conditions of the Atlantis abutments warranty.

About Dentsply Sirona Implants

Dentsply Sirona Implants offers comprehensive solutions for all phases of implant therapy, including Ankylos[®], Astra Tech Implant System[®] and Xive[®] implant lines, digital technologies, such as Atlantis[®] patient-specific solutions and Simplant[®] guided surgery, Symbios[®] regenerative solutions, and professional and business development programs, such as STEPPS[™]. Dentsply Sirona Implants creates value for dental professionals and allows for predictable and lasting implant treatment outcomes, resulting in enhanced quality of life for patients.

About Dentsply Sirona

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies, with a 130-year history of innovation and service to the dental industry and patients worldwide. Dentsply Sirona develops, manufactures, and markets a comprehensive solutions offering including dental and oral health products as well as other consumable medical devices under a strong portfolio of world class brands. As The Dental Solutions Company[™], Dentsply Sirona's products provide innovative, high-quality and effective solutions to advance patient care and deliver better, safer and faster dentistry. Dentsply Sirona's global headquarters is located in York, Pennsylvania, and the international headquarters is based in Salzburg, Austria. The company's shares are listed in the United States on NASDAQ under the symbol XRAY.

Visit www.dentsplysirona.com for more information about Dentsply Sirona and its products.

THE DENTAL
SOLUTIONS
COMPANY[™]

