

Your Hygiene Handpiece Just Got Smarter...Introducing the Nupro Freedom™ Cordless Prophylaxis System with SmartMode™ Technology

York/Salzburg, December 1, 2016

Dentsply Sirona Preventive announces the launch of the Nupro Freedom™ Cordless Prophylaxis System with SmartMode™ Technology. SmartMode Technology provides clinicians with an innovative method of controlling the speed and power of the prophylaxis cup, while also providing freedom from cords and foot pedals. In addition, Nupro Freedom provides 25 percent more speed and power than the previous generation Freedom handpiece.

Nupro Freedom Cordless Prophylaxis System Package includes a cordless handpiece and three autoclavable outer sheaths that allow a dental practice to comply with CDC Infection Prevention Guidelines without needing to purchase three separate handpieces. Also, included are Disposa-Shield® disposable polyethylene barriers to facilitate cleaning and disposable prophylaxis angles.

SmartMode Technology uses the gentle pressure of the prophylaxis cup against the tooth to control its speed. Adding pressure increases the speed of the prophylaxis cup while removing pressure decreases the speed. An optional cordless foot pedal is available for those who prefer the foot pedal operation.

“The Nupro Freedom Cordless Prophylaxis System is quieter than traditional air driven hygiene or low speed handpieces, providing patients with a more pleasant polishing experience,” said Gail Malone, Global Professional Services Senior Manager. “Additionally, when combined with Nupro Freedom® Slim Disposable Prophylaxis Angles, they provide less chatter and vibration for the patient as well as better visibility and access for the clinician. It’s a win-win for the patient and the clinician.”

Customers can request a demo and learn more about the Nupro Freedom with SmartMode Technology Prophylaxis System by contacting Dentsply Sirona Customer Service at 800-989-8825 or by visiting nuprofreedom.com.

Press Contact

Nick Harrison
Marketing Communications Manager
Dentsply Sirona Preventive
1301 Smile Way
York, PA 17404
Phone: 717-767-8538
Nick.Harrison@dentsplysirona.com

ABOUT DENTSPLY SIRONA:

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies, with a 130-year history of innovation and service to the dental industry and patients worldwide. Dentsply Sirona develops, manufactures, and markets a comprehensive solutions offering including dental and oral health products as well as other consumable medical devices under a strong portfolio of world class brands. As The Dental Solutions Company™, Dentsply Sirona's products provide innovative, high-quality and effective solutions to advance patient care and deliver better, safer and faster dentistry. Dentsply Sirona's global headquarters is located in York, Pennsylvania, and the international headquarters is based in Salzburg, Austria. The company's shares are listed in the United States on NASDAQ under the symbol XRAY. Visit www.dentsplysirona.com for more information about Dentsply Sirona and its products.

IMAGES

are available for > **Download** on the website.



Fig. 1: Nupro Freedom™ Cordless Prophy Handpiece




Fig. 2: 1. Cup starts spinning, 2. Cup speed increases as you press against the tooth; reduce and increase speed as you wish with less or more pressure, 3. Cup speed decreases automatically as it is pulled away from the tooth