













## HAND VS. ULTRASONIC INSTRUMENTATION

The chart below provides you with an easy snapshot of how DENTSPLY Cavitron® Inserts can replace all of your hand instrumentation needs.

| DENTSPLY CAVITRON INSERT   | HAND SCALER EQUIVALENT  | DENTSPLY CAVITRON INSERT   | HAND SCALER EQUIVALENT  |
|--|---|--|---|
|  <p>FSI-SLI-10S<br/>SLI-10S<br/>P-10</p>        | <ul style="list-style-type: none"> <li>• Columbia 13/14 Universal curette</li> <li>• 4R/4L Universal curette</li> <li>• Younger Good Universal curette</li> <li>• Any Universal curette</li> <li>• 5/6, 7/8 Gracey</li> </ul> |  <p>FSI-10/PowerLINE™ 10<br/>TFI-10</p>                 | <ul style="list-style-type: none"> <li>• 4R/4L Rigid Universal Curette</li> <li>• 13/14/ Rigid Columbia</li> <li>• 5/6, 7/8 Rigid Gracey</li> </ul> |
|  <p>FSI-SLI-10R<br/>SLI-10R</p>               | <ul style="list-style-type: none"> <li>• 11/12 or 15/16 Gracey or mini or After 5 Gracey</li> <li>• 13/14 or 17/18 Gracey or mini or After 5 Gracey</li> </ul>  |  <p>FSI-3/PowerLINE™<br/>TFI-3</p>                    | <ul style="list-style-type: none"> <li>• Beaver Tail</li> <li>• Younger good</li> </ul>   |
|  <p>FSI-SLI-10L<br/>SLI-10L</p>               | <ul style="list-style-type: none"> <li>• 11/12 or 15/16 Gracey or mini or After 5 Gracey</li> <li>• 13/14 or 17/18 Gracey or mini or After 5 Gracey</li> </ul>  |  <p>CAVITRON THINSERT® INSERT</p>                     | <ul style="list-style-type: none"> <li>• Sickles</li> <li>• Nevi's</li> <li>• Jaquettes</li> </ul>  |
|  <p>FSI SLI 1000</p>                          | <ul style="list-style-type: none"> <li>• Younger good</li> <li>• 13/14 Rigid Columbia</li> <li>• 4R/4L Rigid Universal Curette</li> </ul>   |  <p>SOFTIP ULTRASONIC IMPLANT INSERT</p>              | <ul style="list-style-type: none"> <li>• Implant Scalers</li> </ul>   |
|  <p>FSI-1000/POWERLINE™ 1000<br/>TFI-1000</p> | <ul style="list-style-type: none"> <li>• 4R/4L Rigid Universal Curette</li> <li>• 13/14 Rigid Columbia</li> <li>• 5/6, 7/8 Rigid Gracey</li> </ul>  |  <p>DIAMOND COAT®</p>                                 | <ul style="list-style-type: none"> <li>• Diamond Coated hoes, files and chisels</li> </ul>  |
|  <p>FSI-100/PowerLINE™ 100</p>                | <ul style="list-style-type: none"> <li>• 4R/4L Rigid Universal Curette</li> <li>• 13/14/ Rigid Columbia</li> <li>• 5/6, 7/8 Rigid Gracey</li> </ul>   |  <p>AIR POLISHING INSERT<br/>AIR POLISHING NOZZLE</p> | <ul style="list-style-type: none"> <li>• Prophy</li> </ul>  |

Note: This comparison chart provides the hand scaling equivalent to DENTSPLY Cavitron® ultrasonic inserts for mechanical removal of deposits only. It is important to note that ultrasonic scaling provides further benefits to hand scaling (ie. lavage, efficiency, comfort, ergonomics, etc).