

Surf the canal
with confidence



wave • one[®]
GOLD



DENTSPLY

Small things. Big difference

THE MAKERS OF

PROTAPER[®] | CAVITRON[®] | WAVEONE[®] | SDR[®] | CERAM-X[®] | NUPRO[®]

Keeping the WaveOne® simplicity, with a golden host of benefits


wave • one®
GOLD

1

Reinforces patient safety

- Primary WaveOne GOLD file is 50% more resistant to cyclic fatigue than WaveOne Primary file
- Reduced screwing effect compared to standard rotary systems

Cyclic fatigue resistance




2

Covers a wider range of canal morphologies

- Enhanced file flexibility thanks to the GOLD technology
- Extended size range: small, primary, medium, large

Flexibility




3

Shortens the shaping time

- Takes the cutting efficiency to a higher level
- A single file per treatment translates into faster shaping time and thus more time for irrigation

Shaping working time



Engineered with the latest reciprocating technology

DENTSPLY has teamed up with some of the best endodontic professionals to develop **wave•one[®]** GOLD

The design team speaks:


Dr C. Ruddle, Dr J. Webber, Dr W. Pertot, Dr S. Kuttler

Shaping canals has never been easier


Since its launch in 2011, DENTSPLY's WaveOne has set a new standard in root canal shaping. All around the world general dentists and endodontists have voiced their appreciation of the simplicity WaveOne brings to their efforts to provide quality endodontic treatment to their patients.

We have been working on a new and exciting WaveOne project to bring you WaveOne GOLD.

WaveOne GOLD brings to the endodontic procedure a metallurgically advanced single-file technique for shaping canals. Further, we optimised the tip diameters, tapers, and the cross-section to produce a file that really improves safety, efficiency and flexibility when preparing canals. WaveOne GOLD takes simplicity and treatment success to another level.





Courtesy of Dr J. Webber



Watch the video: www.dentsply.co.uk/gold

Simplicity that facilitates treatment focus

The majority of root canals can be shaped with a single wave•one[®] Primary file
GOLD



Ready-to-use + Single use




Most dental practitioners have the twofold aim of reducing overall treatment costs while improving safety and efficiency. With this in mind, WaveOne GOLD files are supplied in sterilised blister packs. Dentists will benefit from:

- Simplicity, thanks to the ready-to-use sterile packaging
- Increased safety margin, by always using files with maximum fatigue resistance potential
- Maximised cutting efficiency

| One complete solution

With **wave•one^{GOLD}** DENTSPLY promotes a comprehensive treatment solution and optimises predictably successful outcomes



Product references

wave•one[®]
GOLD



Files



Small reciprocating file:
3x small (020 07) – 21mm – A 0753 221 G0S 03
3x small (020 07) – 25mm – A 0753 225 G0S 03
3x small (020 07) – 31mm – A 0753 231 G0S 03



Primary reciprocating file:
3x small (025 07) – 21mm – A 0753 221 G0P 03
3x small (025 07) – 25mm – A 0753 225 G0P 03
3x small (025 07) – 31mm – A 0753 231 G0P 03



Medium reciprocating file:
3x small (035 06) – 21mm – A 0753 221 G0M 03
3x small (035 06) – 25mm – A 0753 225 G0M 03
3x small (035 06) – 31mm – A 0753 231 G0M 03



Large reciprocating file:
3x small (045 05) – 21mm – A 0753 221 G0L 03
3x small (045 05) – 25mm – A 0753 225 G0L 03
3x small (045 05) – 31mm – A 0753 231 G0L 03

Assortment reciprocating file:
1x small/primary/medium/large – 21mm
A 0754 221 G9A 03
1x small/primary/medium/large – 25mm
A 0754 225 G9A 03
1x small/primary/medium/large – 31mm
A 0754 231 G9A 03



WaveOne GOLD paper points



Box of 180x sterile paper points:
small - A 175W 000 00S 03
primary - A 175W 000 00P 03
medium - A 175W 000 00M 03
large - A 175W 000 00L 03



WaveOne GOLD gutta-percha points



Box of 60x variable tapered gutta-percha points:
Small - A 175X 000 00S 00
Primary - A 175X 000 00P 00
Medium - A 175X 000 00M 00
Large - A 175X 000 00L 00



Thermafil[®] for WaveOne GOLD



Blister of 6x obturators:
Small obturators - A 1750 B00 00S 00
Primary obturators - A 1750 B00 00P 00
Medium obturators - A 1750 B00 00M 00
Large obturators - A 1750 B00 00L 00

Box of 30x obturators (5 blisters):
Small obturators - A 1750 P00 00S 00
Primary obturators - A 1750 P00 00P 00
Medium obturators - A 1750 P00 00M 00
Large obturators - A 1750 P00 00L 00

Motors



X-Smart iQ™ Motor - A1054
X-Smart® Plus Motor - A1032

Oven ThermaPrep^{®2}



A 0176 230 002 33 – UK (230V)

Kits



A10541 X-Smart iQ™ Starter Kit:
1 Motor with contra-angle
1 Pack of 15 Primary WaveOne GOLD
1 Pack of 16 Assorted WaveOne GOLD
1 Blister of 3 ProGilder[®]
1 iPad Protective Case

A10358 X-Smart® Plus Starter Kit:
1 Motor with contra-angle
1 Pack of 15 Primary WaveOne GOLD
1 Pack of 16 Assorted WaveOne GOLD
1 Blister of 3 ProGilder[®]

DENTSPLY

Small things. Big difference

THE MAKERS OF

PROTAPER[®] | CAVITRON[®] | WAVEONE[®] | SDR[®] | CERAM-X[®] | NUPRO[®]

CE
0459