

# ProSystem GT<sup>®</sup> Endodontic System



With the ProSystem GT endodontic system, you can handle almost any case requiring root canal therapy, and do it predictably and efficiently. The system contains all the accessories you'll need to get started:

Torque-control motor  
(Aseptico Endo DTC or Endo ITR)

Mini-head push-button contra-angle

GT<sup>®</sup> obturator oven

Two sets of ProFile<sup>®</sup> GT<sup>®</sup> rotary files and  
GT<sup>®</sup> obturators

ProLube<sup>®</sup> root canal conditioner

ProRinse<sup>®</sup> endo  
irrigation needles

Lexicon<sup>®</sup> K-Files

Lexicon<sup>®</sup> FlexSSK<sup>®</sup> files

Lexicon<sup>®</sup> barbed broaches

Lexicon<sup>®</sup> Amalgam burs

Lexicon<sup>®</sup> Diamond Access burs

Lexicon<sup>®</sup> Multi-burs

Lexicon<sup>®</sup> Pulp Shaper burs

Prepi<sup>®</sup> burs

GT<sup>®</sup> Fiber Post Sample Pack

Endo measuring device

ThermaSeal<sup>®</sup> Plus root canal sealer

Complete instruction manual

Tip cards

Introduction DVD




# ProFile GT Rotary Files

## Know How to Handle the Curves.

The ProFile GT file's radial lands and U-shaped flutes lift debris coronally. These design features help keep the instrument centered in the canal preventing stripping and zipping of the canal. They also reduce the risk of transporting, ledging and perforating the canal.



Flat radial lands



U-shaped flutes

ProFile GT rotary files are available in 20 Series, 30 Series, 40 Series and Accessory Sizes. Shaping root canals with GT rotary files generally requires one to three files to cut the initial shape and, after gauging, one additional file to adjust the shape, if necessary. The robust line of files allows for a wide variety of shaping options.

### 20 SERIES



### 30 SERIES



### 40 SERIES



### ACCESSORY FILES



To find out more about our comprehensive line of endodontic products, contact us today.

**DENTSPLY**  
TULSA DENTAL  
SPECIALTIES

Endodontics & Implants

Our Quality. Your Success.  
Their Satisfaction.

#### DENTSPLY Tulsa Dental Specialties

DENTSPLY International, Inc.  
5100 E. Skelly Drive, Suite 300  
Tulsa, OK 74135  
1-800-662-1202  
1-800-597-2779 (fax)  
[www.tulsadentalspecialties.com](http://www.tulsadentalspecialties.com)



See Instructions For Use

Rx only

For Dental Use Only

© 2009 DENTSPLY International, Inc.  
Manufactured in the U.S.

SIZED TO FIT EVERY CASE.