



Biosint® Extra

Steady as a rock

DeguDent

A Dentsply International Company

Biosint® Extra

34001/0703/H

the plaster free phosphate bonded investment material
for partial dental work

- high accuracy of fit
- optimal grain distribution
- exact reproduction of detail
- smooth casting surfaces
- expansion control with Biosol® expansion liquid

Technical data

Mixing ratio (powder : liquid)	100 g: 15 ml
Total expansion (linear)	0.95 – 1.65 %
Mixing time (under vacuum)	60 sec
Processing time span	3 – 5 min
Compressive strength (depending on concentration of the mixing liquid)	15 – 20 N/mm ²

Supply form

Box of 45 sachets, each containing 400 g
with Biosol® E Liquid

Product No.

2554 0001

Represented by:

DeguDent GmbH
P.O.Box 13 64 · 63403 Hanau
GERMANY
www.degudent.com

DeguDent
A Dentsply International Company