

Multi-Vest[®]

High Heat Casting Investment

Contents: 33 lbs. (14.85 kg.)

For Dental Use Only

Precaution: Inhalation of powder may be hazardous to your health.

For casting alloys of the cobalt-chromium system and other high fusing precious and non-precious alloys.

Storage: Investment materials are subject to deterioration when exposed to the atmosphere. Store this product at 59°F to 86°F. Keep container tightly closed once seal is broken to protect contents from moisture and high humidity.

DENTSPLY Trubyte
DENTSPLY International Inc.
York, PA 17405-0872
1-800-786-0085

www.trubyte.dentsply.com
or www.dentsply.com

Made & Printed in U.S.A.
4294-B Rev. B (2/03)

ReOrder
#99010

DENTSPLY
TRUBYTE

WARNING

Contains respirable crystalline silica (RCS). Do not breathe dust. May cause delayed lung injury (silicosis, pneumoconiosis). The I.A.R.C. (International Agency for Research on Cancer) reports (I.A.R.C. Monograph 68) there is sufficient evidence in humans for the carcinogenicity of inhaled crystalline silica in the forms of quartz or cristobalite from occupational sources. The NTP (National Toxicology Program) reports (Ninth annual report on Carcinogens) that RCS is known to be a carcinogen based on sufficient evidence from studies in humans indicating a causal relationship between exposure to RCS and increased lung cancer rates in workers exposed to crystalline silica dust. Follow OSHA Safety and Health Standards for crystalline silica. See Material Safety Data

NPCAHMIS CODE		NOTE
Health	3	Prolonged unprotected overexposure may result in nonreversible silicosis. Sufficient evidence shows crystalline silica causes cancer in humans.
Flammability	0	Non-flammable
Reactivity	0	Non-reactive
Personal Protection	*	Use NIOSH/MSHA/OSHA approved respirator.



See Material Safety Data Sheet (MSDS) for Detailed Information.