



Trubyte® Waxes

DENTSPLY Prosthetics offers a wide variety of high-quality waxes that are easy to handle and provide consistent quality. Trubyte® waxes are specifically formulated and designed to provide laboratories the right material for various indications.

TruWax® Baseplate Wax



TruWax Baseplate Wax is a high-quality wax that carves easily at room temperature and boils-out cleanly. Color: Pink

7709002	Regular – 1 lb. Pkg.
7708902	Extra Hard – 1 lb. Pkg.
7709003	Regular – 5 lb. Pkg.
7708903	Extra Hard – 5 lb. Pkg.

NeoWax® Flexible Baseplate Wax



NeoWax Flexible Baseplate Wax is a popular all-season, all-purpose, high-grade synthetic wax that is tough, pliable and boils-out cleanly. Color: Bright Pink

774140101	0.050" Thick – 5 lb. Pkg.
774140102	0.070" Thick – 5 lb. Pkg.

DENTSPLY® Bite Rim Sticks



Preformed wax for establishing vertical dimension, and arranging artificial denture teeth. Color: Pink

77315	45 per Pkg. - 0.45" Ridge Form x 0.45" Flat x 4.5"
-------	--

DENTSPLY® Set-Up Wax



DENTSPLY Set-Up Wax is specially formulated for arranging and articulating artificial teeth. It carves easily and boils-out cleanly. Available in sheets and sticks. Color: Salmon

7707901	Sheets – 5 lb. Pkg.
77309	Sticks – 32 per Pkg. - 1/4" x 1/4" x 4 3/8"

DENTSPLY® Sticky Wax



A highly tenacious wax used to hold material or pieces in position. Flows easily, sets firmly, and is rigid at room temperature. It flushes clean without residue. Color: Creamy Yellow

7709302	7 1/4 oz. Can
77310	1 lb. Pkg. - Sticks - 0.25" Diameter x 4.5"


Success® Sprue Wax Sticks



Specifically designed to make sprues for the Success Injection System. Color: White

904584	96 sticks – "D" shape 7mm x 180mm
--------	-----------------------------------

DENTSPLY® Utility Wax



A tacky, flexible wax for border moulding and beading impressions. Available in red and white in small or large strips and in square ropes. Colors: Red and White

STRIPS

77301	Small Strips – 198 Strips Red 0.125" Diameter x 11"
77303	Large Strips – 90 Strips Red 0.1875" Diameter x 11"

SQUARE ROPES

77305	Red – 55 Ropes 3/16" x 3/16" x 11"
77306	White – 55 Ropes 3/16" x 3/16" x 11"

DENTSPLY® Boxing Wax



A tacky, flexible wax for boxing impressions. Color: Red

77307	Regular Red – 1 lb. Pkg. – 1.5" x 0.1" x 11"
77308	Extra Thin Red – 1 lb. Pkg. – 1.5" x 0.05" x 11"

DENTSPLY

**Trubyte®
Waxes**

DENTSPLY International, **DENTSPLY Trubyte**,
570 West College Avenue, York PA 17405 1-800-786-0085
Distributed By: DENTSPLY Canada,
161 Vinyl Court, Woodbridge, Ontario L4L 4A3 1-800-263-1437
prosthetics.dentsply.com