



Porcelain Brushes with remarkable PERMALON™ fibers

PERMALON™* is a remarkable man-made fiber that has captured the performance of sable and made it long-lasting. Compared to PERMALON fibers, sable bristles rapidly lose their effectiveness for porcelain handling. When this happens, the frayed worn fibers no longer provide the precise control you need for smooth application and shaping.

In contrast, BIOBOND Brushes with PERMALON fibers resist fraying and separating. The PERMALON fibers retain their “vitality” for an amazingly long time, keeping you always in precise control of your porcelain. It’s almost like working with a brand new brush on every case. And the natural moisture absorbency of BIOBOND Brushes actually draws water from your build-up. The result is easier condensing, less chance of steam fissures and reduced shrinkage.

Try one BIOBOND Brush and you’ll want to try them all.

Fifteen different hand shaped tips on balanced birch shafts assure you of finding a BIOBOND Brush that is perfect for opaque application, porcelain build-up, staining, shaping or whatever your craft demands. In metal departments, BIOBOND Brushes are ideal for applying investment to delicate wax patterns. Denture technicians will also find BIOBOND Brushes useful in painting separator on casts and molds and many other applications.

Request Biobond for lasting consistency stroke after stroke. That’s LabSmarter®.

	Item No.	Brush No.	Use	
ROUNDS	9482601	000	Staining/Opaquing	
	9482606	00	Staining/Opaquing	
	9482611	0	Staining/Opaquing	
	9482616	1	Staining	
	9482621	2	Staining/Porcelain build-up	
	9482626	3	Porcelain build-up	
	9482631	4	Porcelain build-up	
	9482636	6	Porcelain build-up/ Mould, Cast painting/ Cleaning	
	9482661	8	Porcelain build-up/ Mould, Cast painting/ Cleaning	
	9482666	11	Porcelain build-up/ Mould, Cast painting/ Cleaning	
WHIPPERS	9482671	14	Porcelain smoothing/ Mould, Cast painting/ Cleaning	
	9482676	16	Porcelain smoothing/ Mould, Cast painting/ Cleaning (Red Sable)	
FLATS	9482641	2	Porcelain build-up	
	9482646	4	Porcelain build-up/ Mould, Cast painting	
	9482651	6	Porcelain build-up/ Cleaning	

*PERMALON™ is not a registered trademark of DENTSPLY International.

DENTSPLY International, Inc., 570 West College Avenue, York, Pennsylvania 17405-0872 www.dentsply.com