

# Faskut Abrasive Wheels and Discs



For grinding, cutting, or polishing porcelain and metal

- Available in graded grits and color-coded for easy identification.
- Provide quick grinding, cutting and smooth polishing with low heat build-up.

Item No.	Description
9890201	#2 Wheel (Dark Gray) Silicon Carbide Wheels: 5/8" x 1/16" (25/PK) For Porcelain Grinding
9890301	#3 Wheel (Green) Silicon Carbide Wheels: 7/8" x 1/8" (25/PK) For Fine Polishing of Porcelain or Metal
9890401	#4 Wheel (Green) Silicon Carbide Wheels: 5/8" x 1/16" (25/PK) For Fine Polishing of Porcelain or Metal
9890501	#5 Wheel (Blue) Silicon Carbide Wheels: 7/8" x 1/8" (25/PK) For Extra Fine Polishing of Porcelain or Metal

Item No.	Description
9890601	#6 Wheel (Blue) Silicon Carbide: 5/8" x 1/16" (25/PK) For Extra Fine Polishing of Porcelain or Metal
9890701	#7 Wheel (Blue) Silicon Carbide: 5/8" Knife Edge (25/PK) For Extra Fine Polishing of Porcelain or Metal
9890801	#8 Disc (Grayish-Brown) Aluminum Oxide: 7/8" x .025" (100/PK) For Metal Cut Off
9890811	#8C Cup Disc (Grayish-Brown) Aluminum Oxide: 7/8" x .025" (50/PK) For Metal Shaping
9890901	#9 Extra-Thin Disc (Dark Grayish-Brown) Aluminum Oxide: 7/8" x .010" (25/PK) For Cutting Metal
9891001	#10 Wheel (Brown) Aluminum Oxide: 7/8" x 1/8" (25/PK) For Metal Reduction
9891101	#11 Wheel (Copper) Aluminum Oxide (non-contaminating) 7/8" Knife Edge (20/PK) For Metal Finishing
9891201	#12 Wheel (Copper) Aluminum Oxide (non-contaminating) 5/8" Knife Edge (20/PK) For Metal Finishing
9891301	#13 Cup Disc (Dark Gray) Silicon Carbide: 7/8" x .025" (100/PK) For Cutting Porcelain
9891401	#14 Extra-Thin Disc (Gray) Silicon Carbide: 7/8" X .010" (25/PK) For Cutting Porcelain

Non-smearing, no grabbing, less vibration...



DENTSPLY International, Inc.  
DENTSPLY Trubyte  
570 West College Avenue  
York, Pennsylvania 17405  
1.800.786.0085  
prosthetics.dentsply.com



# Faskut® Carbide Cutters

## Faskut Luci-Sof® Burs



	Head /Size 1/10mm	Head Length (mm)	Max RPM (000)
904763 <b>SMALL</b>	40	12.0	50
904764 <b>MEDIUM</b>	60	13.3	50
904765 <b>LARGE</b>	60	13.4	30

**Luci-Sof Burs:** for the finishing of silicone denture base materials, particularly Luci-Sof® Soft Denture Liner.

- Smooth-cutting honeycomb pattern inhibits clogging, grabbing and vibration, yet provides excellent durability.
- Uniform pattern over entire cutting surface lets Luci-Sof Burs be used at any angle.
- High-quality carbide metal provides long-lasting sharpness of cutting facets.

## Faskut Series 200 Carbide Cutters

**Series 200:** for all types of acrylic, precious or non-precious crown and bridge alloy finishing; trimming stone dies and casts.

**Specifications:** 3/32” diameter shank, H.P., 18–20,000 RPM operating capacity, two “diamond cut” blade designs: FINE (F) and COARSE (C)

Item No.	Cutter No.	Shape
9893101	212F	Pear
9893201	212C	Pear
9893103	214F	Taper
9893203	214C	Taper
9893105	216F	Cone
9893205	216C	Cone
9893107	218F	Flame
9893207	218C	Flame
9893109	220F	Egg
9893209	220C	Egg
9893111	222F	Round
9893211	222C	Round

9893113	1 each: 212F, 214F, 216F, 218F, 220F, 222F	FINE Assort.
9893213	1 each: 212C, 214C, 216C, 218C, 220C, 222C	COARSE Assort.



# Faskut Super Band™

## Arbor Bands and Mandrels



For all types of acrylic finishing

- Resin-bonded, fabric-backed abrasive for longer life.
- New layers of aluminum oxide/silicon carbide abrasive are exposed as the arbor band wears, creating greater efficiency through less time changing bands.
- Available in two diameters to meet the varying needs of technicians.
- Faskut Mandrels are precision ground from high-speed steel for dependability and performance. Head screw holds arbor band securely in position.

## Faskut Arbor Bands

Item No.	Description
99122	3/8” bands; tub of 100
99124	3/4” bands; tub of 100
99125	3/4” bands; tub of 500



## Faskut Mandrels

Item No.	Description
9890001	#00 Disc Mandrel (4/PK)
98915	#00L Disc Mandrel Left Hand Thread (4/PK)
9890011	#15 Tapered Shank Mandrel (6/PK)



# Faskut Aluminum Oxide

## Abrasive Points



For all types of acrylic finishing

- Contamination-free finishing.
- Self-cleaning and self-cooling, Faskut Points can withstand high RPMs over extended grinding times.

Item No.	Point No.	Shape
99110A	#2	Small Pear
99111A	#4	Flame
99112A	#5	Round
99113A	#6	Cone
99114A	#7	Large Pear
99115A	#8	Bud
99116A	Assort.	Contains one of each above

# Faskut Large Abrasive Wheels



- Faskut Large Abrasive Wheels are ideal for trimming hard or soft denture base materials, grinding stone casts, and shaping baseplate material.

Item No.	Description
904557	Faskut Beveled 3” Abrasive Wheel
99117A	Faskut Large 3” Abrasive Wheel
904582	Beveled Faskut Wheel Mandrel