



## Increased Productivity in three Translucent Shades.

TruTone is a fast-cure denture base material that provides desirable handling characteristics, increased productivity, and esthetic shades.

TruTone can be packed about 10 minutes after mixing, yet the long working time of the material provides up to 60-minutes of packing time.

For fast turnaround times in a high-production environment, TruTone's fast-cure cycle with a 20-minute boil can be used. When using the fast-cure technique, the denture can be packed, cured and de-flasked in less than one hour. The material can also be processed using conventional curing cycles with a processing time of 2 to 5 hours or overnight, depending on the processing method.

Available in three esthetic translucent shades – Light Pink, Medium Pink and Dark Pink – the fibered material enables laboratories to provide a natural looking denture.

Desirable handling characteristics, increased productivity, and esthetic shades are the TruTone advantages that enable laboratories to LabSmarter®.

### TruTone® Denture Base Resin provides a number of advantages:

- *Rapid gelation action*
- *Long working time*
- *Fast-cure cycle for increased productivity*
- *A choice of three esthetic translucent shades: Light Pink, Medium Pink and Dark Pink*
- *Fibered resins for a natural appearance*

DENTSPLY International  
DENTSPLY Trubyte  
570 West College Avenue  
York PA 17405-2558  
1-800-786-0085  
prosthetics.dentsply.com

Product/Package	Item No.
<b>105-Unit Powder (5 lbs.)</b>	
Light Pink	89252
Medium Pink	89253
Dark Pink	89254
<b>540-Unit Powder (25 lbs.)</b>	
Light Pink	89255
Medium Pink	89256
Dark Pink	89257
<b>300 mL Liquid</b>	89258
<b>950 mL Liquid</b>	89259
<b>3784 mL Liquid</b> (Four 1-qt. containers)	89260