





Trubyte® Biotone® Plastic Anteriors & Posteriors

Plastic Upper Anteriors

MOULD NUMBER	Width of Central (mm)	Length of Central Without Collar (mm)	Width of Anteriors on Curve (mm)	Articulates with Lower Mould
133	7.5	8.75	45.0	2C 2D
134	8.0	9.25	48.5	3D
135	8.75	10.0	52.5	46
136	9.0	10.25	54.5	4H 3N
137	9.5	11.0	57.0	47
1B	7.0	9.0	45.00	2C
1H	8.5	12.0	53.0	1H
2D	7.75	9.5	48.0	2D
3D	8.0	9.25	49.0	3D
4H	8.75	10.25	53.5	4H
A84	8.25	11.5	48.5	1N
1N	8.25	9.75	49.0	1N 46
2N	7.5	10.25	47.0	2N
2P	8.5	11.25	52.0	2P
A23	7.25	9.0	46.0	2C

Plastic Lower Anteriors

MOULD NUMBER	Width of Central (mm)	Length of Central Without Collar (mm)	Width of Anteriors on Curve (mm)	Articulates with Lower Mould
2C	4.50	8.25	33.50	
3M	5.00	8.25	33.50	
2D	5.00	9.25	34.00	
2N	5.00	9.0	35.00	
1N	5.00	9.25	36.50	
3D	5.00	8.75	36.50	
3N	5.00	9.5	36.50	
35	5.50	10.0	37.50	
2E	5.50	9.5	37.50	
2P	5.50	9.75	38.50	
46	5.50	8.75	39.00	
4H	5.50	9.25	39.00	
3P	6.00	10.0	39.50	
26	5.50	11.5	39.50	
47	5.75	9.5	41.50	
1H	5.75	10.75	41.50	
3R	5.75	11.75	42.50	

Plastic Posteriors

	MOULD NUMBER	X	Y	U	L
Blended 33°	30M	30.0	32.4	8.2	8.7
	30L	30.0	32.8	8.0	9.2
	30LS	30.0	32.0	8.1	8.8
	32M	32.0	34.1	8.5	9.0
	32LS	32.0	34.8	8.6	8.8
	32L	32.0	35.0	8.2	8.7
Dentron 20°	34M	34.0	35.9	9.1	9.8
	34L	34.0	36.4	9.0	9.8
	34LS	34.0	36.5	8.7	9.8
	29M	29.0	32.9	7.7	9.2
	29L	29.0	33.2	8.7	9.1
	31M	31.0	35.2	8.1	9.9
Functional 10°	31L	31.0	35.5	9.0	9.8
	33M	33.0	36.9	8.9	10.2
	33L	33.0	37.5	9.6	10.2
	F30	30.0	34.0	8.6	8.3
	F32	32.0	33.6	8.7	9.3
	F33	33.0	35.0	9.5	9.7
Rational 0°	F34	34.0	37.1	9.0	10.0
	29M	29.0	31.9	7.0	8.1
	29L	29.0	31.9	7.2	8.0
	31M	31.0	34.6	7.2	8.1
	31L	31.0	33.6	7.5	8.4
	33M	33.0	35.4	7.9	8.6
Rational Block 0°	33L	33.0	35.9	8.4	9.0
	31L	31.0	33.9	7.3	7.8
	33M	33.0	36.4	8.5	7.9

Anatomical Occlusal

Blended 33° Plastic

30M,LS*
32M,LS*
34M,LS*

62P, 65P, 66P, 67P, 77P, 81P

Semi-Anatomical Occlusal

Dentron®
20° Plastic

29M,L
31M,L
33M,L

59, 62, 65, 66, 67, 77, 81

Functional 10° Plastic

F30, F32,
F33, F34

59P, 62P, 65P, 66P, 67P, 77P, 81P

Flat Plane Occlusal

Rational 0° Plastic

29M,L
31M,L
33M,L

62P, 65P, 66P, 67P, 69P, 77P, 81P


**Rational 0° Block

31L
33M

62P, 65P, 66P, 67P, 77P


*LS indicates Long Buccal-Short Lingual moulds.

**Supplied in sets of 1 X 8 only
(consisting of 2 blocks of 4 right and left - upper and lower).



The average width of a 1 x 4 set of upper posteriors = X;
the lower posteriors = Y.
The average depth of the upper left first molar = U;
the lower left first molar = L.

ISO3336: 1993
ANSI/ADA Specification no. 15-1999
Processing instructions are available upon request.



Trubyte® Biotone® Plastic Anteriors & Posteriors

Plastic Mould Chart

Suggested Articulations & Combinations

Upper Mould Number	Articulates with Lower Mould	Combines with Trubyte® Posterior 33° Mould	Combines with Trubyte® Posterior 20° Mould	Combines with Trubyte® Functional® Posterior Mould	Combines with Trubyte® Rational® Posterior Mould
133	2C, 2D	30M	29M	F30	29M
134	3D	30M	31M	F33	31M
135	46	32M	33M	F33	33M
136	4H, 46	34M	33M	F33	33M
137	47	34L	33M	F34	33M
1B	2C	30M	29M	F32	29M
1H	1H	34L	33L	F33	33L
2D	2D	30L	29L	F32	29L
3D	3D	30M	31M	F33	31M
4H	4H	32M	33M	F33	33M
A84	1N	30L	31L	F32	31L
1N	1N, 3N	30L	31M	F32	31M
2N	2N	30L	29L	F30	29L
2P	2P	30L	31L	F33	31L
A23	2C	30M	29M	F30	29M
A24	1N, 3N	30M	29M	F30	29M
A25	2E	30M	31M	F33	31M
A26	46, 2E	30M	31M	F33	31M
3M	3M	30M	29M	F30	29M
3N	3N	30L	31M	F32	31M
3P	3P	34M	33M	F33	33M
4M	3M	30M	29M	F30	29M
4N	2N	30L	31M	F30	31M
4P	2P	32L	33L	F33	33L
263	2N	30M	29L	F32	29L
264	35, 2E	30M	31L	F34	31L
266	26	34L	33L	F34	33L
267	3R	34L	33L	F34	33L
268	3R	34L	33L	F34	33L

Suggested Anterior-Posterior Shade Combinations

Trubyte® Biotone® Anterior Shade	Trubyte® Biotone® Posterior Shade	Trubyte® Bioform® Posterior Shade
59P	59P	B59
	62P	
61P	59P	B62
	62P	B63
62P	62P	B62
		B54
65P	65P	B65
		B56
66P	66P	B66
		B55
67P	67P	B67
		B56
68P	65P	B65
		B83
69P	69P	B69
		B96
77P	77P	B77
		B84
81P	77P	B81
	81P	B85

DENTSPLY
TRUBYTE

The chemistry of a great smile.™