

Try now the new, state-of-the-art  
core build-up material



core·x<sup>TM</sup> flow

Dual Cure Core Build-up Material and Cement for Endodontic Posts


- 1 Convenient automix syringe
- 2 Dual cure mode
- 3 High compressive strength

For better dentistry

**DENTSPLY**

## Product features at a glance

- Two indications in one product (core build-up and post cementation) help to save time and additional expenses
- Automix syringe guarantees consistent mix of material and predictable quality
- Speed of placement with high depth of cure: Light-cures in 20 seconds for 3 mm
- High compressive strength ("Cuts like dentine")
- Dual cure mode
- Radiopaque: X-ray detection for monitoring
- Stackable viscosity for ease of placement
- Compatible with dual-/self-cure adhesives (we recommend the use of XP BOND® with Self Cure Activator)



<sup>1</sup> LuxaCore Smartmix is not a registered trademark of DENTSPLY International Inc

## Packaging Information

### Core-X™ flow 606.67.330

4 syringes, 4.75 g each  
+ 20 Mixing/Intra Oral Tips

### Mixing/Intra Oral Tips for Core-X™ flow 606.20.116

40 pcs.



#### DENTSPLY DeTrey GmbH

De-Trey-Str. 1  
78467 Konstanz  
GERMANY  
www.dentsply.de  
Phone +49 (0) 75 31 5 83-0

#### DENTSPLY Ltd

Building 1  
Aviator Park, Station Road  
ADDLESTONE, KT15 2PG  
UNITED KINGDOM  
Phone +44 (0) 19 32 85 34 22

For better dentistry

**DENTSPLY**