




Easily identify setting stages with Jeltrate® Chroma!

Formulated to change colors for each preparation phase



0-30 sec.
Fuschia during mixing

30-45 sec.
Lavendar during tray
loading

45 sec-1 min
Light blue during
setting. Wait 1 minute
for mouth removal.

"Besides the great consistency, great accuracy, and great stability, there is the fun color changing component. The patients love to see it change from bright pink to lavender blue. The color changing component eliminates the guesswork on the time it needs to be in the mouth, and I believe that this is an asset. All in all Jeltrate Chroma Color is a product that I would recommend to any dental office as a great product. The accuracy, consistency, and fun color changing aspect are all wonderful."

Carol Guerette, C.D.A.
Hammondsport, NY

- 5-day dimensional stability
- 2-minute Fast set
- Easy mixing
- Non-slumping
- Consistent mix provides smooth models

Jeltrate®
Chroma
Dustless Color-Changing
Alginate Impression Material

Ordering Information

605710 **CASE OF 9 CANISTERS**
Fast Set Intro Kit
Reusable Canister contains:
1 Pouch Jeltrate® Chroma Material (454g),
1 Accessory Pack

605700 **CASE OF 12 POUCHES**
Fast Set (454g each)




Jeltrate® Jeltrate® DUSTLESS Jeltrate® Plus DUSTLESS

The high algin content of **Jeltrate® impression materials** provides for a quality impression without the excessive flow you will find in alginate materials with lower algin contents.

By eliminating excessive flow, Jeltrate® impression materials allow your patients a level of comfort that cannot be found in other alginate impression materials.



Jeltrate®



other alginate impression material



HIGH ALGIN CONTENT – Why is it important?

- Consistent mix.....Very smooth models
- No remakesSaves time
- High tear strengthEasily removed from undercuts without tearing
- Non-slump formulaPatient comfort

Ordering Information

Jeltrate®

Alginate Impression Material

CASE OF 9 CANISTERS
608503 Regular Set – Beige (454g ea)
608522 Fast Set – Pink (454g ea)

608005 **GROSS PACKAGE**
12 Twelve-Pouch Packages (21g ea)
1 Water Measure

NOTE: Pouch Packages not available in Fast Set

BULK PACKAGE
608058 Regular Set (22lbs / 10 kilos)
608247 Fast Set (22lbs / 10 kilos)

Jeltrate® Dustless

Dustless Alginate Impression Material

CASE OF 8 POUCHES
1 Accessory Pack
1 Empty Canister
605100 Regular Set (454g ea)
605200 Fast Set (454g ea)

605220 **CASE OF 20 POUCHES**
Fast Set (454g ea)

Jeltrate® Plus

Antimicrobial Dustless Alginate Impression Material

CASE OF 9 CANISTERS
605502 Regular Set – Beige (454g ea)
605602 Fast Set – Pink (454g ea)

605505 **GROSS PACKAGE**
12 Twelve-Pouch Packages (25g ea)
1 Water Measure

NOTE: Pouch Packages not available in Fast Set

BULK PACKAGE
605503 Regular Set (22lbs / 10 kilos)
605603 Fast Set (22lbs / 10 kilos)



DENTSPLY
CAULK

38 West Clarke Avenue / Milford, DE 19963
800.532.2855 / www.dentsply.com / www.caulk.com

All data on file. ©2013 DENTSPLY International. All rights reserved.
Form #505900 (R 12/3/13)

