

digit Power®  
Dispenser

*Moves With You*



Designed for Use with  
Aquasil® Ultra+  
Smart Wetting® Impression Material



## digit Power<sup>®</sup> Dispenser

Designed for Use with  
**Aquasil<sup>®</sup> Ultra+**  
Smart Wetting<sup>®</sup> Impression Material

### IT'S TIME TO ASK: IS YOUR DELIVERY SYSTEM WORKING FOR YOU?

Low viscosity impression delivery systems are built around the material—not the procedure or clinician—requiring extra steps to backfill a syringe or wrestle with a 50mL dispenser. Instead of suffering with wash delivery that doesn't work for you, it's time to find a better way.

**76%**

of clinicians want a process and material to deliver more value to their procedure\*

# Feel the thrill of technology that **moves with you.**

### SAY "HELLO" TO THE REVOLUTIONARY DIGIT POWER<sup>®</sup> DISPENSER.

With the touch of a button, it lets you focus on what matters most: a procedure done with precision. And when used with Aquasil<sup>®</sup> Ultra+ Smart Wetting<sup>®</sup> Impression Material it represents an impressive step forward in material delivery and procedural excellence.

**92%**

of clinicians find precision placement to be "critically important" to their procedure\*




# Designed to work with you.

## IMPRESSIVE ERGONOMICS


Dentists spend as many as 60,000 hours over the course of their career in improper body positions, negatively impacting their work and safety.<sup>1</sup>

Current dispensing options are designed around the material itself and not the clinician or the procedure. These delivery options are often heavy and cause physical strain. Ergonomically awkward body positioning is often required to achieve stabilization and control to ensure adequate delivery results, despite risk of injury from body position.<sup>2</sup>


With a unique pen grip design and power delivery, the digit Power Dispenser allows you to assume a more ergonomic position and focus on precise delivery of the wash material.




Device weighs up to 42% less limiting placement errors and repetitive fine motor fatigue\*



Up to 84% less force required to dispense than 50mL dispenser\*




Get up to 63% closer to the tooth prep for more control\*



Up to 64% less range of motion due to inadvertent wrist movement\*

Side-by-Side Comparison of Wash Material Delivery Devices



## GAME-CHANGING PRECISION

The digit Power® Dispenser shortens the distance between the hand and the tooth prep to just three inches. With the closer proximity, there is considerably less risk of tip movement while dispensing, resulting in a more accurate flow of wash material around the tooth prep.

50mL dispensing guns require the dentist's hand to work as much as eight inches away from the impression site and to apply approximately seven and a half pounds of pressure to pull the trigger. These and other factors reduce control and can compromise the procedure.



1. Gupta A, Bhat M, Mohammed T, Bansal N, Gupta G. Ergonomics in Dentistry. Int J Clin Pediatr Dent 2014;7(1):30-34.  
 2. American Dental Association. 2004. An Introduction to Ergonomics: Risk Factors, MSDs, Approaches and Interventions.

# A radical new approach to dispensing.

The digit Power® Dispenser's impressive set of design features complement the power of Aquasil® Ultra+ material to deliver results from the moment the material is dispensed. A slim, short mixing tip reduces distance to the impression site while reducing material waste. Single-use cartridges simplify the procedure from material loading to waste disposal. A built-in timer with work and set time notifications keeps you informed throughout the procedure. Two convenient speeds put control at your fingertips. This is simply the most advanced dispensing system on the market made for taking the most accurate impressions with utmost ease.

Designed for Use with  
**Aquasil® Ultra+**  
 Smart Wetting® Impression Material



Convenient unit dose delivery in multiple viscosities and set times

Slim mixing tip profile for access, visibility, and reduced waste

Convenient, autoclaveable cartridge cap system

Two-speed dispensing makes impression-taking flexible and convenient

Fine intraoral tip design with 360° rotation for optimal material placement

Slim, lightweight design for improved hand comfort


Three convenient timer modes: fast set, regular set, user set

Battery powered for ultimate portability

Large information display

Timer Mode Memory that remembers your last used timer mode upon restart

Work and set time notifications to improve compliance



Color options to match your operator/office decor






Optional foot control for variable speed via Bluetooth® technology¹

## ORDER INFORMATION

### digit Power Dispenser Complete Kit

1 handpiece, 1 charging base, power cord with power blades,  
3 cartridge caps, 50 disposable barriers, 1 manual and quick start guide.


678530	White	
678531	Orange	
678532	Mint Green	

### digit Power Dispenser Accessories

678536	Disposable Barriers for digit Power Dispenser (50 ct)
678538	Cartridge Cap Refill (3 ct)
678537	Bluetooth® Foot Control with power cord and power blades
678008	digit Intrасulcular mixing tips (25 ct)
678007	digit Standard Intraoral mixing tips (25 ct)

### Aquasil Ultra+ Impression Material for digit Power Dispenser




50 single-use cartridges, 50 intraoral mixing tips

	FAST SET	REGULAR SET
 XLV	678698	678693
 LV	678694	678696

### Aquasil Ultra+ Impression Material Refill Options




50mL Cartridge 4-pack Refills

4 cartridges (50mL ea), 12 mixing tips

	FAST SET	REGULAR SET
 Medium	678674	678673
 Heavy	678676	678675
 Rigid	678678	

DECA™ Cartridge Standard Packs

2 Cartridges (380mL ea), 20 Mixing Tips, 1 Bayonet Locking Ring

	FAST SET	REGULAR SET
 Medium	678697	678682
 Heavy	678684	678686
 Rigid	678688	

*Not All SKUs are available in all markets*

Contact your local Dentsply Sirona Restorative representative or  
visit [dentsplysirona.com/digitPowerDispenser](https://www.dentsplysirona.com/digitPowerDispenser) to learn more.



\*Data on file.

1. Bluetooth is not a registered Trademark of Dentsply Sirona  
©2017 Dentsply Sirona. All rights reserved.  
ML7860000B (8/23/17)

