

ATLANTIS™ abutments implant compatibility chart



Manufacturer	Implant	ATLANTIS™ Abutment		ATLANTIS™ Crown Abutment		ATLANTIS™ Conus Abutment		Screwdriver
						Custom	Over- denture	
		Titanium Gold- shaded ti	Zirconia	Titanium	Zirconia	Titanium Gold- shaded ti	Titanium	
BioHorizons	INTERNAL							
	Internal 3.5, 4.5, 5.7	x	x	x*	x	x	x*	0.050" (1.25 mm)
	Tapered Internal 3.5, 4.5, 5.7	x	x	x*	x	x	x*	Hex screwdriver
	EXTERNAL							
	External implant 3.5	x				x		0.050" (1.25 mm)
Biomet 3i	External implant 4.0, 5.0	x	x		x	x		Hex screwdriver
	INTERNAL							
	Certain™ 4/3, 5/4	x	x	x*	x	x	x*	0.048"
	Certain™ 3.25, 4.0, 5.0, 6.0	x	x	x*	x	x	x*	Hex screwdriver
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	x	x*	x	x	x*	
	Certain™ XP 4/5, 5/6	x	x	x	x	x	x*	
	EXTERNAL							
	MicroMiniplant™ 3.25	x				x		0.048"
	Miniplant™ 3.25	x	x		x	x		Hex screwdriver
Camlog	Osseotite™ 3.75, 4.0, 5.0, 6.0	x	x		x	x		
	Osseotite™ XP 4/5, 5/6	x	x		x	x		
	INTERNAL							
	Screw-line Implant 3.3	x	x	x		x	x	1.28 mm
	Screw-line / Root-line Implant 3.8, 4.3, 5.0, 6.0	x	x	x	x	x	x	Hex screwdriver
DENTSPLY Implants OsseoSpeed implants are part of the ASTRA TECH Implant System.	INTERNAL							
	ANKYLOS® C/X – C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)	x	x			x		
	ANKYLOS® C/X – /X 3.5, 4.5, 5.5, 7.0 (indexed)	x	x			x	x	1.0 mm Hex screwdriver
	OsseoSpeed™ EV 3.0 – Green	x		x		x	x	0.050"
	OsseoSpeed™ EV 3.6 – Purple	x	x	x	x	x	x	Hex screwdriver
	OsseoSpeed™ EV 4.2 – Yellow	x	x	x	x	x	x	
	OsseoSpeed™ EV 4.8 – Blue	x	x	x	x	x	x	
	OsseoSpeed™ EV 5.4 – Brown	x	x	x	x	x	x	
	OsseoSpeed™ Profile EV 4.2 Yellow	x	x	x		x	x	
	OsseoSpeed™ Profile EV 4.8 Blue	x	x	x		x	x	
	OsseoSpeed™/OsseoSpeed™ TX 3.0 – Yellow	x		x		x	x	
	OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 – Aqua	x	x	x	x	x	x	
	OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 – Lilac	x	x	x	x	x	x	
	OsseoSpeed™ TX Profile 4.5, 5.0	x		x		x	x	
	OsseoSpeed™ 4.0 – Blue (Previous)	x				x		
	XiVE® S 3.0	x		x		x	x	1.22 mm
	XiVE® S 3.4, 3.8, 4.5, 5.5	x	x	x	x	x	x	Hex screwdriver
	FRIALIT™ 3.4, 3.8, 4.5, 5.5	x	x	x	x	x	x	
	Keystone Dental	INTERNAL						
PrimaConnex™ SD 3.3, 3.5, Genesis™ 3.8		x	x	x	x	x	x	Keystone
PrimaConnex™ RD 4.0, 4.1, Genesis™ 4.5		x	x	x	x	x	x	Quad Driver
PrimaConnex™ WD 5.0, Genesis™ 5.5, 6.5		x	x	x	x	x	x	
EXTERNAL								
Restore™ SD 3.3, 3.4		x				x		0.048"
Restore™ RD 3.75, 4.0		x	x		x	x		Hex screwdriver
Keystone Dental	Restore™ WD 5.0, 5.5, 6.0	x				x		

All trademarks and company names are the property of their respective owners.
* 2 day extra Turnaround time.

Please turn over ➡

Manufacturer	Implant	ATLANTIS™ Abutment		ATLANTIS™ Crown Abutment		ATLANTIS™ Conus Abutment		Screwdriver
						Custom	Over- denture	
		Titanium Gold- shaded ti	Zirconia	Titanium	Zirconia	Titanium Gold- shaded ti	Titanium	
Nobel Biocare	INTERNAL							
	NobelActive™ NP 3.5	x	x	x	x	x	x	Nobel Biocare Unigrip
	NobelActive™ RP 4.3, 5.0	x	x	x	x	x	x	
	NobelReplace™ NP 3.5	x	x	x	x	x	x	
	NobelReplace™ RP 4.3	x	x	x	x	x	x	
	NobelReplace™ WP 5.0	x	x	x	x	x	x	
	NobelReplace™ 6.0	x	x	x	x	x	x	
	NobelReplace™ Conical Connection NP 3.5	x	x	x	x	x	x	
	NobelReplace™ Conical Connection RP 4.3, 5.0	x	x	x	x	x	x	
	EXTERNAL							
	Brånemark System™ NP 3.3	x				x		Nobel Biocare Unigrip
	Brånemark System™ RP 3.75, 4.0	x	x		x	x		
	Brånemark System™ WP 5.0, 5.5	x	x		x	x		
	SteriOss Replace™ 3.1, 3.5, 4.3, 5.0, 6.0	x				x		0.050” Hex screwdriver
	SteriOss HL™ 3.25, 3.8, 4.5, 5.0, 6.0	x				x		
Straumann	INTERNAL							
	Bone Level 3.3 NC (Narrow CrossFit™)	x	x	x	x	x	x	SCS screwdriver
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	x	x	x	x	x	x	
	Standard/Standard Plus 4.8 RN (Regular Neck)	x		x		x	x	
	Standard/Standard Plus 4.8 WN (Wide Neck)	x		x		x	x	
	EXTERNAL							
	Standard Plus 3.3 NN (Narrow Neck)	x				x		SCS screwdriver
Zimmer Dental	INTERNAL							
	Tapered Screw-Vent™ 3.5, 4.5, 5.7	x	x	x*	x	x	x*	1.25 mm Hex screwdriver
	Tapered Screw-Vent™ 3.5, 4.5, 5.7 friction fit	x*		x*		x	x*	
	Screw-Vent™ 3.3, 3.7, 4.5	x	x	x*	x	x	x*	
	Screw-Vent™ 3.3, 3.7, 4.5 friction fit	x*		x*		x	x*	
	SwissPlus™ 4.8	x		x		x	x*	

All trademarks and company names are the property of their respective owners.

* 2 day extra Turnaround time.