

Your individual options

Here you can view the complete list of individual options, and a checklist of what is included with basic units.

Intego

Patient chair		Ambi-dextrous
2 treatment positions, S and O position	■	■
Integrated junction box	■	■
Chair lift system 140 kg	■	■
Chair lift system 185 kg	□Er	□Er
Flat headrest	■	■
Double articulating headrest	□	□
Premium upholstery	■	■
Lounge upholstery	□	□
ErgoMotion	□Er	□Er
4-way foot switch	□	□
Booster cushion for small patients	□	□
Right arm rest	□Er	□Er
Left arm rest	□Er	□Er
Valve set for external device	□	□
Vacuum pressure reduction kit	□	□
Cover installation kit	□	□
Separate junction box	□	□
Adapter plate	□	□

Dentist element		Ambi-dextrous
Height-adjustable support arm with pneumatic brake	■	■
5 instrument positions	■	■
Additional camera holder for TS dentist element	□	—
Hydro-pneumatic supply for 2 drill drives w/o fiber optics	■	■
Pneumatic instrument foot control	■	■
C+ electronic foot control	□A, I, T	□A, I, T
EasyPad user interface with 2 user-specific settings	■	■
EasyTouch user interface with 4 user-specific settings	□T	□T
Sprayvit E with light	□A, I	□A, I
Sprayvit E with light (water button on the left)	□A, I	□A, I
Standard 3-way syringe	□S	□S
Turbine hose with light	□Sx2, A, I, T	□Sx2, A, I, T
Borden turbine hose	□	□
1 st BL ISO E motor with hose	□A, I, T	□A, I, T
2 nd BL ISO E motor with hose	□	□
SiroSonic L scaler	□I	□I
Third drill drive	□	□
Intra-oral light source kit for turbine and Sprayvit E	□S, A, I, T	□S, A, I, T
Torque control for endodontic treatments	□E	□E
eControl: Speed range of 100-40.000 rpm	□E	□E
Integrated ApexLocator	□E	□E
Endo 6:1 contra-angle handpiece	□E	□E
Intra-oral x-ray viewer	□	□
Purge and manual sanitization function	□	□
Flexible TS tray holder for 2 standard size trays	—	□
Additional TS tray holder for 2 standard size trays	□	—
Silicone pad for TS dentist element	□	□
Additional CS tray holder for 1 standard size tray	□	□
Additional CS tray holder for 2 standard size trays	□	□
Silicon pad for CS dentist element	□	□
Removable silicone sleeves for handles	□	□

Assistant element		Ambi-dextrous
Compact assistant element, swivel mounted directly on the cuspidor arm	■	■
Comfort assistant element, ideal for two-handed dentistry	□	□
4 instrument positions	■	■
Saliva ejector and HVE	■	■
Sprayvit E with light	□	□
Sprayvit E with light (water button on the left)	□	□
Standard 3-way syringe	□	□
2 nd HVE kit	□	□
Satelec Mini LED curing light	□	□
Touch pad	□	□
Light source for Sprayvit E	□	□
Silicone pad for Comfort assistant element	□	□

Water unit		Ambi-dextrous
Prepared for fresh water bottle	■	■
Cuspidor with manual swivel function	■	■
White cuspidor	■	■
Water unit without cuspidor	□	□
Prepared for standard wet suction	■	■
Prepared for standard wet suction w/ automatic suction control	□	□
Multi-station valve	□	□
Cuspidor valve for standard wet suction	□	□
Automatic separator	□	□
Amalgam separator	□	□
Fresh water bottle w/ol switch to public water supply	□	□
Fresh water bottle with switch to public water supply	□	□
Water heater	□	□
Water flush suction hose cleaning system	□	□

Operating light/Heliodont Plus options		Ambi-dextrous
LEDlight	□S, A	□S, A
LEDview Plus	□I, T	□I, T
LEDview Plus ceiling light	□	—
Adaptation for LEDview Plus/Heliodont Plus ceiling combination	□	—
Adaptation for Heliodont Plus unit model	□	—

Patient communication		Ambi-dextrous
USB interface in the dentist or assistant element	□SiR, SiS	□SiR, SiS
SiroCam F	□SiS	□SiS
SiroCam AF	□	□
SiroCam AF+	□	□
Monitor mount on the light support column	□SiS	□SiS
22" AC monitor (with integrated power supply unit)	□SiS*	□SiS*
Network interface	□SiR, SiS	□SiR, SiS

- Included in basic unit
- Option
- Not possible/applicable
- SiR Sivation Ready package
- SiS Sivation Start package
- SiS* Sivation Start package with monitor
- S Start Pro package
- A Advance Pro package
- I Innovation Pro package
- T Touch Plus package
- E Endo Pro package
- Er Ergo package