

THE DENTAL
SOLUTIONS
COMPANY™



Hugo or Carl

Which one will you choose?

dentsplysirona.com



Just have a seat

Hugo - the ergonomics expert

Hugo is a key element in the Dentsply Sirona ergonomics program and combines unique technologies. It promotes ergonomic and intuitive sitting and enables you to maintain a healthy posture throughout the treatment process. Dynamic, comfortable, hygienic, and versatile – the perfect complement to your treatment center.

- Innovative design and intelligent shape for more comfort when sitting
- OptiMotion supports your natural movements
- Thermo upholstery ensures your comfort when sitting
- Choice of two gas spring lengths for different heights



Hugo Manual

Manual height adjustment



Hugo Manual Plus

Manual height adjustment, with foot ring



Hugo Freehand

Height adjustment by foot pressure

Carl - the all-rounder

The Carl round stool ensures maximum freedom of movement and meets the highest demands in comfort and quality. The round seat features comfortable upholstery and is „seamlessly“ hygienic.

- Round seat with comfortable upholstery
- Height can be adjusted from any sitting position
- The backrest can be rotated 360°
- Choice of two gas spring lengths



Carl Manual

Manual height adjustment



Carl Manual Plus

Manual height adjustment, with foot ring





The color schemes

The elegant world

Elegant upholstery colors for a classic, timeless look.



Carbon

Aubergine

Platinum

Pacific

The natural world

Soft, warm shades for a soothing treatment atmosphere.



Basalt

Rosewood

Cherry

Sapphire

The vital world

Vivid colors and bright tones for a modern, lively atmosphere.



Spring

Sunset

Orchid

Lagoon

Optimized ergonomics – Adapted to you from head to foot

Hugo is a key element in the Dentsply Sirona ergonomics program and combines unique technologies. By promoting an intuitive, ergonomic sitting position, Hugo enables you to maintain a healthy posture throughout the treatment process. Dynamic, comfortable and versatile - the perfect complement to your treatment center.



The backrest – for more dynamic support

The dynamic Hugo backrest gives you more flexibility during treatment. The backrest follows your movements and provides support – while you work and when you relax. Your spine is automatically straightened, effectively preventing muscle strain.



The seat – for more comfort

The innovative OptiMotion mechanism with flexible tilt of the front of the left or right half of the seat ensures more freedom of movement to allow you optimal access to the patient. This simultaneously reduces pressure points on the thighs considerably. The new thermo upholstery provides increased comfort during prolonged treatment sessions.



The stool base – for more flexibility

The Hugo stool can be flexibly adjusted to fit the height of each member of your practice team. On request, the seat height can be adjusted easily and hygienically via a foot switch.



1 The intelligent backrest ensures good sitting comfort and supports every vertebra in your back in each treatment position.

2 Innovative high tech from the aeronautics industry: Reduces heat build-up in the seat. This provides a pleasant effect and ensures that you can sit comfortably even during long treatments.

3 The seat with the unique OptiMotion mechanism supports your natural movements, relieves pressure points, and ensures more freedom of movement.

4 The use of new technologies dispenses with the operating lever. Hugo adjusts automatically to each sitting position. This helps prevent faulty posture and health risks. The functions and design also support hygienic work.



5 At no extra charge you can choose between two different gas springs. This allows you to select the ergonomically correct stool height for your height.

Height of user* (cm)	Seat height
approx. 150–180	Short gas spring 405–540 mm
approx. 170–205	Long gas spring 505–645 mm

* These figures vary from user to user.

6 The replaceable step element on the foot control increases the service life even for frequent use.

7 The optional foot ring allows uncomplicated interaction, even if the dentist and assistant are of different heights. Both work without stress in an ergonomically correct position.

8 The new foot release allows hands-free, hygienic adjustment of the seat height.

Dentsply Sirona

Sirona Dental Systems GmbH
Fabrikstraße 31, 64625 Bensheim, Germany
dentsplysirona.com

Procedural Solutions

Preventive
Restorative
Orthodontics
Endodontics
Implants
Prosthetics

Enabling Technologies

CAD/CAM
Imaging
Treatment Centers
Instruments