

## SiroLaser Xtend – Technical data

Technical data	SiroLaser Xtend
Wavelength	970 +/- 15 nm
Operation modes	Continuous Wave, Chopped Mode, Peak-Pulse Mode
Power	5 W*
Peak pulse	~ 10 W*
Frequency	1-1.000 Hz*
Duty cycle	50 % (fixed)*
Weight	~ 1.0 kg (incl. handpiece)
Dimensions	~ 19.7 x 18.2 x 18.9 cm
Scope of supply software	
User profiles	2 x (fixed)*
History file	-*
Favorites	-*
Scope of supply hardware	
Optical fiber	320 µm (4 pcs.) and 200 µm (1 pcs.)
Single-use tips incl. bending tool	50 x
Options	
Upgrade packages	+
Rechargeable battery	+

\* Upgrade possible.

Accessories	REF
Handpiece sleeve for alternating operation	62 56 767
Set of 320 µm fibers (5 pcs.)	62 56 759
Set of 200 µm fibers (5 pcs.)	62 56 742
Fiber cutter	60 91 669
Package of single-use tips (100 pcs.) incl. bending tool	62 56 692
Laser safety goggles (transparent glass filter)	61 81 585
Laser safety goggles for lens wearers (transparent glass filter)	61 81 593
Laser safety goggles for patients (transparent glass filter)	60 47 745
Rechargeable battery	62 56 833

SiroLaser Xtend REF-No for the following countries:

Italy REF 63 13 824; Great Britain REF 63 13 832; Spain REF 63 13 758; Portugal REF 63 13 840.  
Further countries on request.

This product is supplied without a paper version of the instruction manual. This can be easily downloaded in pdf format from the internet at [www.sirona.com/manuals](http://www.sirona.com/manuals).

Please note the following guidelines:

