



 **SINGLE-UNIT CROWN**  
WORKFLOW SOLUTIONS

FINAL IMPRESSION SOLUTION

**CAPTURE  
EVERY DETAIL  
QUICKLY AND  
RELIABLY**

# Every detail captured quickly and reliably

**1 in 3** dentists have problems with crown fit due to impression errors.<sup>1</sup>

## A solution that works for your workflow.


Our Final Impression Solution creates a more efficient workflow by effectively stopping minor intraoral bleeding to create a clear field for both conventional and digital workflows. Our products help you capture the data you need so the lab can create the restoration you want.



### Hemoban™ Gel 25% Aluminum Chloride Hemostatic Solution


Effectively stop intraoral bleeding—precisely and efficiently—using a transparent gel that stays where you place it and rinses off easily with no residue. Improve your productivity with:

- A clear gel that won't stain
- Performance that won't interfere with bonding
- Hemostasis without coagulation, reducing the risk of secondary bleeding when using cord
- Viscous yet spreadable non-drip gel for easy placement and removal



Cord soaked with liquid hemostatic agent

- Hemostatic agent absorbed fully into cord
- Textured surface



Cord coated with Hemoban™ Gel

- Hemostatic agent available on surface of cord for optimal delivery
- Smooth surface for easy insertion and removal without re-occurring bleeding
- Coated and soaked for easy packing

1. Data on file

### Aquasil® Ultra+ Smart Wetting® Impression Material

Capture the right details every time

Only Aquasil® Ultra+ impression material has an optimized set of intraoral performance attributes—hydrophilicity, tear strength, and long work time—coupled with a precision delivery system for improved impression accuracy. Every time.

- Leading intraoral hydrophilicity for excellent wettability to minimize voids and bubbles while maximizing detail capture
- Leading Intraoral tear strength for an added layer of confidence in even the most complex cases

	INTRAORAL HYDROPHILICITY (2" @80% RH)	CURED-FILM HYDROPHILICITY (5" @50% RH)	INTRAORAL TEAR STRENGTH (PSI @ .265mm thickness)	24 HOUR TEAR STRENGTH (PSI @ .265mm thickness)
AQUASIL® ULTRA+	15°	5°	607	678
IMPRINT™ 4	19°	6°	441	539
TAKE 1™ ADVANCED™	72°	88°	196	338
EXAFAST™ NDS	68°	66°	242	273
FLEXITIME®	55°	58°	291	399
IMPRINT™ 3	70°	43°	505	614
IMPREGUM™	49°	54°	165	259
V-POSIL	52°	42°	414	404

### digit Power® Dispenser

Discover technology that moves with you. The digit Power® Dispenser makes precision easier than ever—at the touch of a button


#### Game-Changing Precision

The digit Power® Dispenser shortens the distance between the hand and the tooth preparation to just three inches. With the closer proximity, there is considerably less risk of tip movement while dispensing, resulting in a more accurate flow of wash material around the tooth preparation.



#### Impressive Ergonomics

Heavy 50mL dispensing guns keep your hands as much as eight inches away from the impression site and require you to apply more than seven pounds of pressure to pull the trigger. The lightweight, pen-grip design of the digit Power® Dispenser allows you to work far more ergonomically. Good for you. Good for your results.



### Primescan™

The latest scanner from Dentsply Sirona provides fast, accurate, and easy to use digital impressing. With large field of view and sharpness up to 20mm in depth, Primescan™ can easily capture full arch impressions in 30 seconds or less. Primescan™ enables restorative workflows with your lab or in-house manufacturing.



**HEMOBAN™ GEL**  
**25% ALUMINUM CHLORIDE HEMOSTATIC SOLUTION**

675110 Hemoban™ Gel Complete Package  
1 Hemoban™ Gel Syringe filled (30mL)  
20 Empty Syringes (1.2mL)  
20 Infusor Tips

**AQUASIL® ULTRA+**  
**SMART WETTING® IMPRESSION MATERIAL**

678695 Introductory Kit  
2 Cartridges Heavy Fast Set (50mL ea)  
1 Cartridge LV Fast Set (50mL)  
1 Cartridge XLV Fast Set (50mL)  
1 Bottle Caulk® Tray Adhesive (14mL)  
1 digit® LV/XLV Fast Set Evaluator Kit  
12 Mixing Tips and 6 Intraoral Tips

**digit Power® Cartridge Refill**

50 Unit Doses, 50 Intraoral Mixing Tips

678698 digit® Power XLV Fast Set (Orange)  
678693 digit® Power XLV Regular Set (Orange)  
678694 digit® Power LV Fast Set (Teal)  
678696 digit® Power LV Regular Set (Teal)

**4-Pack Cartridge Refill**

4 Cartridges (50mL ea), 12 Mixing Tips, 12 Intraoral Tips

678670 XLV Fast Set  
678672 LV Fast Set

**4-Pack Cartridge Refill**

4 Cartridges (50mL ea), 12 Mixing Tips

678676 Heavy Fast Set  
678678 Rigid Fast Set

**Standard Pack**

2 Cartridges (380mL ea), 20 Dynamic Mixing Tips  
1 Bayonet Locking Ring

678684 Heavy Fast Set  
678688 Rigid Fast Set

**DIGIT POWER® DISPENSER**

digit Power® Dispenser Complete Kit  
1 Handpiece, 1 Charging Base and Power Cord,  
1 Power Adapter and Blades, 3 Cartridge Caps,  
50 Disposable Barrier Sleeves

678530 White Handpiece  
678531 Orange Handpiece  
678532 Mint Green Handpiece

678537 Foot Control

**PRIMESCAN™ AC WITH CONNECT SOFTWARE**

For more information about Primescan™, consult your  
Dentsply Sirona Representative.

[For more information and to request a demo,  
contact your Dentsply Sirona Representative.](#)

Dentsply Sirona | 13320-B Ballantyne Corporate Place | Charlotte, NC 28277 | 1-800-532-2855 | [www.singleunitrestorations.com](http://www.singleunitrestorations.com)

Competitor brands are not registered trademarks of Dentsply Sirona. ©2019 Dentsply Sirona Inc. All rights reserved ML661002B (R 3/15/19)

