

SEM picture of SphereTEC* sphe prepolymerized fi Their shape and size distribution create a ball bearing-like ef that enables exc adaptation and sculptability





ructure of SphereTEC* er is able to bind more ee resin compared to ual fillers Therefore Neo ectra™ ST composite es not stick to hand

Neo Spectra[™] ST universal composite with SphereTEC® Technology.

Made with what matters.

Neo Spectra™ ST composite's performance is made possible by SphereTEC® filler technology, our proprietary method of manufacturing microscaled, well-defined spherical superstructures, comprised of submicron glass. The SphereTEC® fillers' morphology, particle size distribution and surface microstructure deliver the benefits that really matter to dentists:



OPTIMIZED HANDLING for easy placement and shaping



BLENDING ABILITY for simplified shade selection and matching



DURABILITY with high polish and stain resistance for lasting esthetics

Neo Spectra[™] ST

Universal Composite Restorative

HIGH VISCOSITY -

Compules[®] Tips Intro Kit (HV) 607.01.970

24 HV Compules® Tips (Gray) Assorted - (0.25 g each) (4 x A1; 4 x A2; 8 x A3; 4 x A3.5; 4 x A4) 1 Compules[®] Tips Gun i-shade label 1 Bottle Prime&Bond Universal[™] Adhesive (2.5 mL) Accessories

607.01.980 Syringe Intro Kit (HV) 5 HV Syringes (Gray) Assorted - (3 g each)

(A1: A2: A3: A3.5: A4) i-shade label 1 Bottle Prime&Bond Universal™ Adhesive (2.5 mL) Accessories

607.01.990 Syringe Eco Kit (HV)

6 HV Syringes (Gray) Assorted - (3 g each) (3 x A2; 3 x A3) i-shade label 1 Bottle Prime&Bond Universal[™] Adhesive (2.5 mL) Accessories

Refills (HV)	Compules [®] Tips	Syringe	
Shades	(16 x 0.25 g)	(1 x 3 g)	
A1 (= VITA ^{®1} shades A1; B1; C1)	607.01.971	607.01.981	
A2 (= VITA ^{®1} shades A2; B2; D2)	607.01.972	607.01.982	
A3 (= VITA ^{®1} shades A3; C2; D3; D4)	607.01.973	607.01.983	
A3.5 (= VITA®1 shades A3.5; B3; B4; C3)	607.01.975	607.01.985	
A4 (= VITA®1 shades A4; C4)	607.01.974	607.01.984	
BW (Bleach White)	607.01.976		

LOW VISCOSITY -

607.01.870 Compules[®] Tips Intro Kit (LV)

24 LV Compules® Tips (Blue) Assorted - (0.25 g each) (4 x A1; 4 x A2; 8 x A3; 4 x A3.5; 4 x A4) 1 Compules® Tips Gun i-shade label 1 Bottle Prime&Bond Universal[™] Adhesive (2.5 mL) Accessories

Syringe Intro Kit (LV) 607.01.880 5 LV Syringes (Blue) Assorted - (3 g each)

(A1; A2; A3; A3.5; A4) i-shade label 1 Bottle Prime&Bond Universal[™] Adhesive (2.5 mL) Accessories

607.01.890 Syringe Eco Kit (LV) 6 LV Syringes (Gray) Assorted - (3 g each) (3 x A2: 3 x A3) i-shade label 1 Bottle Prime&Bond Universal[™] Adhesive (2.5 mL)

Refills (LV)	Compules® Tips	Syringe	
Shades	(16 x 0.25 g)	(1 x 3 g)	
A1 (= VITA ^{®1} shades A1; B1; C1)	607.01.871	607.01.881	
A2 (= VITA ^{*1} shades A2; B2; D2)	607.01.872	607.01.882	
A3 (= VITA®1 shades A3; C2; D3; D4)	607.01.873	607.01.883	
A3.5 (= VITA®1 shades A3.5; B3; B4; C3)	607.01.875	607.01.885	
A4 (= VITA®1 shades A4; C4)	607.01.874	607.01.884	
BW (Bleach White)	607.01.876		

Accessories

Neo Spectra[™] ST Effects Universal Composite Restorative

Compules[®] Tips Intro Kit 607.01.940 24 HV Compules[®] Tips (Gray) Assorted - (0.25g each) (12 x D1; 8 x D3; 4 x E1) Illustrated Technique Guide i-shade label 607.01.950

Syringe Intro Kit (HV) 3 HV Syringes (Gray) Assorted - (3 g each) (1 x D1: 1 x D3: 1 x E1) Illustrated Technique Guide i-shade label

Refills (HV)	Compules [®] Tips	Syringe
Shades	(16 x 0.25 g)	(1 x 3 g)
D1 (use with A1; A2)	607.01.942	607.01.952
D3 (use with A3; A3.5; A4)	607.01.943	607.01.953
E1 (use with A1; A2; A3; A3.5; A4)	607.01.941	607.01.951

Neo Spectra™ ST Shade Guide	607.01.869
1 shade guide holder	
9 shade guide tabs (A1, A2, A3, A3.5, A4, BW, D1	, D3, E1)

Compules® Tips Gun 1 Compules® Tips Gun 606.65.900

— with – **SphereTEC**[®]

NEO SPECTRA™ ST IS MADE WITH MORE OF WHAT MATTERS TO YOU.

Magnification of a 1

To learn more about Neo Spectra™ ST composite, visit

dentsplysirona.com/NeoSpectraST



Neo Spectra[™] ST Universal Composite Restorative



Neo Spectra[™] ST composite isn't just made with SphereTEC® filler technology, it's made with more of what matters to you.

74% of dentists

efficiency matters

Versatile handling, made even easier to work with.

Thanks to SphereTEC[®] filler technology, Neo Spectra™ ST composite offers optimized handling characteristics in both a high and low viscosity to accommodate your preference and technique for placement efficiency.

- Doesn't stick to hand instruments
- Adapts easily to cavity surfaces
- Easy to sculpt and shape
- Excellent slump resistance

Neo Spectra[™] ST Universal Composite Restorative

When we set out to create Neo Spectra™ ST composite, we asked clinicians what matters most to them. Certainly, it needs to be durable. But could the material handle better to save you time and frustration during placement? And what if shade selection was simplified, minimizing stock and creating efficiencies during the procedure? Could this composite offer outstanding stain resistance to keep patients satisfied? Introducing Neo Spectra[™] ST composite made to deliver what matters to you.

28% of Clas

Five

universa

CLOUD

VITA^{®3} Shades dentsplysirona.com/NeoSpectraSTEffects to learn more.

(Al			A2			A	3
A1	B1	C1	A2	B2	D2	A3	C2	D3
				1		1		F
Ι	Ι	Ι	Ι	I	I	I	I	
A1	B1	C1	A2	B2	D2	A3	C2	D
		A1 B1	A1 B1 C1			A1 B1 C1 A2 B2 D2	A1 B1 C1 A2 B2 D2 A3	A1 B1 C1 A2 B2 D2 A3 C2

